


Scale (2016-)



WARNING

To avoid injury, always read this *Instructions for Use* and accompanied documents before using the product.



Mandatory to read the *Instructions for Use*.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2022.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

Content

Content	3
Foreword	4
Intended Use	5
Safety Instructions	5
Preparations	6
Change the Batteries	6
Scale Overview	7
Parts Designation	8
Product Description / Functions	9
Net Weight / Power On-Off / Scroll button	9
Hold / Power On / Select button	9
Power button for foot scale	9
Battery	9
Net weight	10
Hold weight	10
Limit of Indication	10
Center of Zero	11
Dual Range (not for suspended scale)	11
Tilt	11
Change the Scale Settings	12
Main Menu	12
Navigate through the menu / sub-menu	13
1. Operation	13
1.1 UNIT	13
1.2 PWROFF (Power Off)	14
UP	14
Weigh the patient - Gross / Net Weight	15
Gross weight (6 steps)	15
Net weight (8 steps)	16
Cleaning and Disinfection Instructions	17
Troubleshooting / Error Code	17
Care and Preventive Maintenance	18
Caregiver Obligations	18
Technical Specifications	19
Labels on the Scale	22
Accessories	24
NAWI Declaration of Conformity (Class III)	24
Log Book for Class III Scale	24

Foreword

Thank you for purchasing Arjo equipment.

Please contact us if you have any questions about the operation or maintenance of your Arjo equipment.

Please read and fully understand these *Instructions for Use (IFU)* before using your Scale.

Information in this *IFU* is necessary to perform the proper operation and maintenance of the equipment. It will help to protect your product and make sure that the equipment performs to your satisfaction. The information in this *IFU* is important for your safety and must be read and understood to help prevent possible injury.

Unauthorized modifications on any Arjo equipment can affect safety. Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

Routine maintenance must be performed on your *Scale* to maintain safety and reliability of the equipment. See section “*Care and Preventive Maintenance*” on page 18.

If you require further information, please contact Arjo for comprehensive support and maintenance to maximize the long-term safety, reliability and value of the product.

The telephone numbers appear on the last page of this *IFU*.

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Definitions in this *IFU*

WARNING

Means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

CAUTION

Means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

NOTE

Means: This is important information for the correct use of this system or equipment.



Means: The name and address of the manufacturer.

Intended Use

This equipment must be used in accordance with these safety instructions. Anyone using this equipment must also have read and understood the instructions in these *Instructions for Use (IFU)*.

If there is anything you do not completely understand, ask your Arjo representative.

The Scale is intended to be used with patient lifters to measure the weight of a patient who cannot stand on a scale.

It is designed to be used indoors in hospitals, nursing homes, home care or other health care facilities. It must be operated by trained caregivers with adequate knowledge of the care environment, its common practices and procedures, and in accordance with the guidelines in the *Instructions for Use (IFU)*. Caregivers shall continuously attend the resident/patient during the usage.

The Safe Working Load (SWL) or maximum capacity the scale will vary according to the device on which it will be installed. Refer to the *IFU* of the specific product. The Scale system should only be used for the purpose specified in the *IFU*. Any other use is prohibited.

Patient/Resident assessment

We recommend that facilities establish regular assessment routines. Caregivers should assess each resident/patient according to the following criteria prior to use. The patient/resident's weight must not exceed the safe working of the device on which it will be installed.

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

Installation requirements

The equipment must be installed by appropriately trained personnel according to the *IFU*. The scale may also be pre-installed on the lifter.

Expected service life

The expected service life or maximum period of useful life of this equipment is ten (10) years when subjected to preventive maintenance as specified in the "Care and Preventive Maintenance" section of the *IFU*. The scale system has been designed and tested to achieve up to 10000 cycles with its maximum capacity. The expected service life for consumables such as batteries is dependent upon care and usage of the product concerned.

Safety Instructions

NOTE

Do not touch the chair / stretcher / sling or the patient during the weighing procedure, since it can influence the weighting result.

Preparations

Actions Before First Use (5 steps)

1. Visually check the package for damage. If the product looks damaged due to freight, make a complaint to the transport agency.
Do NOT use the product.
2. The packaging should be recycled according to local regulations.
3. Read this *IFU*.
4. For assembly and installation of scales that are detachable, see *Scale Assembly and Installation Instructions*.
For a scale with quick connect attachment, see the medical device equipment *IFU*.
5. Choose a designated area, where this *IFU* should be kept and is easily accessible at all times.

Actions Before Every Use (2 steps)

1. Visually inspect the scale. If any part is damaged - do NOT use the product.
2. Check the battery level.

In between Use


Clean and disinfect the product according to section “*Cleaning and Disinfection Instructions*” on page 17.

Weight Unit

The default unit of the scale is in kg. To change the weight unit to lb, see section “*1.1 UNIT*” on page 13.

Change the Batteries

Fig. 1




This does not apply to scales attached to Alenti, Miranti and Sara Plus, since they are powered by the product’s main battery.

Suspended and Pillar Scale (3 steps) (See Fig. 1)

1. Use a screw driver to unscrew the cover.
2. Use the screw driver to remove old batteries and replace with four new AA batteries.
3. Attach the cover using the screw driver.

Fig. 2



Foot Scale (4 steps) (See Fig. 2)

1. Push down and slightly forward on the cover.
2. Pull out the battery holder.
3. Insert / replace the four AA batteries.
4. Secure the battery compartment cover into place by pushing it down and back into position.

Scale Overview



Suspended Scale

(Operated by separate batteries)

- *Maxi 500*
- *Maxi Twin*
- *Maxi Sky 2*
- *Maxi Sky 440*
- *Maxi Sky 600*
- *Maxi Sky 1000*
- *Ergolift 600*
- *Voyager*
- *Voyager Duo*
- *Minstrel*
- *Tenor*

NOTE

Attachments and color may vary upon medical device.



Pillar Scale

(Operated by product main battery)

- *Alenti*
- *Miranti*



- *Sara Plus*
(Operated by product main battery)



- *Sara 3000*
(Operated by separate batteries)





Foot Scale

(Operated by separate battery)


- *Bolero*
- *Calypso*

Parts Designation



Product Description / Functions

Fig. 3



Net Weight / Power On-Off / Scroll button

This button has the following functions: (See Fig. 3)

- Power on scale - press once.
- Power off the scale - press the button for 5 sec.
- Scroll in menu mode.
- Set Net weight of patient (Tare).

Fig. 4




Hold / Power On / Select button

This button has the following functions. (See Fig. 4)

- Power on scale - press once.
- Hold the displayed weight.
- Select in menu mode.

Fig. 5



Power button for foot scale

For *Calypso* and *Bolero* hygiene lifts, press the foot button or any scale button to power on the scale.





(See Fig. 5)

Fig. 6



Battery

The battery symbols indicates the battery charge:

-  Full
-  Almost full
-  Half full
-  Low

The low battery symbol will flash when the battery charge is below 10%. (See Fig. 6)

Battery charge below 10% has less than 50 readings left.

The battery operation life time is based on the default settings. To change the default setting (power off time function) see section “1.2 PWROFF (Power Off)” on page 14.

The battery symbol will not show up on the Alenti, Miranti and Sara Plus scale display since they are powered by the lift itself.

Fig. 7




Net weight

This is the Net weight of the patient without any accessories etc. (See Fig. 7)

For instruction how to get the Net weight of the patient, see section “Weigh the patient - Gross / Net Weight” on page 15.

Fig. 8



Hold weight

Memorizes the weight on the scale. (See Fig. 8)


Press the Hold button to lock the weight. The lock symbol and weight will display until the unit powers off.

- If the scale powers off, restart the scale and press the Hold button once and the memorized weight will display.
- If there is more than 10 kg (22 lb) on the scale when pressing the Hold button, the memorized weight will be erased.

NOTE

The patient needs to be still for the Hold weight function to work.

Fig. 9




Limit of Indication

Max weight

The maximum weight symbols appears when the patient weight exceeds the maximum weight capacity of the scale depending on product. (See Fig. 9)


Fig. 10



Minus weight (Below zero)

Indication of weight below zero. This can occur if there is a weight on the scale before it's powered on and then removed after the weighing. (See Fig. 10)

Fig. 11



Center of Zero

Means there is zero weight on the scale / zero on the scale. (See Fig. 11)

The deviation from zero is no more than $\pm 0.25\%$ of the precision.

Fig. 12




Dual Range (not for suspended scale)

The dual range symbol indicates that the weight precision has been changed. This symbol will only appear on the pillar and foot scale. (See Fig. 12)

The precision range starts at 100 g (3.5 oz) and changes to 200 g (7 oz) at a certain weight depending on the product.

Once the 200 g (7 oz) precision becomes active it will remain until the weight is back at 0.0 kg /lb.
The weight precision will restart at 100 g (3.5 oz).

Fig. 13




Tilt

When the scale is tilted, the tilt symbol and text will appear. (See Fig. 13)

When tilted, no weight reading will be available.

Change the Scale Settings


Fig. 14



Main Menu

Enter the main menu to change the scale settings.
To access the main menu, press BOTH the Net Weight and Hold buttons for 3 seconds. (See Fig. 14)

Fig. 15




Here are the menus and sub-menus to choose from:
(See Fig. 15)

- » **1. Operation**
 - 1.1 Unit**
 - 1.2 Power Off**
 - UP** (return to Main menu)
- » **2. Configuration** (only for service technician)
- » **3. Service** (only for service technician)
- » **4. Special** (only for service technician)
- » **EXIT** (exit out of Main menu)

NOTE

Only menu 1 is available for the caregiver. Remaining menus are for service and are password protected.

Fig. 16



Navigate through the menu / sub-menu

- Scroll: press the Net Weight button
- Select: press the Hold button
(See Fig. 16)

Fig. 17



1. Operation

Operation menu has two sub-menus and a return:
(See Fig. 17)

- **1.1 UNIT** - to change unit display.
- **1.2 PWROFF** - to change time delay for power off time.
- **UP** - return to the main menu.

Fig. 18





Fig. 19



1.1 UNIT

NOTE

For a Class III Scale, "kg" is the only option.

Once in the 1.1 UNIT menu, follow the steps below to change the weight unit for the scale.

1. Press Hold button to select and edit current unit setting (weight unit will blink).
2. Press Net Weight button to scroll and switch between kilograms (kg) and pounds (lb).
3. Press Hold button again to confirm selection.
4. Press Net Weight button again to exit the unit menu.

The unit menu gives a display choice between kg and lb. (See Fig. 18) and (See Fig. 19)

Fig. 20



1.2 PWROFF (Power Off)

Once in the 1.2 PWROFF menu, follow the steps below to change the delay time for the scale.


1. Press Hold button to select and edit current time setting (number will blink).
2. Press Net Weight button to scroll through the number of minutes.
3. Press Hold button again to confirm selection.
4. Press Net Weight button again to exit the power off menu.

The time range is from 1-7 minutes. (See Fig. 20)

NOTE

A longer power off time will affect the battery life.

Fig. 21




UP

Scroll and select the UP menu to return to the main menu. (See Fig. 21)

Weigh the patient - Gross / Net Weight

Fig. 22




Gross weight (6 steps)

The weight of the patient including the sling, linens etc.

For the Net weight of the patient see section “Net weight (8 steps)” on page 16.


1. Press any button to power on the scale.
(See Fig. 22)

Fig. 23




2. The display will start of by showing all symbols, followed by the calibration, configuration counter and software version. (See Fig. 23)

Fig. 24



3. 0.0 kg or 0.0 lb will be shown on the display.
(See Fig. 24)

Fig. 25




4. Place the patient safely onto the medical device, equipment Added weight will be the sling, spreader bar etc. (See Fig. 25)

NOTE

By default, there is a 4 minute time window to get the patient into the medical device equipment.

5. Lift the patient.
Suspended scale only: Make sure the patient is in an upright position, to avoid the tilt message on the display..

Fig. 26



6. The scale will now display the Gross weight of the patient. (See Fig. 26)

NOTE

When the weight is flashing it means the weight is not stable.

Fig. 27



Fig. 28




Fig. 29




Fig. 30





Fig. 31



Net weight (8 steps)

The weight of the patient excluding the weight of the sling, spreader bar, linen etc.

For the Gross weight of the patient see section “Gross weight (6 steps)” on page 15.

1. Press any button to power on the scale.
(See Fig. 27)
2. The display will start of by showing all symbols, followed by the calibration, configuration counter and software version. (See Fig. 23)
3. 0.0 kg or 0.0 lb will be shown on the display.
(See Fig. 28)
4. Weigh the equipment / accessories, for example a spreader bar, sling, linen etc.
5. Press the Net weight button. This is the weight at 0.0 kg or lb with the accessories added.
(See Fig. 29)
6. Place the patient safely onto the medical device equipment. (See Fig. 30)
7. Lift the patient.
Suspended scale only: Make sure the patient is in an upright position, to avoid the tilt message on the display.
8. The Net weight of the patient is now displayed.
(See Fig. 31)

NOTE

By default, there is a 4 minute time window to get the patient into the medical device equipment.

NOTE

When the weight is flashing it means the weight is not stable.

Cleaning and Disinfection Instructions

Removing visible residues (2 steps)

1. Use a cloth soaked with water.
2. Remove visible residues from the scale from top to bottom.

General cleaning (5 steps)

1. To clean the scale assembly, spray cleaning fluid on a dampened cloth.
2. Rub the product vigorously to promote effective cleaning.
3. Use a clean dampened cloth soaked in water to wipe off all of the cleaning chemicals.
4. Repeat steps as necessary.
5. Let the parts dry.

Disinfection (3 steps)

1. Rub the scale vigorously using a disinfectant wipe or a cloth soaked in disinfectant.
2. Wipe off all remains of the disinfectant with a wet cloth soaked in water.
3. Wipe the scale dry with a towel or paper.

Allowed cleaning chemicals / Accessories

- Soap
- Cloth
- Towel or paper

CAUTION

Do not use methanol-based or any other type of solvent that may damage the finish or plastic parts.


Allowed Disinfectants

- *Arjo General Purpose Disinfectant*
- *Arjo Disinfectant Cleanser IV*
- *Arjo All Purpose Disinfectant*
- Disinfectant wipes soaked in a 70% v/v solution of isopropyl alcohol

CAUTION

Do not use disinfectant with more than 0.5% w/w hydrogen peroxide as it may damage product labelling.

Troubleshooting / Error Code

Problem	Solution
Scale does not turn on when pressing any button.	<ul style="list-style-type: none">• Replace battery.• Check for external damage. If any, call Arjo for service.
Can't access the main menu.	<ul style="list-style-type: none">• Press both buttons for 3 seconds.
Display shows “-” with no load.	<ol style="list-style-type: none">1. Power off the scale.2. Press the Net Weight button  to reset to zero.
Blinking battery symbol. (Not applicable to Alenti, Miranti, and Sara Plus) .	<ul style="list-style-type: none">• Replace the scale batteries.
Loose bolts or loose top / bottom attachments.	<ul style="list-style-type: none">• Call Arjo for service.
Display shows an error code.	<ul style="list-style-type: none">• Power off the scale. Remove the battery for scales powered by product main battery. See section “Scale Overview” on page 7. If error code won't go away, contact a service technician.
The weight is flashing.	<ul style="list-style-type: none">• Make sure the scale is still.

Care and Preventive Maintenance

SCHEDULE Caregiver Obligations Action/Check	Before Each Use / Between Each Patient Use	Every Month	Every Year
Visually check all parts	X		
Visually check mechanical attachments	X		
Clean/Disinfect	X		
Check battery level		X	
Periodic inspection of attachments		X	
Check for battery leakage			X
Calibration: Class III Scales			X

WARNING

To avoid injury to both the patient and the caregiver, never modify the equipment or use incompatible parts. Only use Arjo designated parts.

Caregiver Obligations

Caregiver obligations shall be carried out by personnel with sufficient knowledge and by following the instructions in this *IFU*.

Before Every Use

Visually check all parts

- **Suspended scale:** Make sure that all the split rings, bolts and Quick-Connect attachments that attach the spreader bar or the strap are secured.
- **All type of scales:** If any part is damaged, contact qualified personnel.

Visually check mechanical attachments

- Inspect the enclosure for any signs of cracking or damages. If cracked or damaged, do NOT use the scale. Contact your Arjo representative for service.

Between Every Patient Use

Clean/Disinfect

- See section “*Cleaning and Disinfection Instructions*” on page 17.

Every Month

Check battery level (Not applicable for *Alenti*, *Miranti* and *Sara Plus*)

- Make sure to replace the batteries if low charge

battery is displayed. See section “*Change the Batteries*” on page 6

Periodic inspection of attachments

- A visual inspections of the scale attachment bolts and brackets is required to make sure they are firmly tightened and with no sign of wear.

Every Year

Check for battery leakage (Not applicable for *Alenti*, *Miranti* and *Sara Plus*)

- Check for corrosion and oxide on the connectors.

Calibration: Class III Scales:

CAUTION

Setting the calibration check interval is the responsibility of the end user and will vary upon state regulations. It's recommended every year if no state regulation rule exists. The Class III scale must be verified, calibrated and sealed by certified personnel only. Contact your local dealer to get more details about this.

- Log Book section retains all the information relevant to the inspection, repair or modification See section “*Log Book for Class III Scale*” on page 24.

Technical Specifications

General	
Safe working load (SWL)	See the lowest SWL of the medical device equipment that the scale is attached to, including sling and spreader bar, if applicable.
Scale weight	0.7 kg / 1.5 lb - Suspended scales except Maxi Sky 1000 and Tenor 0.9 kg / 2 lb - Maxi Sky 1000 and Tenor 5 kg / 11 lb - Alenti and Miranti 8 kg / 18 lb - Sara Plus 9 kg / 20 lb - Bolero, Calypso and Sara 3000
Scale display range - depending on weight and product	100 g (3.5 oz) accuracy <ul style="list-style-type: none"> • 2-100 kg (4-220 lb) range - Sara 3000 and Sara Plus • 2-120 kg (4-265 lb) range - Alenti Miranti, Bolero and Calypso • 2-182 kg (4-400 lb) range - Maxi Twin • 2-190 kg (4-418 lb) range - Minstrel 200 g (7 oz) accuracy <ul style="list-style-type: none"> • 100-190 kg (220-418 lb) range - Sara 3000 and Sara Plus • 120-136 kg (265-300 lb) range - Bolero and Calypso • 120-160 kg (265-352 lb) range - Miranti • 120-182 kg (265-400 lb) range - Alenti • 2-272 kg (4-600 lb) range - Maxi 500, Maxi Sky 2, Maxi Sky 440, Maxi Sky 600, Ergolift 600, Voyager and Voyager Duo • 4-320 kg (9-705 lb) range - Tenor 500 g (1.1 lb) accuracy <ul style="list-style-type: none"> • 10-454 kg (9-1000 lb) - Maxi Sky 1000
Accuracy - depending on product and weight	+/- 100 g (3.5 oz), +/- 200 g (7 oz) and +/- 500 g (1.1 lb)
Resolution	0.1 kg (0.2 lb), 0.2 kg (0.4 lb) or 0.5 kg (1.1 lb)
Battery	Alkaline -1.5 Volt DC (4XAA) (suspended, foot scale and Sara 3000) Lead-acid - 24 Volt (Lifter battery - Alenti, Miranti and Sara Plus)
Protection class	IPx5 - Alenti and Miranti IPx4 - Sara 3000, Sara Plus, Bolero, Calypso and suspended scales
Operating forces on buttons	5N (finger activation) 10N (foot activation only foot scale)
Material	Load bearing components: Stainless steel, Zinc plated steel, aluminum. Enclosure: plastic.

Operating conditions	
Ambient temperature	+10°C to 40°C (50°F to 104°F)
Relative humidity	Max 90% at 20°C (68°F)
Atmospheric pressure	800 hPa - 1060 hPa

Transport and Storage	
Ambient temperature	-40°C to + 70°C (-40°F to 158°F)
Relative humidity	0% - 93% at 20°C (68°F)
Atmospheric pressure	500 hPa - 1060 hPa






Safe Disposal at End of Life	
Battery	All batteries in the product must be recycled separately. Batteries are to be disposed in accordance with national or local regulations.
Package	Corrugated cardboard: recyclable Expanded polystyrene (ESP): recyclable
Product	Electric, aluminium / metal and plastic parts shall be separated and recycled according to marking on the unit.
Electrical and Electronic Components	Lift systems having electrical and electronic components or an electrical cord should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.










Allowed Combinations	
Pillar scale on product:	• Alenti, Miranti, Sara 3000 and Sara Plus
Foot Scale on product:	• Bolero and Calypso

Allowed Combinations		
Suspended scale on products, detachable parts and accessories.	• Tenor	
	• Maxi Twin	• S8559973-24 (Powered DPS)
	• Minstrel	• HMX449-1 (4 hook spreader bar)
	• Maxi 500	<ul style="list-style-type: none"> • 700-19341-BOX (DPS with scale) • 700-05505 (2 hook with scale spreader bar) • 700-05502 (2-point spreader bar w anti-sway & scale) • 700-05503 (2-point spreader bar w anti-sway & scale class III) • 700-19344 (manual DPS w anti-sway & scale) • 700-19345 (manual DPS w anti-sway & scale class III)
	• Maxi Sky 440 and Voyager	<ul style="list-style-type: none"> • 700-05725 (Portable scale adaptor) • 700-05461-BOX (2 hook spreader bar)
	• Maxi Sky 2	<ul style="list-style-type: none"> • 700-19415 (2 hook spreader bar) • 700-19420 (2 hook spreader bar) • 700-19430 (4 hook spreader bar) • 700-19475 (Flat manual DPS) • 700-19480 (Open manual DPS) • 700-19350 (Power DPS Medium) • 700-19355 (Power DPS Large) • 700-15695 (Stretcher frame) • 700.05380 (Morgue spreader bar) • 700.05385 (Morgue spreader bar)
	• Maxi Sky 600	<ul style="list-style-type: none"> • 700-05461-BOX (2 hook spreader bar) • 700-19200 (Flat manual DPS) • 700-19340-BOX (Open manual DPS) • 700-19522 (Stretcher frame) • 700.05380 (Morgue spreader bar) • 700.05385 (Morgue spreader bar)
	• Ergolift 600	• A5465 (Carry Bar)
	• Maxi Sky 1000	<ul style="list-style-type: none"> • 700-05441 (4 hook spreader bar) • 700-05401-BOX (2 hook bar)
	• Voyager Duo	• 700.05491 (2 hook spreader bar)
No other combinations are allowed		

Labels on the Scale

On the labels you can find some of the following symbols / information:

Symbol explanation	
XX-XXXX	Model number.
IPXX	Classification of degrees of protection provided by enclosures.
CALCNT 10XXX	Calibration counter for all other parameters.
CFGCNT 10XXX	Configuration counter mainly for the gravity code.
	CE marking indicating conformity with European Community harmonised legislation. Figures indicate Notified Body supervision.
W1	Weight one in 100g / 3.5oz range.
W2	Weight two in 200g / 7oz range.
d=e e	Weight accuracy based on weight range.
	Indicates that the scale is a non-automatic weighing instrument of accuracy class III.
	Reference number.
	Serial number.
	Maximum and minimum temperature.

	Safe Working Load (SWL).
	Direct Current symbol (Located inside the battery cover).
	Waste Electrical and Electronic Equipment (WEEE) - do not dispose of this product in general household or commercial waste.
	Manufacturer name and address.
	Read the <i>IFU</i> before use.
	Metrology mark, indicating compliance with Directive 2014/31/EU (NAWI) - for scales manufactured after April 20th, 2016. (For Class III scales only) yy = year; XXXX = Notified Body number
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745.
	Unique device identifier
	Non-automatic weighing instrument label Class III; yy = year; XXXX= Approved Body Number.

UKCA Symbol explanation

This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.

**UK
CA**
0086

UK marking indicating conformity with UK Medical Devices Regulations 2002 (SI2002 No 618, as amended)
Figures indicate UK Approved Body supervision.


UK Responsible Person & UK Importer:

Arjo (UK) Ltd, ArjoHuntleigh House, Houghton Regis. LU5 5XF

Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended)

For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.

Accessories



Crash Bracket for Foot scale - Calypso and Bolero

- Part number: CFA0010-90
(Accessory Kit includes bracket, screws, mounting plate and Loctite)

NAWI Declaration of Conformity (Class III)

NOTE

This applies only to Class III scales.

To obtain a copy of the NAWI Declaration of Conformity, contact your local Arjo representative listed on the address page at the end of this *IFU*.

Log Book for Class III Scale

The purpose of the log Book:

- This Log Book section retains all the information relevant to the inspection, repair or modification of the instrument.
- The logbook can be provided by Arjo, the approved inspector or the approved repair engineer.
- The logbook must accompany the scale at the place of use.
- The logbook must always accompany the scale, even if the scale is mobile.

Responsibility of the Holder of the Equipment:

- The holder (end user) must make this Log Book available at all times for the local regulatory authority, the approved inspector or the approved engineer.
- The responsibility for producing the Log Book lies with the holder (end user) of the equipment.
- The holder (end user) is responsible for looking after the integrity of the Log Book.
- The absence of or the damage to the Log Book will result in refusal to sign for the periodic verification or the refusal of an approved engineer to work on the scale.
- If a Log Book replaces another one, the approved inspector or the approved engineer must note this on the new Log Book.
- The holder (end user) is responsible to ensure that the approved bodies for verification fill in the Log Book and should sign the Log Book after the signature of the approved inspector or the approved engineer.

Verification interval (years):	
Name and address of the customer:	
Scale Model	
Serial number of product:	
Identification number of the scale	
Classification of scale	

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	____ ____ ____ ____ / ____ ____ / ____ ____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	____ ____ ____ ____ / ____ ____ / ____ ____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	 ____ ____ ____ ____ / ____ ____ / ____ ____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	 ____ ____ ____ ____ / ____ ____ / ____ ____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

<input type="checkbox"/> Scheduled verification <input type="checkbox"/> Repair (maintenance) - calibration <input type="checkbox"/> Modification	_____ / _____ / _____ Date work completed (YYYY ,MM, DD)
Verification Body Identification of the Organisation – Service provider (name and address)	
Represented by	_____
Signature of Organization	_____
<input type="checkbox"/> Approved	<input type="checkbox"/> Refused

AUSTRALIA

Arjo Australia
Building B, Level 3
11 Talavera Road
Macquarie Park, NSW, 2113,
Australia
Phone: 1800 072 040

BELGIQUE / BELGIË

Arjo Belgium
Evenbroekveld 16
9420 Erpe-Mere
Phone: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.com

BRASIL

Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelice, 329 PB02
Galpão - Lapa
São Paulo – SP – Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA

Arjo Canada Inc.
90 Matheson Boulevard West
Suite 350
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 (905) 238-7880
Free: +1 (800) 665-4831
Fax: +1 (905) 238-7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA

Arjo Czech Republic s.r.o.
Na Strži 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
E-mail: info.cz@arjo.com

DANMARK

Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail: dk_kundeservice@arjo.com

DEUTSCHLAND

Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA

ARJO IBERIA S.L.
Poligono Can Salvatella
c/ Cabanyes 1-7
08210 Barberà del Valles
Barcelona - Spain
Telefono 1: +34 900 921 850
Telefono 2: +34 931 315 999

FRANCE

Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG

Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building, North Tower,
Dubai Science Park,
Al Barsha South
P.O. Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND

Arjo Nederland BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND

Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@Arjo.com

NORGE

Arjo Norway AS
Olaf Helsets vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

ÖSTERREICH

Arjo Austria GmbH
Lemböckgasse 49 / Stiege A / 4.OG
A-1230 Wien
Tel: +43 1 8 66 56
Fax: +43 1 866 56 7000

POLSKA

Arjo Polska Sp. z o.o.
ul. Ks Piotra Wawrzyniaka 2
PL-62-052 KOMORNIKI (Poznań)
Tel: +48 61 662 15 50
Fax: +48 61 662 15 90
E-mail: arjo@arjo.com

PORTUGAL

Arjo em Portugal
MAQUET Portugal, Lda.
(Distribuidor Exclusivo)
Rua Poeta Bocage n.º 2 - 2G
PT-1600-233 Lisboa
Tel: +351 214 189 815
Fax: +351 214 177 413
E-mail: Portugal@arjo.com

SUISSE / SCHWEIZ

Arjo Switzerland AG
Fabrikstrasse 8
Postfach
CH-4614 HÄGENDORF
Tél/Tel: +41 (0) 61 337 97 77
Fax: +41 (0) 61 311 97 42

SUOMI

Arjo Scandinavia AB
Riihitontuntie 7 C
02200 Espoo
Finland
Puh: +358 9 6824 1260
E-mail: Asiakaspalvelu.finland@arjo.com

SVERIGE

Arjo International HQ
Hans Michelsensgatan 10
SE-211 20 MALMÖ
Tel: +46 (0) 10 494 7760
Fax: +46 (0) 10 494 7761
E-mail: kundservice@arjo.com

UNITED KINGDOM

Arjo UK and Ireland
Houghton Hall Park
Houghton Regis
UK-DUNSTABLE LU5 5XF
Tel: +44 (0) 1582 745 700
Fax: +44 (0) 1582 745 745
E-mail: sales.admin@arjo.com

USA

Arjo Inc.
2349 W Lake Street Suite 250
US-Addison, IL 60101
Tel: +1 (630) 307-2756
Free: +1 (800) 323-1245
Fax: +1 (630) 307 6195
E-mail: us.info@arjo.com

JAPAN

Arjo Japan K.K.
東京都港区虎ノ門三丁目7 番8 号
ランディック第2 虎ノ門ビル9 階
Tel: +81 (0)3-6435-6401
Fax: +81 (0)3-6435-6402
E-mail: info.japan@arjo.com

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.



ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo



04.CG.00_EN

CE
2797