	DVT5	≤ 36 cm / 14 inch
	DVT10	≤ 43 cm / 17 inch
	DVT20	≤ 58 cm / 23 inch
	DVT60L	≤ 81 cm / 32 inch
	L501-M	≤ 43 cm / 17 inch
	TRP10	≤ 43 cm / 17 inch
	DVT30	≤ 71 cm / 28 inch
	DVT40	≤ 89 cm / 35 inch
	L503-M	≤ 71 cm / 28 inch
	TRP30	≤ 71 cm / 28 inch
	TRP40	≤ 89 cm / 35 inch

EN To be used only under the direction of a physician • Only for use with Flowtron® DVT prevention pumps manufactured by ArjoHuntleigh • Not made with natural rubber latex • Non-Sterile

Indications

- To help prevent Deep Vein Thrombosis (DVT).

Contraindications

IPC should not be used in the following conditions

- Severe arteriosclerosis or other ischaemic vascular diseases.
- Known or suspected acute Deep Vein Thrombosis (DVT) or phlebitis.
- Severe congestive cardiac failure or any condition where an increase of fluid to the heart may be detrimental.
- Pulmonary Embolism.
- Any local condition in which the garments would interfere, including gangrene, recent skin graft, dermatitis or untreated, infected leg wounds.
- If you are unsure of any contraindications refer to the patient's physician before using the device.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor.

In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Recommendations

- Continuous intermittent pneumatic compression is recommended until the patient is fully ambulatory. Uninterrupted use of the system is recommended.
- In the non-surgical patient, the system should be initiated immediately once the risk of DVT formation is identified.
- Garments should be removed regularly to inspect the skin for signs of redness or pressure points.

Cautions

- Garments should be removed immediately if patient experiences tingling, numbness or pain.
- When used for DVT prophylaxis, continuous use is recommended and any interruption of therapy for a substantial length of time should be at the discretion of physician.
- Lower limb positioning in relation to the garment and tubing should also be considered particularly when a patient is unconscious, or has reduced sensation/ability to move their legs.
- Additional care should be taken when placing the garments on any deformed leg or foot, or on legs with significant oedema.
- This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse and is therefore intended for single patient use. Attempts to clean or sterilize these devices may result in a biocompatibility, infection or product failure risk to the patient.

End of Life Disposal

Garment material or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.

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UK SYMBOL EXPLANATION	
This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.	
	UK marking indicating conformity with UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended)
UK Responsible Person & UK Importer:	
Arjo (UK) Ltd, ArjoHuntleigh House, Houghton Regis. LU5 5XF	
Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended). For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.	