
INSTRUCTIONS FOR USE

Sara Plus Active Slings

EN · ZH · DE

使用说明书 · Bedienungsanleitung

Contents

Foreword	2	Cord Attachment & Detachment	9
Customer Support	2	Attach the Cord (3 steps)	9
Definitions in this IFU	2	Detach the Cord (3 steps)	9
Intended Use	3	Applying the Sling	10
Patient/Resident Assessment	3	In Chair/Wheelchair (17 steps)	10
Expected Service Life	3	Apply the Walking and Standing Sling	11
Pre-sale Information	4	In Chair/Wheelchair (24 steps)	11
Safety Instructions	5	Removing the Sling	13
Safety Practises	5	In Chair/Wheelchair (16 steps)	13
Parts Designation	5	Cleaning and Disinfection	14
Preparations	6	Cleaning Instruction (5 steps)	14
Before First Use (6 steps)	6	Disposable Sling – Do NOT wash	14
Before Every Use (7 steps)	6	Cleaning Chemicals	14
After every use	6	Disinfection	14
Select Sling Size	7	Care and Preventive Maintenance	15
Sling model KKA5120, KKA5370, KKA5130M and		Before and After Every Use	15
MFA4000	7	When Soiled, Stained and Between Patients	15
Sling Selection	8	Storage	15
Transfer	8	Service and Maintenance	15
Allowed Combinations	8	Troubleshooting	16
Safe Working Load (SWL)	8	Technical Specifications	17
Clip Attachment & Detachment	9	Parts and Accessories	17
Attach the Clips (5 steps)	9	Label on the Sling	18
Detach the Clips (3 steps)	9		

EN

WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product. Mandatory to read the Instructions for Use.

Design Policy and Copyright

© and ™ are trademarks belonging to the Arjo group of companies. © Arjo 2022.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

The content of this publication may not be copied either whole or in part without the consent of Arjo.

Foreword

Thank you for purchasing Arjo equipment. Please read this *Instructions For Use (IFU)* thoroughly!

Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

Customer Support

If you require further information, please contact your local Arjo representative. The contact information is located at the end of this *IFU*.

Definitions in this IFU

WARNING

Means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

CAUTION

Means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

NOTE

Means: This is important information for the correct use of this system or equipment.

Intended Use

The Active sling is a product intended for assisted transfer and rehabilitation of patients/residents with limited ability to move.

The Active disposable sling has been designed to support hospital or care facility patients/residents while using a patient/resident transfer lifter, under the supervision of trained nursing staff.

Active disposable slings are to be used for a limited period only, and, by nature of their design, must be treated as a disposable and patient/resident specific product.

The Active disposable sling must not be used to support a patient/resident while bathing or showering.

The Active disposable sling cannot be washed.

The Active sling / Active disposable sling should be used together with Arjo lift devices in accordance with the Allowed Combinations specified in the instructions for use (IFU).

The Active sling / Active disposable sling shall only be used by appropriately trained caregivers with adequate knowledge of the care environment, its common practices and procedures and in accordance with the instructions outlined in the Instructions for use.

The Active sling / Active disposable sling should only be used for the purpose specified in this Instructions for Use. Any other use is prohibited.

Patient/Resident Assessment

We recommend that facilities establish regular assessment routines. Caregivers should assess each patient/resident according to the following criteria prior to use:

- Patient/resident sits in wheelchair.
- Is able to partially bear weight on at least one leg.
- Has some trunk stability.
- Dependent on caregiver in most situations.
- Physically demanding for caregiver.
- Stimulation of remaining abilities is very important.

Safe Working Load (SWL) for Active sling / Active disposable sling:

- KKA5130M (Transfer: 140 kg (309 lbs))
- KKA5130M (Walking and standing), KKA5120, KKA5370: 190 kg (418 lbs);
- MFA4000 (Sara plus disposable standing sling): 190 kg (418 lbs)

The right type and size of slings should be used after proper assessment of each patient/resident's size, condition and the type of lifting situation.

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

Expected Service Life

The expected service life of Active sling / Active disposable sling is the maximum period of useful life.

The expected life of the Active sling / Active disposable sling is dependent on the actual use conditions. Therefore, before use, always make sure that the product does not show any signs of fraying, tearing or other damage and that there is no damage (i.e. cracking, bending, breaking), in which case do not use the sling.

The expected service life:

Service time

- 2 years, Active sling
- 2 weeks (based on 4 transfers per day), Active disposable sling

Shelf time

- 5 years, Active sling / Active disposable sling

Pre-sale Information

General

Intended use, the intended user group and the intended environment; including any contraindications	See <i>"Intended Use"</i> on page 3
Expected lifetime of the product	See <i>"Intended Use"</i> on page 3
Maximum user mass	See <i>"Sling Selection"</i> on page 8
Sizing	See <i>"Select Sling Size"</i> on page 7
Compatible products	See <i>"Allowed Combinations"</i> on page 8
Name and address of the manufacturer	See back page

NOTE

Any organization that purchases hoists and body-support units shall make sure, that the combination(s) are safe either by requiring compatibility documentation for the combination(s) from the manufacturer(s) or by performing compatibility testing themselves hereby transferring the responsibility for a safe combination to the organization.

Safety Instructions

WARNING

To avoid injury, make sure that the patient is not left unattended at any time.

WARNING

To avoid falling, make sure that the user weight is lower than the safe working load for all products or accessories being used.

WARNING

To avoid injury, never wash a disposable sling. A disposable sling is only intended for "single patient use".

WARNING

To avoid injury, never let the patient smoke when using the sling. The sling can catch on fire.

WARNING

To avoid injury, only store the equipment for a short period of time. If stored longer than stated in the *IFU*, it can cause weakness and breakage in the material.

WARNING

To avoid injury, keep the equipment away from sun/UV-light. Exposure to sun/UV-light can weaken the material.

WARNING

To avoid injury, always assess the patient prior to use.

EN


Safety Practises

At any time, if the patient becomes agitated, stop transferring/transporting and safely lower the patient.

Serious incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Parts Designation

1. Waist part
2. Sling label (located on the outside)
3. Attachment loop for cord (mid section)
4. Support belt (hook and loop)
5. Strap (between the legs)
6. Attachment strap
7. Attachment clip
8. Buckle (body)
9. Attachment loop for cord (located on the outside of the leg flap)
10. Buckle (leg)
11. Leg flap
12. Disposable sling identification label (located on the outside)

Example of slings. Not all models are represented.

Preparations

Before First Use (6 steps)

1. Check all parts of the sling, see section “Parts Designation” on page 5. If any part is missing or damaged - Do NOT use the sling.
2. Read this IFU thoroughly.
3. Check that the sling is clean.
4. Choose a designated area where the IFU should be kept and easily accessible at all times.
5. Make sure to have a rescue plan ready in case of a patient emergency.

EN

Before Every Use (7 steps)

1. Patient assessment must always determine the methods used in everyday care, see “Intended Use” on page 3.

WARNING

To avoid the patient from falling, make sure to select the correct sling size according to the IFU.

2. Check the allowed combinations against the sling and lift, see “Allowed Combinations” on page 8.

WARNING

To avoid injury, always make sure to inspect the equipment prior to use.

3. Check all parts of the sling, see “Parts Designation” on page 5. If any part is missing or damaged - Do NOT use the sling. Check for:
 - fraying
 - loose stitching
 - tears
 - fabric holes

6. For questions, contact your local Arjo representative for support.

NOTE

Arjo recommends the name of the patient is written on the sling care label to avoid cross contamination between patients.

- soiled fabric
- damaged clips/buckles
- unreadable or damaged label

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

4. Check that the sling is clean. If the sling is not clean see “Cleaning and Disinfection” on page 14.
5. Always check the washing symbol on the disposable sling’s label. If the disposable sling has been washed, throw the disposable sling away. See “Disposable Sling – Do NOT wash” on page 14.
6.  If this symbol appears, do NOT use the disposable sling. The label indicates that the disposable sling has been washed.
7. Before using the sling, read the lift IFU for transfer, transport and bathing.

After every use

When the patient leaves the hospital, throw away the patient specific disposable sling.


Select Sling Size

Sling model KKA5120, KKA5370, KKA5130M and MFA4000

1. Measure the patient's waist circumference.
2. Follow the sizing chart below to pick the correct size.

EN

		 Waist circumference					
		cm (in)					
		S	M	L	XL	XXL	One size fits most
Sling models	KKA5120, KKA5370, KKA5130M	70-97 (27.5-38)	98-117 (38.5-46)	118-137 (46.5-54)	138-169 (54.5-66.5)	/	/
	MFA4000	/	/	/	/	/	80-125 (31.5-49)

Sling Selection

Transfer

Article number	Safe Working Load	Product Description	Size
KKA5120	190kg (418 lbs)	Sara Plus standing sling	S, M, L, XL
KKA5370	190kg (418 lbs)	Sara Plus standing sling – extra postural support	S, M, L, XL
MFA4000	190kg (418 lbs)	Sara Plus disposable standing sling	One Size
KKA5130M	Transfer: 140 kg (309 lbs) Walking and Standing: 190 kg (418 lbs)	Sara Plus walking sling	S, M, L, XL

EN

Allowed Combinations

WARNING

To avoid injury, always follow the allowed combinations listed in this *IFU*. No other combinations are allowed.

Safe Working Load (SWL)

- Only use the earlier listed slings with the *Sara Plus* standing and raising aid.
- The SWL of the *Sara Plus* is 190 kg (418 lbs)
- Always follow the lowest SWL of the total system. E.g. the *Sara Plus* lift has a SWL of 190 kg (418 lbs) and the KKA5130M sling for transfer only has a SWL of 140 kg (309 lbs). This means that the KKA5130 sling (for transfer only) has the lowest SWL. The patient is not allowed to weigh more than the lowest SWL.

	Standing/Raising Aid	Sara Plus
	SWL	190 kg (418 lbs)
Active Sling	SWL	Size
KKA5120	190 kg (418 lbs)	S, M, L, XL
KKA5130M	140 kg (308 lbs)	S, M, L, XL
KKA5370	190 kg (418 lbs)	S, M, L, XL
MFA4000	190 kg (418 lbs)	One size

Clip Attachment & Detachment


Attach the Clips (5 steps)

1. Place the clip on the spreader bar lug. (See Fig. 1)
2. Pull the strap down.
3. Make sure the lug is locked at the top end of the clip. (See Fig. 2)
4. Make sure the strap is not squeezed in between the clip and the spreader bar.
5. Make sure the straps are not twisted.

Detach the Clips (3 steps)

Make sure that the weight of the patient is taken up by the receiving surface before removing the clip.

1. Pull the strap up. (See Fig. 3)
2. Make sure the lug is unlocked at the bottom of the clip.
3. Remove the clip.


Cord Attachment & Detachment

Attach the Cord (3 steps)

1. Pull the cord with the cone through the attachment loop. (See Fig. 4)
2. Fit the cord into the opening of the cup section (A). (See Fig. 5)
3. Make sure to pull the other end of the cord so that the cone (B) fits into the cup (A) tightly. (See Fig. 6)

Detach the Cord (3 steps)

1. Release the tension in the cord by pushing it towards the cup.
2. Remove the cup from the cone.
3. Pull out the cord from the attachment loop.

Applying the Sling

In Chair/Wheelchair (17 steps)

1. Place the *Sara Plus* lift nearby.
2. Apply the brakes on the wheelchair.
3. Lean the patient forward and place the sling around the patient's lower back. (See Fig. 7)
4. Position the bottom of the sling, aligned horizontally, about five centimetres (two inches) above the patient's waistline.
5. Make sure that the:
 - sling attachment points are forward of the body and have enough reach to attach to the spreader bar,
 - patients arms are outside of the sling,
 - support belt is positioned loosely around the patient's body,
 - support belt is not twisted or trapped behind the patient's back,
 - support belt fits tightly but comfortable for the patient.
6. Fasten the support belt securely by pressing the buckles or hook and loop together. The strap should be tight but comfortable for the patient. (See Fig. 8)

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

Position the *Sara Plus* lift in front of the patient with its protective pad just below the patient's knees. Open up the chassis legs, if needed. See the *Sara Plus IFU*.

7. Apply the brakes on the *Sara Plus*.
8. Place the patient's feet on the *Sara Plus* foot support.
9. Lower the *Sara Plus* Arc-Rest as low as possible, allowing for obstructions e.g chair arms. See the *Sara Plus IFU*.

10. Find the attachment loops on the side of the sling and attach each cord to the loop across from it. The cords shall not cross each other. (See Fig. 9) *Sara Plus disposable sling*: before the transfer of the patient, make sure to use the relevant pair of loop attachment in accordance with the patients position/distance to the lift. It's recommended to use either the innermost or the outermost pair.

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

11. Once the cords are attached, pull and lock them in the retaining cleats (A) on the lift. (See Fig. 10)
12. Pull the cords so that any slack is taken up by each cord and the back section supports the patient comfortably and securely. (See Fig. 11)
13. Ask the patient to hold on to the hand grips with their arms resting on the *Sara Plus Arc-Rest*. See the *Sara Plus IFU*. If the patient is too far from the hand grips, adjust the patient's position by choosing relevant loop combination, see step 11.
14. Slightly lift the patient up to create some tension. See the *Sara Plus IFU*.
15. Make sure the sling is attached correctly and the cords are locked. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clip/cord.
16. Ask the patient to assist during the lifting exercise.

Fig. 7

Fig. 8

Fig. 9

Apply the Walking and Standing Sling

In Chair/Wheelchair (24 steps)


(Only applies to sling KKA5130M)

For transfer/transport refer to *Sara Plus* IFU.

1. Place the *Sara Plus* lift nearby.
2. Apply the brakes on the wheelchair, if any.
3. Lean the patient forward and place the sling around the lower back of the patient.
4. Make sure that:
 - the bottom edge of the sling is levelled with the base of the spine (coccyx),
 - the sling attachment points are forward of the body and have enough reach to attach to the spreader bar,
 - the patients arms are outside of the sling,
 - the support belt is positioned loosely around the patient's body,
 - the support belt is not twisted or trapped behind the patient's back,
 - the support belt fits tightly but comfortable for the patient,
 - each leg flap is placed under the patient's thigh and connected with the support belts and
 - no pieces are twisted under the patient's legs and back.
5. Fasten the support belt securely, by pressing the hook and loop together. The strap should be tight, but comfortable for the patient. (See Fig. 12)
6. Pull the leg flap under the patient's legs. (See Fig. 13)
7. Only for patient's up to 190 kg (420 lbs): Connect each leg buckle on the leg flap to the corresponding body buckle. (See Fig. 14)
8. Position the *Sara Plus* in front of the patient with the front leg support just below the patient's knees. Open up the chassis if needed. See the *Sara Plus* IFU.
9. Apply the brakes on the *Sara Plus*.
10. Place the patient's feet on the *Sara Plus* foot support.
11. Lower the *Sara Plus* Arc-Rest as low as possible, allowing for obstructions e.g chair arms. See the *Sara Plus* IFU.
 - **For patient's up to 140 kg (308 lbs):** Attach each cord opposite from the patient, through the sling attachment loops on the leg flaps (See Fig. 15)
 - **For patient's up to 190 kg (420 lbs):** Attach each cord opposite from the patient, through the sling attachment loops on the mid section. (See Fig. 16)

EN

Continue with the steps on the next page

12. Once the cords are attached, pull and lock the cords in the cleats (A) on the lift. There should be no slack in the cords. (See Fig. 10)

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

13. Attach the clips of the sling to the *Sara Plus* lugs, located under the Arc-Rest. (See Fig. 17) & (See Fig. 2)
14. Make adjustments on the cords so that any slack is taken up by each cord and the back section supports the patient comfortably and securely. (See Fig. 18)
15. Allow the patient to hold on to the hand grips with their arms resting on the *Sara Plus* Arc-Rest. See the *Sara Plus IFU*. (See Fig. 19)
16. Slightly lift the patient up to create some tension. See the *Sara Plus IFU*.

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

17. Make sure the sling is attached correctly and the cords are locked. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clip.
18. Ask the patient to assist during the lifting exercise and to actively stand up.
19. Continue lifting until the patient is in a comfortably supported standing position. See the *Sara Plus IFU*.
20. If any discomfort is experienced by the patient, return to seated position and re-adjust the sling.
21. If a walking practice is to be carried out make sure that the patient is correctly and comfortably supported.
 - Readjust straps and cords so these equally take up any slack
 - Readjust Arc-Rest, if necessary.
22. Release the brakes of the *Sara Plus*.
23. Move the *Sara Plus* away from the chair and close the chassis legs if opened. See the *Sara Plus IFU*.
24. Transport the patient or continue with walking practice by removing the foot stand. See the *Sara Plus IFU*.

EN

Fig. 17


Fig. 18



Fig. 19

Removing the Sling

In Chair/Wheelchair (16 steps)

1. Apply the brakes on the wheelchair.
2. Add the foot support for the patient if he/she has been practising walking.
3. Position the patient with his/her back against of the chair/wheelchair. If needed, open up the chassis legs on the *Sara Plus* lift.
4. Apply the brakes on the *Sara Plus* lift.
5. Lower the patient down into the chair/wheelchair. See the *Sara Plus* IFU. (See Fig. 20)
6. Make sure the patient's lower back is positioned all the way in against the chair/wheelchair and that the weight of the patient is taken up by the receiving surface before removing the clip/loop.
7. Detach the sling clips from the *Sara Plus* lift, if any. (See Fig. 3)
8. Unlock the cords from the cleats (A). (See Fig. 21)
9. Unlock and remove the cord from the sling loop. (See Fig. 22)
10. Remove the patient's feet from the foot support.
11. Release the brake on the *Sara Plus* lift.
12. Move the *Sara Plus* lift away from the patient.
13. Unbuckle the leg straps from the body straps, if any. (See Fig. 23)
14. Open up the support belt, if any. (See Fig. 24)
15. Remove the leg flaps, if any, from underneath the patient's legs.
16. Lean the patient forward and remove the sling from behind the patient.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

EN

Fig. 20


Fig. 21

Fig. 22


Fig. 23

Fig. 24

Cleaning and Disinfection

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

WARNING

To avoid material damage and injury, clean and disinfect according to this IFU:

- No other chemicals are allowed.
- Never clean with chlorine.
- Chlorine will deteriorate the surface of the material.

WARNING

To avoid injury, always remove the sling before disinfecting the lift.

WARNING

To avoid injury, never wash a disposable sling. A disposable sling is only intended for “single patient use”.

All Sara Plus active slings, except the disposable slings, should be cleaned when it is soiled or stained and between patients.

Cleaning Instruction (5 steps)

Sara Plus active slings (NOT disposable slings)

1. Detach the sling from the lift.
2. Close all buckles and hook and loop attachment on the slings and belts prior to washing.
3. Check the symbols on the product label for cleaning recommendations.
4. Machine wash the sling at 70°C (158°F). In order to obtain sufficient level of disinfection. If allowed according to the product label, tumble dry in low temperature, max 60°C (140°F).
5. **Do NOT**
 - wash with other items that have rough surfaces or sharp objects
 - use any mechanical pressure, pressing or rolling during the washing and drying process
 - use bleach
 - use gas sterilization
 - use autoclave
 - dry clean
 - steam
 - ironing

Disposable Sling – Do NOT wash

The disposable sling is for single patient use. Do NOT clean, wash, disinfect, wipe or sterilize the disposable sling. If the disposable sling has been subjected to any such treatment, it shall be discarded. The disposable sling is marked with the “Do NOT wash” symbol. (See Fig. 25)
If the disposable sling is washed, the “Do NOT use/lift” symbol will appear. (See Fig. 26)

Cleaning Chemicals

For all slings (NOT the disposable slings) use common commercial detergent without optical brightener. No other chemicals are allowed e.g. chlorine, softener, iodine based disinfectants, bromine and ozone.

Disinfection

The only allowed disinfection is disinfection by washing.

This does not apply to the disposable slings. Do NOT clean, wash or disinfect a disposable sling.

Fig. 25

Fig. 26

Care and Preventive Maintenance

WARNING

To avoid injury to both patient and caregiver, never modify the equipment or use incompatible parts.

Before and After Every Use

Visually check all exposed parts

The caregiver shall inspect the sling before and after every use. The complete sling should be checked for all deviations listed below. If any of these deviations are visible, replace the sling immediately.

- Fraying
- Loose stitching
- Tears
- Holes
- Discolouration or stains from bleaching
- Sling soiled or stained
- Unreadable or damaged label
- For all slings (not disposable slings): If the sling is soiled and stained – wash the sling according to section Cleaning and Disinfection. If the soil and stains are still visible, replace the sling immediately.
- For disposable slings: Make sure the “Do NOT wash” symbol on the head support is intact. If damaged or missing and the underlying red “Do NOT use/lift” symbol appears, replace the sling immediately. See “*Label on the Sling*” on page 18.

When Soiled, Stained and Between Patients

Clean/Disinfect, all slings except the disposable slings

The caregiver should make sure that the sling is cleaned according to “*Cleaning and Disinfection*” on page 14, when it is soiled or stained and between patients.

Disposable slings

Washing or disinfection is not allowed. Do not wipe, sterilize or soil the disposable slings. If a disposable sling has been subjected to any treatment, it shall be discarded.

Storage

When not in use, the slings should be stored away from direct sunlight where they are not subject to unnecessary strain, stress or pressure, or to excessive heat or humidity. The slings should be kept away from sharp edges, corrosives or other things that could cause damage on the sling.

Service and Maintenance

It is recommended that qualified personnel inspect the condition of the sling twice a year (every 6 month), according to ISO 10535.

EN

Troubleshooting

Problem	Action
The patient is not positioned correctly in the sling.	<ul style="list-style-type: none"> • Make sure the clip/lift cords are attached at the right length. • Make sure the straps are not twisted or that there are no knots on the lift cords. • Make sure the patient is positioned on the inside of the sling. The outside has a label for reference. • Make sure the patient is centred in the sling.
The patient feels discomfort in the leg area when seated in the sling.	<ul style="list-style-type: none"> • Make sure that there are no creases on the leg flaps of the sling. • Make sure the leg flaps are attached at the right length.
It is difficult to apply the clips.	Make sure the clip is not caught between the lug on the lift and the clip attachment.
It is difficult to use the sling in combination with the lift.	Use the sling only in accordance with the combinations described in section <i>"Allowed Combinations"</i> on page 8.
The position of the patient in the sling should be more reclined or seated.	Use the clips/lift cords to adjust the patient's position.
The patient wishes to have more support for the back.	Use the inner clips for attachment or tighten the lift cords to adjust the patient's position.
The patient feels discomfort in the waist area.	Adjust the support belt, tighten or loosen.
The sling is sliding up under the patient's arms.	Tighten the support belt for a more comfortable fit.
The patient feels discomfort while using the sling.	<ul style="list-style-type: none"> • Make sure the patient is holding on to the handles on the lift. • Adjust the support belt for a more comfortable fit.
The disposable sling is soiled, stained or wet.	Throw the disposable sling away and replace it with a new one.
The "Do NOT wash" label is damaged or missing and the underlying red label "Do NOT use/lift" appears.	Throw the disposable sling away and replace it with a new one.

EN

Technical Specifications

General	
Safe working load (SWL) = Maximum total load	See "Allowed Combinations" on page 8
Service life - Recommended period of use	Disposable sling: 2 weeks, based on 4 transfers/day* All other slings: 2 years* * Refer to "Expected Service Life" on page 3.
Shelf life - Maximum storage period of new unpacked product	5 years
Model and Type	See "Sling Selection" on page 8

EN

Operating, Transport and Storage Environment	
Temperature	0°C to +40°C (+32°F to +104°F) Operating and Storage -25°C to +60°C (-13°F to +140°F) Transport
Humidity	Max 15–70% at +20°C (68°F) Operating and Storage Max 10–95% at +20°C (68°F) Transport

End of Life Disposal	
Package	The bag consists of soft plastic, recyclable to local legislations.
Sling	Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.

Parts and Accessories

Sling (X is the size)	Part Description	Article number (X is the size)
KKA5120-X	Belt	KKX52810-X
KKA5130-X	Belt	KKX52810-X
KKA5370-X	Belt	KKX52810-X

Label on the Sling

EN

Care and Washing Symbols	
	Machine wash 70°C (158°F)
	No bleaching allowed
	No tumble drying allowed
	Tumble drying Does NOT apply to disposable slings
	No ironing allowed
	Do NOT wash the disposable sling. Only applies to disposable slings.
	No dry cleaning allowed

Certificates/Markings	
	CE marking indicating conformity with European Community harmonised legislation
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745
	Unique device identifier

Disposable Sling Symbols	
	“Do NOT wash” symbol. Located on the outside of the sling.
	“Do NOT use/lift” symbol. The disposable sling has been washed. Located on the outside of the sling.

Symbol for Clip/Cord	
	Use a clip and/or cord spreader bar.

Article Number	
REF XXXXXX-X	Article number with -X refers to the sling size.
REF XXXXXXX	If the article number is stated without the size letter at the end, it means the sling is a One Size.

Fibre Content	
PES	Polyester
PE	Polyethylene
PU	Polyurethane
PA	Polyamide

Misc. Symbols	
	Safe Working Load (SWL)
	Patient's name symbol
	Single patient multiple use
	Record symbol
	Read the <i>IFU</i> before use
	Article number
	Manufacturing date and year
	Manufacturer name and address

UK Symbol explanation	
This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.	
	UK marking indicating conformity with UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended)
<p>UK Responsible Person & UK Importer: Arjo (UK) Ltd, ArjoHuntleigh House, Houghton Regis. LU5 5XF Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended). For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.</p>	

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Intentionally left blank

AUSTRALIA

Arjo Australia
 Building B, Level 3
 11 Talavera Road
 Macquarie Park, NSW, 2113, Australia
 Phone: 1800 072 040

BELGIQUE / BELGIË

Arjo Belgium
 Evenbroekveld 16
 9420 Erpe-Mere
 Phone: +32 (0) 53 60 73 80
 Fax: +32 (0) 53 60 73 81
 E-mail: info.belgium@arjo.com

BRASIL

Arjo Brasil Equipamentos Médicos Ltda
 Rua Marina Ciufuli Zanfelice, 329 PB02
 Galpão - Lapa
 São Paulo - SP - Brasil
 CEP: 05040-000
 Phone: 55-11-3588-5088
 E-mail: vendas.latam@arjo.com
 E-mail: servicios.latam@arjo.com

CANADA

Arjo Canada Inc.
 90 Matheson Boulevard West
 Suite 350
 CA-MISSISSAUGA, ON, L5R 3R3
 Tel/Tél: +1 (905) 238-7880
 Free: +1 (800) 665-4831
 Fax: +1 (905) 238-7881
 E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA

Arjo Czech Republic s.r.o.
 Na Strži 1702/65
 140 00 Praha
 Czech Republic
 Phone No: +420225092307
 E-mail: info.cz@arjo.com

DANMARK

Arjo A/S
 Vassingerødvej 52
 DK-3540 LYNGE
 Tel: +45 49 13 84 86
 Fax: +45 49 13 84 87
 E-mail: dk_kundeservice@arjo.com

DEUTSCHLAND

Arjo GmbH
 Peter-Sander-Strasse 10
 DE-55252 MAINZ-KASTEL
 Tel: +49 (0) 6134 186 0
 Fax: +49 (0) 6134 186 160
 E-mail: info-de@arjo.com

ESPAÑA

ARJO IBERIA S.L.
 Poligono Can Salvatella
 c/ Cabanyes 1-7
 08210 Barberà del Valles
 Barcelona - Spain
 Telefono 1: +34 900 921 850
 Telefono 2: +34 931 315 999

FRANCE

Arjo SAS
 2 Avenue Alcide de Gasperi
 CS 70133
 FR-59436 RONCQ CEDEX
 Tél: +33 (0) 3 20 28 13 13
 Fax: +33 (0) 3 20 28 13 14
 E-mail: info.france@arjo.com

HONG KONG

Arjo Hong Kong Limited
 Room 411-414, 4/F, Manhattan Centre,
 8 Kwai Cheong Road, Kwai Chung,
 N.T.,

HONG KONG

Arjo Hong Kong
 Tel: +852 2960 7600
 Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
 Via Giacomo Peroni 400-402
 IT-00131 ROMA
 Tel: +39 (0) 6 87426211
 Fax: +39 (0) 6 87426222
 E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
 Office 908, 9th Floor,
 HQ Building, North Tower,
 Dubai Science Park,
 Al Barsha South
 P.O. Box 11488, Dubai,
 United Arab Emirates
 Direct +971 487 48053
 Fax +971 487 48072
 Email: Info.ME@arjo.com

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At Arjo, we believe that empowering movement within healthcare environments is essential to quality care. Our products and solutions are designed to promote a safe and dignified experience through patient handling, medical beds, personal hygiene, disinfection, diagnostics, and the prevention of pressure injuries and venous thromboembolism. With over 6500 people worldwide and 65 years caring for patients and healthcare professionals, we are committed to driving healthier outcomes for people facing mobility challenges.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

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Arjo Inc.