INSTRUCTIONS FOR USE

Sara Plus wipeable standing sling

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使用说明书 · Bedienungsanleitung · Mode d'emploi · Istruzioni per l'uso · Gebruiksaanwijzing

Contents

Foreword	3	Cord Attachment & Detachment.....	7
Customer Support.....	3	Attach the Cord (3 steps)	7
Definitions in this IFU	3	Detach the Cord (3 steps)	7
Intended Use	3	Applying the Sling.....	8
Patient/Resident Assessment	4	In Chair/Wheelchair (17 steps).....	8
Expected Service Life	4	Removing the Sling	9
Safety Instructions	4	In Chair or Wheelchair (11 steps).....	9
Safety Practises	4	Cleaning and Disinfection.....	9
Parts Designation	5	Disinfection Agents.....	9
Preparations	5	Allowed Disinfectant Agents for Wiping.....	10
Before First Use (6 steps)	5	Care and Preventive Maintenance	11
Before Every Use (5 steps)	6	Before and After Every Use.....	11
After every use	6	When Soiled and Stained as well as	
Sling Model Selection	6	Between Patients	11
Select Sling Size.....	6	Storage.....	11
Measure (2 steps)	6	Service and Maintenance.....	11
Allowed Combinations	7	Troubleshooting	11
Safe Working Load (SWL).	7	Technical Specifications	12
		Labels on the Sling	13

WARNING: To avoid injury, always read this Instructions For Use and accompanied documents before using the product.

Mandatory to read the Instructions For Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2023.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

Foreword

Thank you for purchasing Arjo equipment. Please read this Instructions For Use (IFU) thoroughly! Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Customer Support

If you require further information, please contact your local Arjo representative. The contact information is located at the end of this IFU.

Definitions in this IFU

WARNING

Warning means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

CAUTION

Caution means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

NOTE

Note means: This is important information for the correct use of this system or equipment.

Intended Use

This sling is intended to be used together with a lifting device to provide body support during raising to a standing position and short transfers of adult patients/residents. It is intended to be cleaned by means of wiping.

The Sara Plus wipeable standing sling should be used together with Arjo lift devices (Sara Plus) in accordance with the Allowed Combinations which are specified in the Instructions for Use.

The Sara Plus wipeable standing sling shall only be used by appropriately trained caregivers with adequate knowledge of the care environment, its common practices and procedures, and in accordance with the guidelines in the Instructions for Use.

The Sara Plus wipeable standing sling is intended to be used in hospitals, nursing homes, assisted living and group homes, and for home care.

The Sara Plus wipeable standing sling should only be used for the purpose specified in the Instructions for Use. Any other use is prohibited.

Patient/Resident Assessment

Arjo recommends that facilities establish regular assessment routines.

Prior to the transfer, caregivers should assess each patient/resident according to the following criteria:

- Size
- Mobility level
- Medical conditions
- Ability to understand and respond to instructions

The Sara Plus wipeable standing sling is intended for patients/residents:

- who are able to bear weight on at least one leg and have some trunk stability.
- whose weight is no more than 200 kg (440 lb).
- has a length in the range of 140 to 200 cm (55 1/8" - 78 3/4").

If the patient/resident does not meet these criteria, an alternative equipment/system shall be used.

Expected Service Life

The expected service life of Sara Plus wipeable standing sling is the maximum period of useful life.

The expected service life of this sling is two (2) years, subject to preventive maintenance being carried out in accordance with the instructions for care and maintenance found in the Instructions for use.

The shelf life, i.e., the period before the first use of the device is five (5) years.

Article Number	Product Description	Service Life	Shelf Life
KKA6000	Sara Plus wipeable standing sling	2 years, based on 12 transfers per day	5 years

Safety Instructions

WARNING

To avoid injury, make sure that the patient is not left unattended at any time when in sling.

WARNING

To avoid falling, make sure that the patient weight is lower than the safe working load for all products or accessories being used.

WARNING

To avoid injury, never let the patient smoke when using the sling. The sling can catch on fire.

WARNING

To avoid injury, only store the sling for a short period of time. If stored longer than stated in the *IFU*, it may cause weakness and breakage in the material.

WARNING

To avoid injury, keep the sling away from sun/UV-light. Exposure to sun/UV-light may weaken the material.

WARNING


To avoid injury, always assess the patient prior to use.

Safety Practises


At any time, if the patient becomes agitated, stop transferring/transporting and safely lower the patient.

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Parts Designation

1. Safety belt
2. Buckle
3. Attachment hook (C-Hook) for Sara Plus cord
4. Flap fold
5. Sling label/print

Preparations

Before First Use (6 steps)

1. Check all parts of the sling, see section "Parts Designation" on page 5. If any part is missing or damaged - Do NOT use the sling.
2. Read this IFU thoroughly.
3. Check that the sling is clean.
4. Choose a designated area where the IFU should be kept and easily accessible at all times.

5. Make sure to have a rescue plan ready in case of a patient emergency.
6. For questions, contact your local Arjo representative for support.

NOTE

To keep track of the service life of the sling, document the serial number and the date it was first used.

Before Every Use (5 steps)

- Patient assessment must always determine the methods used in everyday care, see "Intended Use" on page 3.

WARNING

To avoid the patient from falling, make sure to select the correct sling size according to the IFU.

- Check the allowed combinations against the sling and lift, see "Allowed Combinations" on page 7.

WARNING

To avoid injury, always make sure to inspect the sling prior to use. If any part is missing or damaged - Do NOT use the sling.

- Check all parts of the sling, see "Parts Designation" on page 5. If any part is missing or damaged - Do NOT use the sling. Check for:
 - fraying
 - loose stitching/delamination

- cracks
- tears
- fabric holes
- soiled fabric
- damaged buckles/c-hooks
- unreadable or damaged label/print

WARNING

To prevent cross-contamination, always follow the cleaning and disinfection instructions in this IFU.

- Check that the sling is clean. If the sling is not clean, clean and disinfect it according to the "Cleaning and Disinfection" section on p. 9.
- Before using the sling, read the lift IFU for transfer and transport.

After every use

- Clean and disinfect the sling according to the instruction in this IFU on page 9.
- Check all parts of the sling, see "Parts Designation" on page 5. If any part is missing or damaged, dispose the sling.

Sling Model Selection

Article Number	Safe Working Load (SWL)	Product Description	Size
KKA6000	200 kg (440 lb)	Sara Plus wipeable standing sling	S, M, L, XL

Select Sling Size

Measure (2 steps)

The patient's physical disabilities, weight distribution and general physique need to be taken into consideration when selecting a sling.

- Measure the patient's waist in seated position.
- Select a sling size that is closest to that measurement. In case of doubt choose the smallest one. Each sling size has its own colour, which is found on the print outside of the sling.

For reference, align the flap fold with patient's side (See Fig. 1). If sling is too small, this flap fold will be too far back (See Fig. 2). If sling is too large, the flap fold will be further in the front (See Fig. 3).

Article Number	Product Description	Waist Size in cm (in)			
		S	M	L	XL
KKA6000	Sara Plus wipeable standing sling	64-84 cm (25" - 33")	76-104 cm (30" - 41")	88-124 cm (34½" - 49")	108-160 cm (42½" - 63")

Flap fold

Patient's side

Allowed Combinations

⚠️ WARNING

To avoid injury, always follow the allowed combinations listed in this *IFU*. No other combinations are allowed.

E.g. the Sara Plus lift has a SWL of 190 kg (420 lb) and the KKA6000 sling has a SWL of 200 kg (440 lb). This means that the Sara Plus lift has the lowest SWL. The patient must not weigh more than the lowest SWL.

Safe Working Load (SWL)

The SWL of the Sara Plus is 190 kg (420 lb). Always follow the lowest SWL of the total system.

⚠️ WARNING

To avoid injury, always follow the lowest SWL in the complete system (lift or sling).

Sara Plus

Article Number	Sling SWL	Sara Plus SWL	Combination SWL	Compatible Size
KKA6000	200 kg (440 lb)	190 kg (420 lb)	190 kg (420 lb)	S, M, L, XL




Cord Attachment & Detachment

Attach the Cord (3 steps)

1. Pull the cord with the cone through the C-Hook. (See Fig. 4)
2. Fit the cord into the cup section opening (A). (See Fig. 5)
3. Pull the other end of the cord so that the cone (B) fits into the cup (A) tightly. (See Fig. 6)

Detach the Cord (3 steps)

1. Release the tension in the cord by pushing it towards the cup.
2. Remove the cup from the cone.
3. Pull out the cord from the C-Hook.

Applying the Sling



In Chair/Wheelchair (17 steps)

1. Place the Sara Plus lift nearby.
2. Apply the brakes on the wheelchair.
3. Ask or assist the patient to lean forward and place sling around the patient's lower back. (See Fig. 7) Make sure that the sling is not upside down and that the patient is positioned on the inside of the sling (soft lilac). For reference, the outward side has a print label.
4. Center the sling around the patient's waistline.
5. Ensure that:
 - patient's arms are outside the sling,
 - safety belt is positioned loosely around the patient's body and
 - safety belt is not twisted or trapped behind the patient's back.
6. Fasten the safety belt securely by fastening the buckle. The safety belt should be tight but comfortable for the patient. Use both hands to adjust the safety belt, one to feed the belt through the buckle and the second one to pull on the belt end. (See Fig. 8)

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the Sara Plus Arc-Rest.

7. Position the Sara Plus lift in front of the patient with its protective pad just below the patient's knees. Open up the chassis legs, if needed. See the Sara Plus IFU.
8. Place patient's feet on Sara Plus foot support.
9. Apply the brakes on the Sara plus.
10. Lower the Sara Plus Arc-Rest as low as

- possible, allowing for obstructions e.g chair arms. See the Sara Plus IFU.
11. Find the C-Hook on the side of the sling and attach each cord to the C-Hook across from it. The cords should not cross each other. (See Fig. 9)

WARNING

To avoid the patient from falling, ensure that the sling attachments are attached securely before and during the lifting process.

12. Once the cords are attached, pull and lock them in the retaining cleats (A) on the lift. (See Fig. 10)
13. Ask or assist the patient to lean forward and hold on to the hand grips with their arms resting on the Sara Plus Arc-Rest. See the Sara Plus IFU.
14. Pull the cords so that any slack is taken up by each cord and the back section supports the patient comfortably and securely. (See Fig. 11)
15. Slightly lift the patient up to create some tension. See the Sara Plus IFU.
16. Make sure the sling is attached correctly and the cords are locked. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the cord.
17. Encourage the patient to assist during the raising exercise by holding on to the support grips with both hands. Continue raising until the patient is in a comfortable standing position. Make sure the patient leans back against the sling at all times. See the Sara Plus IFU.

Removing the Sling

In Chair or Wheelchair (11 steps)

1. Apply the brakes on the wheelchair.
2. Position the patient with his/her back against the backrest of the chair/wheelchair. If needed, open up the chassis legs on the Sara Plus lift.
3. Apply the brakes on the Sara Plus lift.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the Sara Plus Arc-Rest.

4. Lower the patient down into the chair/wheelchair. See the Sara Plus IFU. (See Fig. 12)

Fig. 12

Fig. 13

5. Make sure the patient's lower back is positioned all the way in against the chair/wheelchair and that the weight of the patient is taken up by the receiving surface before unlocking the cords from the cleats (A). (See Fig. 13)
6. Unlock and remove the cord from the C-Hook. (See Fig. 14)
7. Remove the patient's feet from the foot support.
8. Release the brake on the Sara Plus lift.
9. Move the Sara Plus lift away from the patient.
10. Open up the safety belt. (See Fig. 15)
11. Ask or assist the patient to lean forward and remove the sling from behind the patient.

Fig. 14


Fig. 15

Cleaning and Disinfection

WARNING

To prevent cross-contamination between patients, it is the responsibility of the caregivers to clean and disinfect the sling according to the cleaning and disinfection instructions in this IFU.

WARNING

To avoid material damage and injury, clean and disinfect according to this IFU:

- No other methods or chemicals are allowed.

WARNING

To avoid injury, always detach the sling from the lift before cleaning, disinfecting.

Disinfection Agents

When soiled or stained, the Sara Plus wipeable standing sling should be cleaned and disinfected using the methods prescribed in the following table (on the next page). The sling should always be cleaned and disinfected between patients. To avoid injury, only use the allowed cleaning and disinfecting methods.

(continued on next page)

Article Number	Product Description	Cleaning/Disinfection	Drying
KKA6000	<p>Sara Plus wipeable standing sling</p> 	<p>The only allowed disinfection method for wipeable slings is disinfection by wiping. DO NOT WASH.</p> <ol style="list-style-type: none"> 1. Detach the sling from the lift and patient. 2. If needed, rinse to remove any deposits or visible dirt from the sling. 3. Wipe all parts of the sling using a solution containing one of the compatible disinfectants from the "Allowed Disinfection Agents for wiping" table below. 4. Disinfection methods should comply with the local or national guidelines (Decontamination of Medical Devices), depending on the Healthcare Facility or country of use. If uncertain, consult local Infection Control specialist. 5. Allow contact time according to the cleaning agent manufacturer's instruction for use. 6. If stated in the cleaning agent manufacturer's instruction for use, rinse to remove chemical residue. 7. Remember: <ul style="list-style-type: none"> • DO NOT launder • DO NOT use gas sterilization • DO NOT use any mechanical pressure • DO NOT tumble dry • DO NOT use autoclave • DO NOT dry clean • DO NOT steam • DO NOT iron 	Air dry at room temperature

Allowed Disinfectant Agents for Wiping	Sara Plus wipeable standing sling (KKA6000)
Quaternary Ammonium	$\leq 0.28\%$ (2,800 ppm) X
Isopropyl Alcohol	$\leq 70\%$ X
Phenolic	$\leq 1.56\%$ (15,600 ppm) X
Bleach (Sodium Hypochlorite)	$\leq 1\%$ (10,000 ppm) X*
Peracetic Acid	$\leq 0.25\%$ (2,500 ppm) X
Hydrogen Peroxide	$\leq 1.5\%$ X

* Sling may show some signs of discolouration and become slightly sticky after repeated contact with bleach.

Care and Preventive Maintenance

WARNING

To avoid injury to both patient and caregiver, never modify the equipment or use incompatible parts.

Before and After Every Use

Visually check all exposed parts

The caregiver shall inspect the sling before and after every use. The complete sling should be checked for all deviations listed below. If any of these deviations are visible, replace the sling immediately.

- Fraying
- Loose stitching/delamination
- Tears
- Cracks
- Holes
- Major discolouration or stains from bleaching
- Unreadable or damaged label/print
- If the sling is soiled and stained - wash or disinfect the sling according to section Cleaning and Disinfection. If any stains are still visible, replace the sling immediately.

When Soiled and Stained as well as Between Patients

Clean/Disinfect

The caregiver should make sure that the sling is cleaned and disinfected when soiled or stained. The caregiver should make sure that the sling is cleaned and disinfected between patients. The sling should be cleaned and disinfected according to the "Cleaning and Disinfection" section on page 9.

Storage

When not in use, the slings should be stored away from direct sunlight where they are not subject to unnecessary strain, stress or pressure, or excessive heat or humidity. The slings should be kept away from sharp edges, corrosives or other things that could cause damage on the sling.

Service and Maintenance

It is recommended that qualified personnel inspect the condition of the sling twice a year (every 6 months), according to ISO 10535.

Disposal of product and components should be compliant with local regulations.

Troubleshooting

Problem	Action
Patient is not correctly positioned in the sling.	<ul style="list-style-type: none"> • Make sure the lift cords are attached at the right length. • Make sure that there are no knots on the lift cords. • Make sure the patient is positioned on the inside of the sling. The outside has a label or print for reference. • Make sure the patient is centred in the sling. • Make sure that the sling is not upside down. The label or print should be readable.
It is difficult to use the sling in combination with the lift.	Use the sling only in accordance with the combinations described in section " <i>Allowed Combinations</i> " on page 7.
The position of the patient in the sling should be more reclined or seated.	Use the lift cords to adjust the patient's position.
Patient wishes to have more back support.	Tighten the lift cords to adjust the patient's position.

(continued on next page)

Patient feels discomfort in the waist area.	Adjust the safety belt, tighten or loosen.
Sling is sliding up under the patient's arms.	<ul style="list-style-type: none"> • Ensure the proper sling size is used • Tighten the safety belt for a more comfortable fit • Ensure proper clinical assessment.
Patient feels discomfort while using the sling.	<ul style="list-style-type: none"> • Make sure the patient is holding on to the handles on the lift. • Adjust the safety belt for a more comfortable fit.
Patient cannot fully reach standing position.	Select another sling size or adjust the cord length and Sara Plus arc-rest height.

Technical Specifications

General	
Safe working load (SWL) = Maximum total load	See <i>Allowed Combinations</i> on page 7
Service life - Recommended period of use	2 years, based on 12 transfers per day* * refer to <i>Expected Service Life</i> on page 4
Shelf life - Maximum storage period of new unpacked product, starting from the manufacturing date	5 years
Model and Type	See <i>Sling Model Selection</i> on page 6
Operating, Transport and Storage Environment	
Temperature	0°C to +40°C (+32°F to +104°F) Operating and Storage -25°C to +60°C (-13°F to +140°F) Transport
Humidity	Max 15–70% at +20°C (68°F) Operating and Storage Max 10–95% at +20°C (68°F) Transport
End of Life Disposal (device should be recycled according to national regulations)	
Package	The box consists of cardboard, recyclable to local legislations. The bag consists of soft plastic, recyclable to local legislations.
Sling	Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.

Labels on the Sling

Care and Washing Symbols		Article Number
	No tumble drying allowed	REF XXXXXX-X Article number with -X refers to the sling size.
	No ironing allowed	
	Do NOT wash	
	No dry cleaning allowed	
	Wipe clean	

Certificates/Markings		Misc. Symbols
	CE marking indicating conformity with European Community harmonised legislation	Safe Working Load (SWL)
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745	Read the Instructions for Use before using
	Unique device identifier	REF Article number
		SN Serial number
		Manufacturing date and year
		Manufacturer name and address
		Size of the sling illustrated by colour and letter. S for Small, M for Medium, L for large, and XL for extra large.

Symbol for C-Hook	
	Use a C-Hook attachment mean.

UK Symbol explanation	
This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.	
	UK marking indicating conformity with UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended)
UK Responsible Person & UK Importer:	

Arjo (UK) Ltd, ArjoHuntleigh House, Houghton Regis, LU5 5XF

Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended).

For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.

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Intentionally left blank

AUSTRALIA
Arjo Australia
Building B, Level 3
11 Talavera Road
Macquarie Park, NSW, 2113, Australia
Phone: 1800 072 040

BELGIQUE / BELGIË
Arjo Belgium
Evenbroekveld 16
9420 Erpe-Mere
Phone: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.com

BRASIL
Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufoli Zanfelice, 329 PB02
Galpão - Lapa
Sao Paulo - SP - Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA
Arjo Canada Inc.
90 Matheson Boulevard West
Suite 350
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 (905) 238-7880
Free: +1 (800) 665-4831
Fax: +1 (905) 238-7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA
Arjo Czech Republic s.r.o.
Na Strži 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
E-mail: info.cz@arjo.com

DANMARK
Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail: dk_kundeservice@arjo.com

DEUTSCHLAND
Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA
ARJO IBERIA S.L.
Polígono Can Salvatella
c/ Cabanyes 1-7
08210 Barberà del Valles
Barcelona - Spain
Telefono 1: +34 900 921 850
Telefono 2: +34 931 315 999

FRANCE
Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG
Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung,
N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA
Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST
Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building,North Tower,
Dubai Science Park,
Al Barsha South
P.O. Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND
Arjo Nederland BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND
Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@arjo.com

NORGE
Arjo Norway AS
Olaf Helssets vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

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At Arjo, we believe that empowering movement within healthcare environments is essential to quality care. Our products and solutions are designed to promote a safe and dignified experience through patient handling, medical beds, personal hygiene, disinfection, diagnostics, and the prevention of pressure injuries and venous thromboembolism. With over 6500 people worldwide and 65 years caring for patients and healthcare professionals, we are committed to driving healthier outcomes for people facing mobility challenges.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

CE

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Arjo Inc.