INSTRUCTIONS FOR USE

TheraRest

Mattress Replacement System

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Bedienungsanleitung · Gebruiksaanwijzing · Istruzioni per l'uso ·
Mode d'emploi · Instrucciones de uso

**WARNING**

To avoid injury, always read this Instructions for use and accompanied documents before using the product.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2022.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

Contents

Introduction	4
Intended use.....	4
Indications	4
Contraindications.....	4
Risks and Precautions.....	4
Safety Information	5
Serious Incident	7
End of life disposal	8
Preparation for Use	9
Mattress Installation.....	10
Inspection	10
Patient placement and nursing care	11
CPR.....	11
Skin Care.....	11
Incontinence / Drainage.....	12
General Operation.....	12
Care and Cleaning	13
Fire Barrier.....	13
Cover Cleaning.....	13
Preventive Maintenance Schedule	15
Daily Cleaning	15
Inspection / System Check-Out.....	15
Troubleshooting	15
Technical specifications	16
Labels	17
Replacement parts	18

Introduction

Intended use

The TheraRest™ Mattress Replacement System (MRS) family of products provide a range of pressure reduction foam mattresses.

Perimeter Plus™ - An advanced technology foam mattress that incorporates a firm side perimeter and a sloped heel section.

Perimeter Plus™ VE - Foam mattress incorporates visco-elastic foam as a top comfort layer. By combining the standard features of the Perimeter Plus mattress with a visco-elastic topper the Perimeter Plus VE mattress is a premier choice for high-performance, pressure-reducing mattresses.

Indications

TheraRest MRS family of products is indicated for the early intervention and prevention of pressure ulcers on low-risk patients.

or

For patients seeking improved comfort.

Contraindications

None.

Risks and Precautions

Transfer – Precaution should be taken during patient transfer, including the locking of caster wheels.

Side Rails and Patient Restraints

WARNING

Use or non-use of restraints, including Side Rails, can be critical to patient safety. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of Side Rails or other restraints.

See related Safety Information.

Patient Migration – Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and /or migration into hazardous positions of entrapment and / or inadvertent bed exit.

NOTE

Monitor patients frequently to guard against patient entrapment.

Skeletal Traction or Unstable Fracture (if not contraindicated) – With skeletal traction, unstable pelvic fracture, or any other unstable fracture (to the extent not contraindicated), maintain physician directed angle of articulation and guard against risks of patient migration.

Safety Information

Patient Entrance / Exit – Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

Bed Frame - Always use a standard healthcare bed frame with this TheraRest mattress, with any safeguards or protocols that may be appropriate. Bed frame and side rails (if used) must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body. It is recommended that bed and side rails (if used) comply with all applicable regulations and protocols.

Brakes – Caster brakes should always be locked once the bed is in position. Verify that the wheels are locked before any patient transfer to or from the bed.

Side Rails / Patient Restraints – Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from the bed) in conjunction with individual patient needs, and should discuss use or non-use with patient

and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, vulnerable patient profile and guidance to further reduce entrapment risks, refer to FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. Monitor patients frequently to guard against patient entrapment.

WARNING

When selecting a standard mattress, ensure the distance between the top of side rails (if used) and top of mattress (without compression) is at least 8.66 in (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

Bed Height – To minimize the risk of falls or injury the bed should always be in the lowest possible position when the patient is unattended.

Head of Bed Elevation – Keep head of bed as low as possible to help prevent patient migration.

Skin Care – Monitor skin conditions regularly, especially at bony prominences and areas where moisture or incontinence may occur or collect and consider adjunct or alternative therapies for high acuity patients. Early intervention may be essential to preventing serious skin breakdown.

Avoid Fire Hazards – To minimize risk of fire, review and follow FDA’s Safety Tips for Preventing Hospital Bed Fires.

No Smoking in Bed – Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

General Protocols – Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

CPR - Level the bed. Lower side rails and initiate CPR per facility protocols. Consider use of a backboard if indicated. After CPR remove the backboard, if used, raise siderails and reconfigure bed and accessories as in initial placement.

Fire Barrier Contact - Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. If the fire barrier becomes damaged or soiled it should be replaced.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

NOTE

For safety information concerning the bed frame refer to the manufacturers’ User Guide.

End of life disposal

Fabric material used on the mattresses or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.

Mattresses at the end of life should be disposed of as waste according to the national or local requirements which may be landfill or combustion.

Preparation for Use

1. Open shipping container.

CAUTION

Do not use sharp instruments to open boxes. Damage to mattress could result.

NOTE

For the TheraRest MRS Perimeter Plus VE models; if stored below room temperature make sure mattress is at room temperature before proceeding to the next step.

NOTE

The mattress cover may appear wrinkled when unpacked. To remove wrinkles, allow mattress up to 24 hours to accommodate; see Troubleshooting for more information. Wrinkles will not affect inflation or function, mattress may be used immediately if needed.

2. Remove the mattress from the plastic protective cover.
3. Check the mattress surface for tears or cracking; do not use if tears or cracks are present.
4. If re-installing mattress onto a new frame or for a new patient, check mattress surface for staining and soiling; clean and / or disinfect as required (see Care and Cleaning).
5. Level bed and lock brakes.
6. Remove the existing mattress from the bed frame, if applicable.

Mattress Installation

1. Position the mattress on the bed frame with the logo facing up and product information tags at the foot end of the bed.
2. Ensure there are no gaps between the mattress and the bed frame or the side rails.

CAUTION

Always use a standard healthcare bed frame with safeguards or protocols that may be appropriate. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.

Inspection

- Check the product surface for staining and soiling, clean and / or disinfect as required.
- Check the product surface for tears or cracking, do not use if significant tears or cracks are present.
- If the mattress cover was recently replaced, ensure that the foam is installed correctly in the new cover, with the top of the foam (i.e. contoured surface or 6.4 cm (2.5 inches) topper foam) towards the teal blue top cover and the sloped foam surface toward the foot end.

Patient placement and nursing care

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications**, **Safety Information** and **Risks and Precautions sections** prior to placing a patient on any **TheraRest MRS**.

1. Transfer the patient following all applicable safety rules and institution protocols.
2. Center patient side-to-side and head-to-foot on TheraRest surface.
3. Ensure that all sections of the mattress fully support the patient.

CAUTION

Mattress handles are for mattress transport only.

CPR

1. Level the bed.
2. Lower or remove the side rail on the caregiver's side if necessary.
3. Begin CPR as per facility protocol. Consider use of the backboard if indicated.
4. After CPR is performed:
 - Remove the backboard if used.
 - Raise or install the side rail as necessary.
 - Reconfigure the bed and accessories as in initial placement.

Skin Care

- Remove excess moisture and keep skin dry and clean.
- Check the patient's skin regularly, particularly in areas where incontinence and drainage occur.
- Ensure linens under the patient are not wrinkled.

Incontinence / Drainage

- Use moisture-impermeable underpads for incontinent patients.
- Wipe surface clean and replace the bed linens as required (see “Care and Cleaning” on page 13).

General Operation

Avoid contact of sharp instruments with TheraRest MRS.

Care and Cleaning

The following process is recommended, but should be adapted to comply with local institution protocols. If you are uncertain, you should seek advice from your local Infection Control Specialist. TheraRest MRS should be routinely decontaminated between patients and at regular intervals while in use.

CAUTION

Do not use Phenol-based solutions or abrasive compounds or pads on the cover during decontamination process as these will damage the surface coating. Do not boil or autoclave the cover. Covers should be wiped down immediately after it has been exposed to fluids or spills.

Fire Barrier

Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over handling of the fire barrier and handle with care. The fire barrier should be replaced if soiled or damaged.

Cover Cleaning

Sewn TheraRest MRS (Top cover cannot be removed from base)

1. Remove or push the bed linens to the center of the mattress.

CAUTION

The sewn TheraRest MRS cover is cleaned by wipe down methods only. Do not launder as damage to the base may occur.

2. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.

Continue with the steps on next page

3. Rinse with a wet towel to remove chemicals.
4. After wipe down, dry the surface with a towel.
5. Ensure the bed linens are refitted and not wrinkled under the patient.
6. Clean pump and tubes (if needed) by wiping with a damp cloth.

Preventive Maintenance Schedule

Preventive maintenance for TheraRest consists of regular cleaning (see Care and Cleaning) and an overall system check-out to be performed at the intervals described below.

All components must be cleaned, disinfected and inspected after each patient's use and before use by a new patient. Always use standard precautions, treating all used equipment as contaminated. Institutions should follow local protocols for cleaning and disinfection.

Daily Cleaning

The cover should be wiped daily with a mild soap and water solution. **After wipe down, dry the surface with a towel.**

Inspection / System Check-Out

Check each of the following before placing TheraRest MRS with a new patient:

1. Check the mattress surface for tears or cracking; do not use if tears or cracks are present.
2. Ensure the mattress is free of stains and is not overly faded.

Troubleshooting

Do not attempt troubleshooting outside this guide or where the solution recommends to contact Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

Technical specifications

WEIGHT	
Maximum recommended patient weight*	158,8 kg (350lb)
TheraRest Perimeter Plus	15 kg (33lb) based on 35 in x 80 in x 7 in MRS
TheraRest Perimeter Plus VE	15.9 kg (35lb) based on 35 in x 80 in x 7 in MRS
*Patient weight capacity may vary depending on frame use. Please consult frame manufacturer.	

THERAREST MATTRESS	CATALOGUE NUMBER	DIMENSIONS (mm/in)		
		Width	Length	Height
TheraRest MRS Perimeter Plus™	TRMRLNB35080TEG	880/35	2020/80	180/7
TheraRest MRS Perimeter Plus™	TRMRLNB35084TEG	880/35	2140/84	180/7
TheraRest MRS Perimeter Plus™ VE	TVMRLNB35080TEG	880/35	2020/80	180/7
TheraRest MRS Perimeter Plus™	TRMRLNB34078TEG	860/34	1980/78	180/7
FOR SWISS MARKET ONLY				
TheraRest MRS Perimeter Plus™	TRMRLNB35080CHG	880/35	2020/80	180/7
TheraRest MRS Perimeter Plus™	TRMRLNB35084CHG	880/35	2140/84	180/7
TheraRest MRS Perimeter Plus™ VE	TVMRLNB35080CHG	880/35	2020/80	180/7
TheraRest MRS Perimeter Plus™	TRMRLNB34078CHG	860/34	1980/78	180/7

Labels

SYMBOL EXPLANATION			
	Do not wash		Do not tumble dry
	Do not iron		Wipeable surface
	Maximum Patient Weight		Do not use Phenol-based cleaning solutions
	Operating instructions - Consult Instructions for use		Foot end
	CE marking indicating conformity with European Community harmonised legislation		Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745
	Unique device identifier		

UK SYMBOL EXPLANATION	
This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.	
	UK marking indicating conformity with UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended)
<p>UK Responsible Person & UK Importer:</p> <p>Arjo (UK) Ltd, ArjoHuntleigh House, Houghton Regis. LU5 5XF</p> <p>Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended).</p> <p>For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.</p>	

Replacement parts

Replaceable TheraRest components are listed below. For more information such as pricing or additional spare parts that are not on this list, please contact your local Arjo representative.

TheraRest mattress components

FOAM SHELLS	CATALOGUE NUMBER
TheraRest Perimeter Plus foam shell	TRFXXX35080XXX TRFXXX35084XXX
TheraRest Perimeter Plus VE foam shell	TVFXXX35080XXX
TheraRest Perimeter Plus foam shell	TRFXXX34078XXX

FIRE BARRIER LINING	CATALOGUE NUMBER
Fire barrier lining	633180

THERAREST MATTRESS COVERS	CATALOGUE NUMBER
Complete Cover, Reliant Top, Non-Skid Black Bottom, with handles	TRCRLNB35080TXG
Complete Cover, Reliant Top, Non-Skid Black Bottom, with handles	TRCRLNB35084TXG
Complete Cover, Reliant Top, Non-Skid Black Bottom, with handles	TVCRLNB35080TXG
Complete Cover, Reliant Top, Non-Skid Black Bottom, with handles	TRCRLNB34078TXG

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Intentionally left blank

AUSTRALIA

Arjo Australia
 Building B, Level 3
 11 Talavera Road
 Macquarie Park, NSW, 2113, Australia
 Phone: 1800 072 040

BELGIQUE / BELGIË

Arjo Belgium
 Evenbroekveld 16
 9420 Erpe-Mere
 Phone: +32 (0) 53 60 73 80
 Fax: +32 (0) 53 60 73 81
 E-mail: info.belgium@arjo.com

BRASIL

Arjo Brasil Equipamentos Médicos Ltda
 Rua Marina Ciufulli Zanfelice, 329 PB02
 Galpão - Lapa
 São Paulo - SP - Brasil
 CEP: 05040-000
 Phone: 55-11-3588-5088
 E-mail: vendas.latam@arjo.com
 E-mail: servicios.latam@arjo.com

CANADA

Arjo Canada Inc.
 90 Matheson Boulevard West
 Suite 350
 CA-MISSISSAUGA, ON, L5R 3R3
 Tel/Tél: +1 (905) 238-7880
 Free: +1 (800) 665-4831
 Fax: +1 (905) 238-7881
 E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA

Arjo Czech Republic s.r.o.
 Na Strizi 1702/65
 140 00 Praha
 Czech Republic
 Phone No: +420225092307
 E-mail: info.cz@arjo.com

DANMARK

Arjo A/S
 Vassingerødvej 52
 DK-3540 LYNGBE
 Tel: +45 49 13 84 86
 Fax: +45 49 13 84 87
 E-mail: dk_kundeservice@arjo.com

DEUTSCHLAND

Arjo GmbH
 Peter-Sander-Strasse 10
 DE-55252 MAINZ-KASTEL
 Tel: +49 (0) 6134 186 0
 Fax: +49 (0) 6134 186 160
 E-mail: info-de@arjo.com

ESPAÑA

ARJO IBERIA S.L.
 Polígono Can Salvatella
 c/ Cabanyes 1-7
 08210 Barberà del Valles
 Barcelona - Spain
 Telefono 1: +34 900 921 850
 Telefono 2: +34 931 315 999

FRANCE

Arjo SAS
 2 Avenue Alcide de Gasperi
 CS 70133
 FR-59436 RONCQ CEDEX
 Tél: +33 (0) 3 20 28 13 13
 Fax: +33 (0) 3 20 28 13 14
 E-mail: info.france@arjo.com

HONG KONG

Arjo Hong Kong Limited
 Room 411-414, 4/F, Manhattan Centre,
 8 Kwai Cheong Road, Kwai Chung,
 N.T.,

HONG KONG

Tel: +852 2960 7600
 Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
 Via Giacomo Peroni 400-402
 IT-00131 ROMA
 Tel: +39 (0) 6 87426211
 Fax: +39 (0) 6 87426222
 E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
 Office 908, 9th Floor,
 HQ Building, North Tower,
 Dubai Science Park,
 Al Barsha South
 P.O. Box 11488, Dubai,
 United Arab Emirates
 Direct +971 487 48053
 Fax +971 487 48072
 Email: Info.ME@arjo.com

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At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

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Arjo Inc.