11004 Bariatric Loop Repositioning Sling

EN Instruction for use
FR Instructions d'utilisation
ES Instrucciones de uso
DE Gebrauchsanweisung
IT Istruzioni per l'uso
NL Gebruiksaanwijzing
CN 使用說明書

Loop strap
Sangle avec boucles
Correas de bucles
Schlaufenbänder
Cinghie con anelli
Bevestigingslussen
錶帶

Body
Pièce de corps
Pieza del cuerpo
Körperstück
Piece corpo
Stoffering
本体

Product label
L'étiquette du produit
Etiqueta del producto
Produktetikett
Etichetta del prodotto
Etiket van het product
產品標籤

Warning! - Avertissement! - Advertencia! - Attenzione! - Warnung! - Waarschuwing - 警告! 請

Always read this instruction and the instruction for the lifter to be used with the sling
Toujours lire cette instruction et l'instruction pour l'élévateur à être utilisés avec la fronde
Siempre lea esta Instrucción y la Instrucción para el elevador que se utilizar n con la Honda Lesen
Sie immer diese Anleitung und die Anleitung für die Lifter, die mit der dem Gurt verwendet werden
Leggere sempre questa Istruzione e l'Istruzione per il sollevatore da utilizzare con la fionda
Lees altijd deze instructie en de instructie van de tillift waar je deze tilband mee gebruikt
务必阅读本使用说明书以及使用吊兜的升降机的说明书。

Technical Specifications

Part no.	Description	Size	SWL	Material	Dimension
11004	Bariatric loop repositioning sling	Onesize	Max. 454 kg / 1000 lbs	Polyester (PES) Cotton (CO)	145 x 210 cm / 57 x 82 in

Definition in this IFU

Warning: Means safety warning. Failure to understand or obey this warning may result in injury to you or to others.

Note means: Important information for the correct use of this system or equipment.

Intended Use

The Sling is intended for lateral transfer of patients between adjacent support surfaces, moving patient up in the bed and turning.

Note: The sling must not be used for moving and lifting in sitting position.

Safety instructions

Always read the label and check that the sling is in good condition without signs of tear and wear such as loose seams, frays and material inconsistency.

Always make sure not to exceed the Service Work Load (Safe Working Load) defined for the product.

Always read the IFU for the lifting equipment to perform the lifting.

Warning! Do not use the product if any signs of deterioration. If sling is dirty, it should be cleaned according to the washing instructions on label.

Before and after every use

1. Always perform a risk assessment before starting the specific lifting session in order to make sure the product model is correct and that the lifting procedure is fully understood.
2. Before using the sling, read this IFU and the IFU for the lifting equipment to be used.
3. When soiled or stained the sling must be replaced immediately.
4. Check the limb sling for any fraying, loose stitching, attachments or any deteriorations.
5. Warning: Remember to mark patient name on the care label in order to avoid cross infection between patients.
6. After use, check sling according to the step 3 and 4 above.
7. Store in a dry area.

Attachments

The loop slings will fit most lifts with standard hangerbars with loop attachment hooks. Check the lift manufacturers documentation for details.

Service & maintenance

The sling shall be checked by authorized personnel every 6 months according to ISO10535

Customer support/trouble shooting

Please refer to the manual for the lifting equipment used
In doubt please contact your dealer who supplied the product

Applying the sling: Moving up in bed or lateral transfer

1. Position the repositioning sling centrally on the bed with the label pointing at the head end of the bed. The label side faced down towards the mattress.
2. Position the hangerbar carefully over the patient (across the patients chest); connect the straps to each side of the hanger bar, starting with the head end straps, then the patients' hips straps and leg straps. For lateral transfer and moving up in bed, use a minimum of five straps on each side.
3. Check patient position and carefully make the required adjustment to achieve the best comfort for the patient, especially, that the head is supported and legs and knee's position is not causing any unattended force on knee joints before lifting.
4. Raise the hanger bar until all straps are tense to make sure the straps are securely attached.
5. Carefully perform the lift, repositioning or lateral transfer, whilst keeping eye on the patient position and comfort.
6. At desired position, slowly lower the hangerbar until all straps can be released from the hangerbar hooks.

Applying the sling: Turning in bed

1. Position the repositioning sling centrally on the bed with the label pointing at the head end of the bed. The label side faced down towards the mattress.
2. Position the hanger bar carefully in parallel position over the patient.
3. Attach a minimum of four straps on one side of the sling to the hangerbar. Use the straps at the patients shoulder and hip area. Use a minimum of two hooks on the hangerbar.
4. Check patient position and carefully make the required adjustment to achieve the best comfort for the patient.
5. Check the support of the head and that arms and legs are in the right position.
6. Raise the hanger bar until all straps are tense to make sure all straps are securely attached.
7. Carefully perform the turning, whilst keeping eye on the patient position and comfort.
8. At desired position, slowly lower the hanger bar and release the straps from the hangerbar hooks.

Expected service life & end of life disposal

- Expected service life varies between 1–5 years, depending on usage and number of washing cycles. When washed once a week the sling has an expected service life of 2 years. Expected service life is based upon a correct usage and washing procedure in accordance with this IFU. The sling is washable according to the marking on the label.
- Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.
- The packaging bag consist of soft plastic and should be recycled according to local legislations.

Storing and shelf life

The slings should be stored in a dry and clean environment and the shelf life is 5 years in original packaging. Store in temperatures between 0 - 40°C (32 - 104°F).

Safe Working Load (SWL)

Max. 454 kg / 1000 lbs.

Product care & maintenance

Always follow the instructions on the product label and check the sling for deteriorations. Mark the sling with the specific patient name and keep in dry and clean environment.

Washing Instruction

Machine wash the sling at 70°C (176°F).
Tumble dry at low heat.
Close the buckle prior to washing.
Do Not heat press slings equipped with buckles.
Do Not Bleach.
Do Not clean dry.
Do Not Iron.

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