
# 22003 Disposable Loop Turning Sling

Patient Specific Disposable  
 Jetable, patient spécifique  
 Paciente desechable específico  
 Patientenspezifische Einweg  
 Specifico paziente monouso  
 Patiënt Specifieke Disposable  
 患者一次性专用产品

EN Instruction for use  
 FR Instructions d'utilisation  
 ES Instrucciones de uso  
 DE Gebrauchsanweisung  
 IT Istruzioni per l'uso  
 NL Gebruiksaanwijzing  
 CN 使用说明书

## Warning! - Avertissement! - Advertencia! - Attenzione! - Warnung! - Waarschuwing - 警告! 請

Always read this instruction and the instruction for the lifter to be used with the sling  
 Toujours lire cette instruction et l'instruction pour l'élévateur à être utilisés avec la fronde  
 Siempre lea esta Instrucción y la Instrucción para el elevador que se utilizarán con la Honda Lesen  
 Sie immer diese Anleitung und die Anleitung für die Lifter, die mit dem Gurt verwendet werden  
 Leggere sempre questa Istruzione e l'Istruzione per il sollevatore da utilizzare con la fionda  
 Lees altijd deze instructie en de instructie van de tillift waar je deze tilband mee gebruikt  
 务必阅读本使用说明书以及使用吊兜的升降机的说明书。

## Technical Specifications

Part no.	Description	Size	SWL	Material	Dimension
22003	Disposable loop turning sling	Onesize	Max. 100 kg / 220 lbs	Polypropylene (PP) Polyethylene (PE)	80 x 45 x 14 cm / 31.5 x 17.7 x 5.5 in

## Definition in this IFU

---

Warning: Means Safety warning. Failure to understand or obey this warning may result in injury to you or to others.

Note: Means important information for the correct use of this system or equipment.

## Intended Use

---

The sling will assist the caregiver in turning/rolling a patient in a supine position to perform various tasks such as, inspection of skin, treatment, and personal hygiene.

Note! The sling is designed for turning and not lifting.

## Safety instructions

---

- Always read the label and check that the sling is in good condition without signs of tear and wear such as loose seams, frays and material inconsistency.
- Always make sure not to exceed the SWL (Safe working Load) defined for the product.
- Always read the IFU for the lifting equipment to perform the turning.

**Warning! Do not use the product if there is any signs of deterioration. Do Not Wash - this product is designed for single patient use only and should be disposed of if in contact with liquid of any kind.**

## Before and after every use

---

1. Always perform a risk assessment before starting the specific turning session in order to make sure the product model is correct and that the turning procedure is fully understood.
2. Before using the sling, read this IFU and the IFU for the lifting equipment to be used.
3. When soiled or stained the sling must be replaced immediately.
4. Check the sling for any fraying, loose stitching, attachments or any deteriorations.
5. Warning: Remember to mark patient name on the care label in order to avoid cross infection between patients.
6. After use, check sling according to the step 3 and 4 above
7. Store in a dry area.

## Attachments

---

The sling will fit most lifts with standard hangerbars with loop attachment hooks. Check the lift manufacturers documentation for details.

## Applying the sling

---

1. Push the velcro strap by hand between the mattress and the patient's hollow at the back.
2. Pull the sling carefully from the other side and attach the strap to the bed frame to which direction you want to turn the patient. Position the sling and make sure the sling is not folded or wrinkled. If needed, adjust the length of velcro strap.
3. Check patient position, especially arm and legs before mounting to the lifter.
4. Place the loop straps in the hanger bar according to the manual for the lifter.
5. Carefully raise the hanger bar and perform the turning, while keeping eye on the patient position.

## Expected service life & end of life disposal

---

- The service life is expected to be 56 performed lifting sessions, however this will be subject to the actual conditions during use. The sling are disposable and designed for single patient use and should be disposed of after patient hospital stay or when soiled or subjected to liquid.
- Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.
- The packaging bag consist of soft plastic and should be recycled according to local legislations.

## Storing and shelf life

---

The slings should be stored in a dry and clean environment and the shelf life is 5 years in original packaging. Store in temperatures between 0 - 40°C (32 - 104°F).

## Safe working load (SWL)

---

Max. 100 kg / 220 lbs.

## Product care & maintenance

---

- Always follow the instructions on the product label and check the sling for deteriorations.
- Mark the sling with the specific patient name and keep in dry and clean environment.
- This is a patient specific disposable sling, that must be disposed of after the patient has been discharged.

## Service & maintenance

---

The sling shall be checked by authorized personnel according to ISO10535.

## Customer support/trouble shooting

---

- Please refer to the manual for the lifting equipment used
- In doubt please contact your dealer who supplied the product.

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