
Flowtron DVT30

arjo

Thigh Garment

To be used only under the direction of a physician • Non-Sterile • Not made with natural rubber latex • For Single Patient Use Only

 For use only with Flowtron® DVT-prevention pumps manufactured by Arjo. Must not be used with Flowtron® Hydroven 3 or Flowtron® Hydroven 12 Intermittent Pneumatic Compression (IPC) pumps.

 More comprehensive information on the garment can be found in the relevant Flowtron DVT pump Instructions For Use document.

Instructions for Use

1. Plug the pump into a suitable electrical outlet. Do not turn the pump on at this time.
2. Remove the garments from the sealed bag. The garments may be used on either leg. Unfold the garment and position the inflatable bladder directly behind the patient's calf.
3. Snugly wrap the garment around the patient's leg and secure the fastener tabs. Repeat for the other leg.
4. Attach the garments to the pump tubing set ensuring a 'click' is heard from each snap-lock connector.
5. If your Flowtron DVT pump has an adjustable pressure regulator, turn the regulator dial to the recommended pressure unless otherwise directed by the physician. For further details, refer to the relevant Flowtron DVT pump Instructions For Use document.
6. Turn the pump on. The green power indicator lights should illuminate.

Intended Use

Inflate the thigh garment by Intermittent Pneumatic Compression system. The garment will squeeze the thigh muscles to facilitate the blood and Lymph circumference.

Indications

To help prevent Deep Vein Thrombosis (DVT).

Recommendations

— General recommendations

- Garments should be removed regularly to inspect the skin for signs of redness or pressure points.
- Where appropriate, patients should be instructed in the proper use of the system, the purpose of the therapy and that any problems should be reported to the nursing staff.

DVT prophylaxis:

- The garments should be applied to the patient pre-operatively, prior to the induction of anaesthesia.
- The system should be used continuously for no less than 72 hours post-operatively or until the patient becomes fully ambulatory.
- If the garment cannot be applied to the operative limb during surgery, it may be applied to the limb once the patient reaches the recovery unit.
- In the non-surgical patient, the system should be initiated immediately once the risk of DVT formation is identified.

Contraindications

IPC should not be used in the following conditions:

- Severe arteriosclerosis or other ischaemic vascular diseases.
- Known or suspected acute Deep Vein Thrombosis (DVT) or phlebitis.

- Severe congestive cardiac failure or any condition where an increase of fluid to the heart may be detrimental.
- Pulmonary embolism.
- Any local condition in which garments would interfere, including gangrene, recent skin graft, dermatitis or untreated, infected leg wounds.

If you are unsure of any contraindications refer to the patient's physician before using the device.

Cautions

- Garments should be removed immediately if the patient experiences tingling, numbness or pain.
- When used for DVT prophylaxis, continuous use is recommended and any interruption of therapy for a substantial length of time should be at the discretion of the physician.

Manufacturing date: see label.

Expired date: 3 years

Symbols

Single Patient Use

CE marking indicating conformity with European Community harmonized legislation

Indicates the product is a Medical Device according to EU Medical Device regulation 2017/745

Operating instructions - Consult Instructions for use

Unique device identifier

Design Policy and copyright

® and ™ are trademarks belonging to the Arjo group of companies.


As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

The content of this publication may not be copied either whole or in part without the consent of Arjo.

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