


INSTRUCTIONS FOR USE

Skin IQ 1000

Microclimate Manager



Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2023.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

Intentionally left blank

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

ARJO HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJO PRODUCT AS DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY ARJO SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL ARJO BE LIABLE FOR ANY INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJO TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in Arjo printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties, except as set forth in the written limited warranty included with this product. Information in this publication may be subject to change at any time. Contact Arjo for updates.

IMPORTANT INFORMATION FOR USERS

In order for Arjo products to perform properly, Arjo recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- **WARNING:** Assembly, operations, adjustments, extensions, modifications, technical maintenance or repairs must be performed only by qualified personnel authorized by Arjo. Contact Arjo for information regarding maintenance and repair.
- Ensure the electrical installation of the room complies with the appropriate national electrical wiring standards. To avoid the risk of electric shock, this product must be connected to a grounded power receptacle.

Specific indications, contraindications, warnings, precautions and safety information exist for Arjo's therapeutic support systems. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

NOTICE

This product has been configured from the manufacturer to meet specific voltage requirements. Refer to the power supply label for specific voltage.

TABLE OF CONTENTS

QUICK REFERENCE.....2

INTRODUCTION4

 Indications4

 Contraindications.....4

 Intended Care Setting4

 Compatibility5

 Risks and Precautions5

 Safety Information6

UNPACKING AND INSPECTION.....7

INSTALLATION.....7

CARE AND CLEANING8

END OF LIFE DISPOSAL8

SPECIFICATIONS9

ELECTROMAGNETIC COMPATIBILITY10

WARRANTY13

SYMBOLS USED.....14



QUICK REFERENCE

Pull coverlet tightly over corners.


For powered air surfaces that have a CPR release feature, please make sure the CPR is visible.




Smooth out wrinkles.




Connect power supply to fan.



Move hand to end of bed and feel fan vibration, listen for humming sound.





Skin IQ 1000

INTRODUCTION

This document should be saved in an easily accessible place for quick reference.

It is recommended that all sections of these instructions be read prior to product use. Carefully review the **Indications, Contraindications, Risks and Precautions** and **Safety Information** prior to placing a patient on the Skin IQ® 1000 Microclimate Manager (MCM).



These instructions do not provide specific safety or operational information for the pressure redistribution surface and / or bed frame provided by the facility for use with the Skin IQ 1000 MCM. Consult product labeling for information.

Caregivers should discuss **Safety Information, Risks and Precautions** and **Contraindications** with the patient (or the patient's legal guardians) and the patient's family.

The *Skin IQ* 1000 MCM is a single patient use, disposable device that provides Negative Airflow Technology (NAT) to manage the microclimate of the skin at the patient surface when fitted over a customer-provided pressure redistribution surface.

It is also designed to reduce friction and improve patient comfort.

The *Skin IQ* 1000 MCM is suitable for use in acute and post acute facilities, is vapor permeable, and has a fluid-resistant nylon taffeta cover.

Indications

The *Skin IQ* 1000 MCM is indicated for use in conjunction with a pressure redistribution surface in order to aid in the prevention and treatment of skin breakdown and pressure ulcers (Stages I-IV) for patients who require microclimate management of the skin.

Contraindications

Although *Skin IQ* 1000 MCM has no associated direct contraindications the caregiver should refer to and follow any contraindications in the product labeling for the pressure redistribution surface and / or bed frame being used with the *Skin IQ* 1000 MCM.

Intended Care Setting

- Acute Care
- Post Acute Care

Compatibility

The *Skin IQ* 1000 MCM is designed to fit on a pressure redistribution surface that is 203 cm - 213.4 cm (80 - 84 in) long by 121.9 cm (48 in) wide by 17.8 cm (7 in) high.

Consult product labeling for the pressure redistribution surface and / or bed frame for compatibility.

Risks and Precautions

Transfer

This product is not intended for use as a transfer device.

Duration of Use

Recommended duration of use for single patient is not more than 60 days for patients up to 453.59 kg (1000 lb).

However, patient specific duration of use may vary. Clinical conditions such as, but not limited to incontinence, skin condition, nutrition status, medications, mobility, weight, or etiology, need to be considered when assessing duration of use for the *Skin IQ* 1000 MCM.

Height

The *Skin IQ* 1000 MCM will increase the height of the pressure redistribution surface it is applied to by approximately 6.35 mm (0.25 in).

Use With Other Devices

All *Skin IQ* 1000 MCM components are designed to be used as a single system device.

The included power supply should only be used with the *Skin IQ* 1000 MCM. The coverlet should only be powered by the *Skin IQ* 1000 Power Supply Arjo part number 4103832.



Any attempt to connect and use the power supply with any other device, or use any other brand or model of power supply other than Arjo part number 4103832 will result in improper operation of equipment, possibly leading to increased risk of patient injury.

Patient Migration

Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit. Monitor patients frequently to guard against patient entrapment.

Safety Information



To avoid serious injury or death, the CPR rapid deflation unit must be visible and accessible at all times.



Please refer to and follow any safety information in the product labeling for the pressure redistribution surface and / or bed frame being used with the Skin IQ 1000 MCM.

Power Cord

Only use a grounded power outlet and the power cord supplied with the *Skin IQ 1000* power supply. The power cord should be positioned to avoid a tripping hazard and / or damage to the cord. The *Skin IQ 1000* MCM should never be operated with a worn or damaged power cord. Should the power cord become worn or damaged, contact Arjo or an Arjo authorized representative to order a replacement.

Coverlet

Use care when handling or transporting. Dropping or other sudden impacts may result in damage to the device.

Skin Care

Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.



The Skin IQ MCM Product Family has an antimicrobial agent formulated into the patient contact layer. Although rare, there is a potential that some patients may experience sensitivity or a reaction during use. Regularly monitor the patient's skin condition. Discontinue use and seek medical treatment if any signs of a reaction are observed.

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

For more information about the Skin IQ MCM Product Family, please visit www.arjo.com.

Patient Weight

The maximum patient weight for this device is 453.59 kg (1000 lb). In addition, consult the specifications for the pressure redistribution surface and / or bed frame being used. Additional weight limitations may apply.

Disposal

The coverlet contains electrical components that may be removed and disposed of separately. See Disposal on page 8 for removal.

General Protocols

- Avoid contact of sharp instruments with the *Skin IQ 1000* MCM. Punctures, cuts and tears will prevent proper operation.
- Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

UNPACKING AND INSPECTION

Unpack the *Skin IQ* 1000 MCM from the shipping box and locate items as listed.

- coverlet with instructions for use booklet (included in complete kit or coverlet only)
- power supply (included in complete kit or available as a separate item)
- power cord (included in complete kit or power supply only)

Inspect all items carefully. If any items are damaged or missing, contact Arjo or an Arjo authorized representative.

INSTALLATION

See the **Quick Reference** on pages 2 and 3 for illustrations of the following procedures.



Failure to properly secure the coverlet to the existing surface may lead to patient or user injury or equipment damage.

1. Remove coverlet from shipping bag.
2. Remove all covers and sheets from the existing pressure redistribution surface.
3. Place the coverlet on top of the existing surface, ensuring the foot graphic on the coverlet is at the foot end of the bed.
4. Connect the power supply to the fan located underneath the foot end of the coverlet.
5. Ensure the cord that runs from the fan to the power supply is placed on the floor under the bed. Improper placement of the cord could cause injury.
6. Pull the coverlet over the pressure redistribution surface by stretching it over each corner securely. Do not trap the power cord between the coverlet and mattress.
7. Smooth any wrinkles on the coverlet.
8. Secure the coverlet to the bed by using the hook-and-loop straps, located on the underside of the coverlet.
9. Ensure strap placement does not interfere with the operation of the bed functions. Failure to do so could result in patient injury or equipment damage.
10. Write the therapy start date on the law tag at the foot end of the coverlet to track use. Failure to track duration of use may void warranty.
11. Confirm there are no sharp objects in the immediate area which may damage the coverlet.
12. Connect the power supply to the power cord.
13. Connect the power cord to a properly grounded electrical outlet and confirm outlet has power. Verify that the electrical outlet can be easily accessed when disconnecting the device from mains power. Ensure power supply and cord are properly stored on the floor beneath bed.
14. Move hand to the end of the bed and press the foot end of coverlet to feel fan vibration. A low humming sound indicates the fan is working.

Connecting Skin IQ MCM to Arjo Surfaces with Skin IQ Integration

The Skin IQ MCM coverlet can connect directly to Arjo surfaces that have an integrated Skin IQ port located on the pump, removing the need to connect the Skin IQ MCM coverlet to the mains power.

To connect the Skin IQ MCM coverlet to Arjo surfaces with Skin IQ integration, a separate accessory is required: Product Code 636377

Skin IQ Integration Coiled Power Cable

1. Connect one end of the Coiled Skin IQ Power Cable to the fan located underneath the foot end of the coverlet.
2. Connect the other end of the Coiled Skin IQ Power Cable to the Skin IQ port located on the pump.
3. The Skin IQ MCM coverlet will be in full operation at all times when the pump is switched on.
4. To confirm operation of the Skin IQ MCM coverlet move hand to the end of the bed and press the foot-end of the coverlet to feel fan vibration. A low humming sound indicates the fan is working.

CARE AND CLEANING

During patient use, clean *Skin IQ 1000* Coverlet by wiping with a mild soap and water solution. Do not launder the *Skin IQ 1000* Coverlet. The *Skin IQ 1000* Coverlet is for single patient use only.

Avoid spilling fluids on any part of the *Skin IQ 1000* power supply. If spills do occur:

- disconnect the power cord from the wall outlet
- clean fluids from the product




Ensure that there is no moisture in or near the power supply and power plug before reconnecting the power cord.

END OF LIFE DISPOSAL

The coverlet itself is a single patient use only product, but some of the items that come with it can be reused if they are handled properly when removed. Follow the steps below to remove the fan and dispose of all components of the *Skin IQ 1000* MCM.

1. Disconnect power cord from wall outlet.
2. Disconnect power supply from fan.



Power supply and cord can be reused. Consider all facility policies and procedures with regard to cleaning, inspection and reuse of electronic equipment. If not reused, dispose of power supply per approved local institutional procedures.

3. Cut fan out of coverlet as shown.
4. Dispose of coverlet and fan according to approved local institutional procedures.



Fan assembly contains electronic components that may require alternate disposal than the soft goods of the coverlet.



Improper disposal of any component may result in regulatory non-compliance.

Fabric material used on the coverlet or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.

Units have electrical and electronic components that should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

SPECIFICATIONS*

Maximum Weight Capacity.....453.59 kg (1000 lb)

Recommended duration of use for single patient is not more than 60 days for patients up to 453.59 kg (1000 lb).



**Consult the specifications for the pressure redistribution surface being used.
Additional weight limitations may apply.**

Electrical:

Voltage 100 - 240 VAC
Frequency 50 / 60 Hz
Voltage (Saudi Arabia only) 230 VAC
Frequency (Saudi Arabia only)..... 60 Hz
Ampere Rating 0.5 A
Maximum Electrical Leakage 100 uA at 115 VAC 60 Hz and 200 uA at 230 VAC 50 Hz
Power Cord Length 6 m (19.69 ft)

Environmental Conditions:

Operating:

Temperature Range.....14°C (57.2°F) to 35°C (95°F)

Transport / Storage:

Temperature Range.....-29°C (20.2°F) to 60°C (140°F)

The *Skin IQ 1000 Coverlet* is classified as a Type B applied part under IEC 60601-1:2005/A1:2012 (3.1 edition).

***Specifications subject to change without notice.**

ELECTROMAGNETIC COMPATIBILITY

Electromagnetic Interference - Although this equipment conforms with the intent of the directive 2014/30/EU in relation to Electromagnetic Compatibility (EMC), all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer.

Portable and mobile RF communications equipment can effect medical electrical equipment.

Radios, cell phones and similar devices may affect this equipment and should be kept at least 2 m (6.5 ft) away from the equipment.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the following tables.

The following tables document compliance levels and guidance from the IEC 60601-1-2:2014 (4th edition) Standard, for the electromagnetic environment in which the *Skin IQ 1000 MCM* should be used in a clinical environment.

Guidance and manufacturer's declaration - electromagnetic emissions		
The <i>Skin IQ 1000 MCM</i> is intended for use in the electromagnetic environment specified below. The customer or user of the <i>Skin IQ 1000 MCM</i> should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic Environment
RF emissions CISPR 11	Group 1	The <i>Skin IQ 1000 MCM</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. This equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations / flicker emissions IEC 61000-3-3	yes	

Guidance and manufacturer's declaration - electromagnetic immunity

The *Skin IQ* 1000 MCM is intended for use in the electromagnetic environment specified below. The customer or user of the *Skin IQ* 1000 MCM should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV Contact ±15 kV Air	±8 kV Contact ±15 kV Air	In accordance with IEC 60601-1-2, if the floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±1 kV lines to ±2 kV for power supply lines	
Surge IEC 61000-4-5	1 kV line(s) to Line(s) 2 kV line(s) to Earth	1 kV line(s) to Line(s) 2 kV line(s) to Earth	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5% half cycle 40% for 5 cycle 70% for 25 cycle 5% for 5 seconds	5% half cycle 40% 5 cycles 70% 25 cycles 5% for 5 seconds	
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/M	30 A/M	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_i is the a.c. mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the *Skin IQ* 1000 MCM


The *Skin IQ* 1000 MCM is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the *Skin IQ* 1000 MCM can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *Skin IQ* 1000 MCM as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter meters		
	150 kHz to 80 MHz not applicable	80 MHz to 800 MHz $d = \left\lceil \frac{3.5}{E_1} \right\rceil \sqrt{P}$	800 MHz to 2.5 GHz $d = \left\lceil \frac{7}{E_1} \right\rceil \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.37	0.74
1	N/A	1.2	2.3
10	N/A	3.7	7.4
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separate distance *d* in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1, At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from surfaces, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity			
The <i>Skin IQ</i> 1000 MCM is intended for use in an electromagnetic environment specified below. The customer or user of the <i>Skin IQ</i> 1000 MCM should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3Vrms 150K - 80 MHz	3 Vrms 150K - 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the <i>Skin IQ</i> 1000 MCM, including cables, than the recommended separation distance calculated from the equation application to the frequency of the transmitter. Recommended Separation Distance Battery Operated Device $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (see note a) should be less than the compliance level in each frequency range (see note b). Interference may occur in the vicinity of equipment marked with the following symbol: </p>
NOTE 1, At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the <i>Skin IQ</i> 1000 MCM is used exceeds the applicable RF compliance level above, the <i>Skin IQ</i> 1000 MCM should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the <i>Skin IQ</i> 1000 MCM.</p> <p>b) Over the frequency range 150kHz, field strengths should be less than 3 V/m.</p>			

WARRANTY

In connection with your purchase of a *Skin IQ 1000 MCM*, Arjo and its affiliates, (collectively referred to herein as "Arjo") warrants, to the original purchaser, its *Skin IQ 1000 MCM* against manufacturer's defects in material and / or workmanship for a period of sixty (60) days from first date of product placement on patient surface; or one year from date of product purchase, whichever comes first under normal usage and so long as all applicable procedures are followed per the instructions for use (IFU).

This Limited Warranty is non-transferable. THIS LIMITED WARRANTY SHALL BE IN LIEU OF ANY AND ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ALL OF WHICH ARE EXPRESSLY DISCLAIMED BY ARJO REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS WARRANTY IS THE CUSTOMER'S SOLE REMEDY. ARJO SHALL NOT BE LIABLE FOR ANY INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, FOR THE BREACH OF ANY WARRANTY ON THIS PRODUCT. SOME STATES AND JURISDICTIONS MAY NOT ALLOW THESE LIMITATIONS ON WARRANTIES. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE OR JURISDICTION TO JURISDICTION.

Conditions and Limitations:

Under this warranty, upon reasonable notice, Arjo will replace defective parts or whole units covered under this warranty at its sole option. No returns or replacements will be allowed without authorization from Arjo. This warranty applies to material/parts replacement costs and labor costs only and does not include shipping and handling costs or disposal fees.

This warranty does not extend to, nor cover:

- Patient wearables, such as, but not limited to, non-integrated cover sheets(s) and disposables used in conjunction with the *Skin IQ 1000 MCM* product; or
- Normal wear and tear; or
- Damage, or product failure due to causes beyond Arjo's control such as, but not limited to, abuse, theft, fire, flood, wind, lightning, freezing, power failure, power reduction, clogging of coverlet pores due to tobacco smoke, unusual atmospheric conditions, or force majeure.

This warranty is VOID if:

- Proof of purchase cannot be made by the original purchaser;
- Adjustment, modification, and / or repair beyond that described in the *Skin IQ 1000 MCM* IFU is carried out by entities not authorized by Arjo;
- The electrical installation of the room does not comply with the appropriate national electrical wiring standards;
- The product is not installed or used in accordance with the instructions and warnings enumerated in the *Skin IQ 1000 MCM* IFU

This warranty is in no way to be construed as an extension of any other agreement entered into by the parties.

Arjo reserves the right to make material changes to the product predicated upon the availability of raw materials of like quality.

SYMBOLS USED



Conforms to AAMI ES60601-1-6.
Certified to CSA Std. C22.2 No. 60601-1.



No Hooks



Important Operational Information



Foot End



CE marking indicating conformity with European Community harmonised legislation. Figures indicate Notified Body supervision.



Direct Current



Consult Instructions for Use



Single Patient Use



Manufacturer



Temperature Low and High Limits



Protected against ingress of liquids



Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745



Hospital name



Maximum patient weight



Do Not Open With Scissors



Content Information



Tripping Hazard



This product or its parts are designated for separate collection at an appropriate collection point. At the end of useful life, dispose of all waste according to local requirements, or contact your local Arjo representative for advice.



Alternating Current



Keep Dry



Warning of possible hazard to system, patient or staff



Date of Manufacture



Do Not Shower



Type B Applied Part



Protective Earth Class 1 Device



Wipe surfaces with cleaning solution, then wipe with a cloth moistened with water and dry thoroughly



Date of first use

AUSTRALIA

Arjo Australia
Building B, Level 3
11 Talavera Road
Macquarie Park, NSW, 2113,
Australia
Phone: 1800 072 040

BELGIQUE / BELGIË

Arjo Belgium
Evenbroekveld 16
9420 Erpe-Mere
Phone: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.com

BRASIL

Arjo Brasil Equipamentos Médicos
Ltda
Rua Marina Ciufuli Zanfelice, 329
PB02 Galpão - Lapa
São Paulo - SP - Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA

Arjo Canada Inc.
90 Matheson Boulevard West
Suite 350
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 (905) 238-7880
Free: +1 (800) 665-4831
Fax: +1 (905) 238-7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA

Arjo Czech Republic s.r.o.
Na Strži 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
E-mail: info.cz@arjo.com

DANMARK

Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail: dk_kundeservice@arjo.com

DEUTSCHLAND

Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA

ARJO IBERIA S.L.
Poligono Can Salvatella
c/ Cabanyes 1-7
08210 Barberà del Valles
Barcelona - Spain
Telefono 1: +34 900 921 850
Telefono 2: +34 931 315 999

FRANCE

Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG

Arjo Hong Kong Limited
Room 411-414, 4/F,
Manhattan Centre,
8 Kwai Cheong Road,
Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building, North Tower,
Dubai Science Park,
Al Barsha South
P.O. Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND

Niederland BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND

Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@Arjo.com

NORGE

Arjo Norway AS
Olaf Helseth vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

ÖSTERREICH

Arjo Austria GmbH
Lemböckgasse 49 / Stiege A / 4.OG
A-1230 Wien
Tel: +43 1 8 66 56
Fax: +43 1 866 56 7000

POLSKA

Arjo Polska Sp. z o.o.
ul. Ks Piotra Wawrzyniaka 2
PL-62-052 KOMORNIKI (Poznań)
Tel: +48 691 119 999
E-mail: arjo@arjo.com

PORTUGAL

Arjo em Portugal
MAQUET Portugal, Lda.
(Distribuidor Exclusivo)
Rua Poeta Bocage n.º 2 - 2G
PT-1600-233 Lisboa
Tel: +351 214 189 815
Fax: +351 214 177 413
E-mail: Portugal@arjo.com

SUISSE / SCHWEIZ

Arjo Switzerland AG
Fabrikstrasse 8
Postfach
CH-4614 HÄGENDORF
Tel/Tél: +41 (0) 61 337 97 77
Fax: +41 (0) 61 311 97 42

SUOMI

Arjo Scandinavia AB
Riihitontuntie 7 C
02200 Espoo
Finland
Puh: +358 9 6824 1260
E-mail:
Asiakaspalvelu.finland@arjo.com

SVERIGE

Arjo International HQ
Hans Michelsensgatan 10
SE-211 20 MALMÖ
Tel: +46 (0) 10 494 7760
Fax: +46 (0) 10 494 7761
E-mail: kundservice@arjo.com

UNITED KINGDOM

Arjo UK and Ireland
Houghton Hall Park
Houghton Regis
UK-DUNSTABLE LU5 5XF
Tel: +44 (0) 1582 745 700
Fax: +44 (0) 1582 745 745
E-mail: sales.admin@arjo.com

USA

Arjo Inc.
2349 W Lake Street Suite 250
US-Addison, IL 60101
Tel: +1 (630) 307-2756
Free: +1 (800) 323-1245
Fax: +1 (630) 307-6195
E-mail: us.info@arjo.com

JAPAN

Arjo Japan K.K.
東京都港区虎ノ門三丁目7番8号
ランドィック第2虎ノ門ビル9階
Tel: +81 (0)3-6435-6401
Fax: +81 (0)3-6435-6402
E-mail: info.japan@arjo.com

At Arjo, we believe that empowering movement within healthcare environments is essential to quality care. Our products and solutions are designed to promote a safe and dignified experience through patient handling, medical beds, personal hygiene, disinfection, diagnostics, and the prevention of pressure injuries and venous thromboembolism. With over 6500 people worldwide and 65 years caring for patients and healthcare professionals, we are committed to driving healthier outcomes for people facing mobility challenges.



ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

CE
2797