
INSTRUCTIONS FOR USE

Citadel Plus ENT-ACC101

Trapeze Mount Assembly

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 **WARNING**

To avoid injury, always read this Instructions For Use and accompanied documents before using the product.

Mandatory to read the Instructions For Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2021.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.
The content of this publication may not be copied either whole or in part without the consent of Arjo.

Model **ENT-ACC101** is a lifting pole intended to assist patients in positioning themselves on the bed. **ENT-ACC101** is supplied with an adjustable chain and handle, and is suitable for use with Arjo Citadel Plus Bariatric Care System.

WARNING

Do not exceed the Safe Working Load of 150kg (330 lbs.)

Do not use the lifting poles to push or pull the bed.

Do not suspend other equipment from the lifting poles.

Only use the recommended chain and handle assembly, i.e. part number 20969.

The lifting poles must be used as a set for proper use of ENT-ACC101.

Refer to the bed instructions for use to confirm compatibility before using this accessory.

The lifting poles must be positioned so that the handle assembly hangs over the middle of the bed.

Do not allow the chain and handle assembly to hang outside the perimeter of the bed.

To reduce the risk of patient entrapment, the head end panel must remain fitted at all times while the lifting poles are installed on the bed.

The bed should not be occupied by a patient while the lifting poles are being fitted.



Before installing the lifting poles, check that the bed's lifting pole sockets are not damaged in any way.

Ensure that the lifting poles and handle assembly are positioned to suit the size and medical needs of the patient.


The chain and handle assembly must not be used by persons under twelve years of age.

Inspect the chain, handle assembly and trapeze mount assembly regularly. If any sign of wear or damage is found, remove it from use immediately and replace the complete unit.

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.


Installation and use


Insert the lifting poles into each socket (1) at the head end of the bed, with the horizontal part of the pole above the center of the bed.

Twist the pole so that the notch (2) locates over the pin at the bottom of the socket.

When the lifting pole is correctly installed, the indicator line (3) will be level with the top of the socket.

Remove the threaded knobs (4) from the trapeze mount assembly (5). From the head end of the bed, slide the trapeze mount assembly (5) over the horizontal part of the lifting poles. Ensure that the holes in the top plate are located over the four projecting pegs located on top of the lifting poles. Secure threaded knobs (4) to the trapeze mount assembly (5).

To adjust the height of the handle above the bed: press desired chain loop onto the "s" hook provided.

WARNING

Do not adjust the chain while a patient is holding the handle.

Do not twist the chain through more than 180° as this may cause premature wear.

After use: remove the trapeze mount assembly, lift the poles out of the sockets and store the accessory away from the bed.

Cleaning

WARNING

Remove the trapeze mount assembly from the lifting poles before cleaning.

Only use approved disinfectants and cleaning methods. Other products may cause degradation of the equipment and possibly premature failure.

Use all cleaners and disinfectants with care. Obey the manufacturer's instructions for use.

Do not disinfect the trapeze mount assembly in a pressure chamber, autoclave or steam cleaner as this may damage the product.


Lifting pole: refer to the instructions for use supplied with the bed.

Product weight (including trapeze mount assembly)



ENT-ACC101 = 20.1kg

Spare Parts

Item	Part Number	Description	QTY
1	AH1001111	Assembly, Trapeze Mount (Threaded Knobs Included)	1
2	KKX51060	Threaded Knob	3
3	AH1001055	Assembly, Fixed Position Lifting Pole	2
4	AH1001056	Assembly, Trapeze Handle	1

Symbol explanation

	Refer to Instruction Manual		Manufacturer
	Reference Number		Date of Manufacturer
	Serial Number		Accessory Weight
	CE marking indicating conformity with European Community harmonised legislation		Safe Working Load
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745		Separate electrical and electronic components for recycling in accordance with 2012/19/EC (WEEE)

End of Life Disposal

Components that are primarily made up of different kinds of metal (containing more than 90% metal by weight) for example sling bars, rails, upright supports, etc., should be recycled as metals.

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AUSTRALIA

Arjo Australia
Building B, Level 3
11 Talavera Road
Macquarie Park, NSW, 2113,
Australia
Phone: 1800 072 040

BELGIQUE / BELGIË

Arjo Belgium
Evenbroekveld 16
9420 Erpe-Mere
Phone: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.com

BRASIL

Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelice, 329 PB02
Galpão - Lapa
São Paulo – SP – Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA

Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSIONSSAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA

Arjo Czech Republic s.r.o.
Na Strzi 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK

Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND

Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA

ARJO IBERIA S.L.
Polígono Can Salvatella
c/ Cabanyes 1-7
08210 Barberà del Valles
Barcelona - Spain
Telefono 1: +34 900 921 850
Telefono 2: +34 931 315 999

FRANCE

Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG

Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building,North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

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At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

CE

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Arjo Inc.