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1. European Pressure Injury Advisory panel, National Pressure Injury Advisory Panel, Pan Pacific Pressure Injury Advisory Alliance. The Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA, 2019.
2. European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel & Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: The International Guideline 2019. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA, 2019. Section 8: Introduction to Repositioning and Early Mobilisation. Page 115.
3. Gefen A. How much time does it take to get a pressure ulcer? Integrated evidence from human, animal, and in vitro studies. *Ostomy Wound Manage* 2008;54(10):26-35.
4. Gefen A, Soppi E. (2021) what is new in our understanding of pressure injuries: the inextricable association between sustained tissue deformations and pain and the role of the support surface. *Wound Practice & Research*, Volume 28, No2 - June 2020
5. Coleman S, Nixon J, Keen J et al, A new pressure ulcer conceptual framework. *J Adv Nurs*. 2014; 70(10):2222-2234
6. International Review. Pressure ulcer prevention: pressure, shear, friction and microclimate in context. A consensus document. London: Wounds International, 2010.
7. Le KM, Madsen BL, Barth PW. An in-depth look at pressure sores using monolithic silicon pressure sensors. *Plast Reconstr Surg* 1984; 74(6):745-56.
8. Kottner, J, Balzer K, Dassen T et al, (2009) Pressure ulcers: A critical review of definitions and classifications. *Ostomy Wound Management*, 2009;55(9):22-29
9. Reger S.I, Ranganathan W.K, Orsted, HL, Ohura T, Gefen A (2010) Shear and friction in context. In International Review. Pressure ulcer prevention: pressure, shear, friction and microclimate in context. A consensus document. London: Wounds International, 2010.
10. Stekelenburg A, Gawlitta D, Bader DL, Oomens CW. Deep tissue injury: how deep is our understanding? *Arch Phys Med Rehabil* 2008;89:1410-3.
11. Gefen A, Brienza D, Cuddigan J, Haesler E, Kottner J (2021) Our contemporary understanding of the aetiology of pressure ulcers/injuries. *International Wound Journal*, 2021; 1-13. Available at: <https://doi.org/10.1111/iwj.13667>. Last accessed December 2021.
12. Linder-Ganz E, Engelberg S, Scheinowitz M, et al. Pressure time cell death threshold for albino rat skeletal muscles as related to pressure sore biomechanics. *J Biomech*. 2006; 39(14):
13. Kottner J, Black J, Call E, Gefen A, Santamaria N (2018). Microclimate: A critical review in the context of pressure ulcer prevention. *Clinical Biomechanics*, Volume 59, p62-70.
14. Lachenbruch C, Tzen YT, Brienza DM et al. The relative contributions of interface pressure, shear stress and temperature on tissue ischaemia: a cross sectional pilot study. *Ostomy Wound Manage*. 2013; 59(3):25-34.
15. European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel, Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: The International Guideline 2019. Emily Haesler (Ed.). EPUAP/NPIAP/PPPIA:2019. Section 2: The Aetiology of Pressure Injuries - Contributors to cell damage and tissue necrosis in pressure injuries, Pages 22-23.
16. Gefen A, Ousey K. Update to device-related pressure ulcers: SECURE prevention. COVID-19, face masks and skin damage. *J Wound Care*. 2020 May 2; 29(5):245-259.
17. Gefen A. The future of pressure ulcer prevention is here: Detecting and targeting inflammation early. *EWMA Journal*. 2018; 19(2):7-13.
18. Gefen, A. The Subepidermal Moisture Scanner: the technology explained. *Journal of Wound Care*. 2020; 29(Sup2c):S10-S16.
19. Gefen A, Ousey K Update to device-related pressure ulcers: SECURE prevention. COVID-19, face masks and skin damage. *J Wound Care*. 2020 May 2; 29(5):245-259.

20. National Pressure Injury Advisory Panel. Support Surface Standards Initiative. Terms and Definitions Related to Support Surfaces. November 2019. https://cdn.ymaws.com/npiap.com/resource/resmgr/s3i/10-23_Terms_and_Defs_2019_We.pdf. Last accessed, December 2021.
21. European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel, Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: The International Guideline 2019. Emily Haesler (Ed.). EPUAP/NPIAP/PPIA:2019. Section 10 Support Surfaces. Pages 155-180.
22. Gefen A, Soppi E (2020) The pathophysiological link between pressure ulcers and pain and the role of the support surface in mitigating both. *Wounds International* volume 11, Issue 4, 2020: Pages 38-44. www.woundsinternational.com. Last accessed December 2021.
23. Thurman K, Deppisch M, Morello S, Moseley R, Slayton S, Sylvia C, Todd J (2021) Overcoming the challenge of support surface selection: utilization of standards. *Advances in Skin and Wound Care*. October 2021 : Pages 526-531.
24. Y-K Jan, Brienza DM, Boninger ML & Brenes (2011). Comparison of skin perfusion response with alternating and constant pressures in people with spinal cord injury. *Spinal Cord* (2011) 49, 136-141.
25. Mawson A.R. et al. Risk factors for early occurring pressure ulcers following spinal cord injury. *American Journal of Physical Medical Rehabilitation*. 1988;123-127.
26. Guther R.A. Clark M. The effect of dynamic pressure-redistributing bed support surface upon systemic lymph flow and composition. *Journal of Tissue Viability*. 2000;10 (3 suppl): 10-15.
27. Goossens R.H. Rithalia S.V.S. Physiological response of the heel tissue on pressure relief between three alternating pressure air mattresses. *Journal of Tissue Viability*. 2008; 17(1):10-14.
28. Rithalia S, Heath G.H. Alternating pressure redistribution mattresses: appearances can be deceptive. *MedTex Conference*. 2007.
29. Laird D.A. Did you sleep well? *Rev. Rev.* 1935; 91 (Feb):23-70.
30. Exton-Smith A.N. Sherwin R.W. The prevention of pressure sores – significance of spontaneous bodily movements. *Lancet*. 1961;2 (7212): 1124-1126.
31. Keane F.X. The minimum physiological mobility requirement (MPMR) for man supported on a soft surface. *Paraplegia* 1978; 16:383-389.
32. Rithalia S.V.S, Comparison of performance characteristics of the Nimbus and Airwave mattresses. *International Journal of Rehabilitation and Research*, 1995;18:182-185.
33. European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel, Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: The International Guideline 2019. Emily Haesler (Ed.). EPUAP/NPIAP/PPIA:2019. Section 9 Heel Pressure Injuries. Pages 145-154.
34. Tissue Viability Society. Laboratory measurement of the interface pressures applied by active therapy support surfaces: A consensus document. *Journal of tissue Viability*. 2010; 19(1):2-6.
35. Teleten O, Kirkland-Kyn, H. Pain T, Ballesteros R.J (2019). The use of pressure mapping: An educational report. *Wounds* 2019; 31(1):E5-E8. Available at: <https://www.hmpgloballearningnetwork.com/site/wounds/brief-communication/use-pressure-mapping-educational-report>. Last accessed November 2021.
36. American National Standards Institute/Rehabilitation Engineering and Assistive Technology Society of North America. SS-1-2019 RESNA Standard for Support Surfaces—Volume 1: Requirements and Test Methods for Full Body Support Surfaces. Arlington, VA: Rehabilitation Engineering and Assistive Technology Society of North America; 2019.
37. National Pressure Injury Advisory Panel (NPIAP). Guidance on Interpretation of Performance Standards for Support Surfaces. <https://npiap.com/page/S3i>.
38. Okonkwo H, Bryant R, Milne J et al. A blinded clinical study using subepidermal moisture biocapacitance measurement device for early detection of pressure injuries. *Wound Repair & Reg* 2020;11-11.

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© Arjo, 2022
Arjo AB · Hans Michelsensgatan 10 · 211 20 Malmö · Sweden · +46 10 335 4500