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Quellen

1. Girard R, Baboi L, et al. The impact of patient positioning on pressure ulcers in patients with severe ARDS: results from a multicentre randomised controlled trial on prone positioning. *Intensive Care Med* (2014) 40:397-403
2. Bunker D, Thomson M, Chin Necrosis as a Consequence of Prone Positioning in the Intensive Care Unit Case Reports in Medicine Volume 2015, Article ID 762956, 3 Seiten
3. Guerin C, Reignier J, Richard JC, Beuret P, Gacouin A, Boulain T, Mercier E, Badet M, Mercat A, Baudin O, Clavel M, Chatellier D, Jaber S, Rosselli S, Mancebo J, Sirodot M, Hilbert G, Bengler C, Richecoeur J, Gainnier M, Bayle F, Bourdin G, Leray V, Girard R, Baboi L, Ayzac L, Group PS (2013) Prone positioning in severe acute respiratory distress syndrome. *N Engl J Med* 368:2159-2168
4. Ward C, Wubbels M, Schembri L (2012). Using Complete Pressure Off-loading and Advanced Wound Care to Treat a Complex Sacral Pressure Ulcer. Posterpräsentation
5. Clancy J (2011). Winning the war against pressure ulcers. Posterpräsentation
6. Malbrain M, Hendriks B, Wijnands P et al (2010). A pilot randomised controlled trial comparing reactive air and active alternating pressure mattresses in the prevention and treatment of pressure ulcers among medical ICU patients. *Journal of Tissue Viability* 19(1): 7-15
7. Finnegan MJ (2008). Comparing the effectiveness of a specialised alternating air pressure mattress replacement system and an air fluidised integrated bed in the management of post-operative flap patients. A randomised controlled study. *Journal of Tissue Viability* 17(1): Februar 2008
8. Russell L, Reynolds TM (2000). Randomised controlled trial of two pressure relieving systems. *Journal of Wound Care* 9(2): 52-55
9. Evans D, Land L, Geary A (2000). A clinical evaluation of the Nimbus 3 alternating pressure mattress replacement system. *Journal of Wound Care* 9(4): 181-186
10. Land L, Evans D, Geary A et al (2000). A clinical evaluation of an alternating pressure mattress replacement system in hospital and residential care settings. *Journal of Tissue Viability* 10(1): 6-11
11. Phillips L (2010). Nimbus range of pressure redistributing mattresses. *Wounds UK*, 2010, 6 (2): 116-122
12. The Impact of Cleaning Chemicals on Polyurethane Mattress Cover Materials and their Propensity for Physical Damage, Jo Milnes MA(Hons), MSc(Dist): Technical Manager, Dartex Coatings Limited, Acton Close, Long Eaton, Nottingham, NG10 1FZ. 2013

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© Arjo, 2020
Arjo AB · Hans Michsensgatan 10 · 211 20 Malmö · Schweden · +46 10 335 4500
Arjo Deutschland GmbH · Peter-Sander-Straße 10 · 55252 Mainz-Kastel · Deutschland · +49 6134-186 0
Arjo AG · Fabrikstrasse 8/Postfach · 4614 Hügendorf · Schweiz · +41 (0) 61 337 97 77
Arjo Austria GmbH · Lemböckgasse 49A · 123 Wien · Österreich · +43 (0) 866 56