

WIPEABLE CLIP SLING FOR ARJO STANDING AND RAISING AIDS*

Sara Flex wipeable standing sling



Wipe to
clean

With no laundering requirement our Wipeable Standing Clip Sling is designed to be able to be wiped clean between patients or residents; hence helping to reduce the cross infection risks associated with shared patient handling devices^{1,2}.



Design benefits

A: Completely wipeable - Avoid the laundry process, i.e. no lost or missing slings, the sling is always available. The body, the straps and attachment points of the sling are made of non-porous material and hence completely wipeable.

B: Welded seams - No stitching on the seams of the sling is designed to reduce potential of fluid ingress and contamination of the sling.

C: Smooth patient contact surface - The inside of the sling in contact with the patients, is free from seams, straps and other components, designed to create a smooth surface towards the patients body.

D: The belt is fixed to the sling - The belt and sling are attached to prevent the two being separated and the belt being lost.

E: Colour coded sizes - All slings are size colour-coded for ease when choosing the right sling. There are four sizes available: S, M, L and XL.

*Compatible lifters:

Sara® Steady, Sara® Flex and Sara® 3000

PRODUCT SPECIFICATIONS

General	
*Compatible lifters	Sara® Steady, Sara® Flex and Sara® 3000
Safe Working Load, SWL	200 kg / 440 lbs
Size (Article Number)	Waist size in cm (in)
Small (TSS.500W-L1)	64-84 cm (25" - 33")
Medium (TSS.501W-L1)	76-104 cm (30" - 41")
Large (TSS.502W-L1)	88-124 cm (34 ½" - 49")
Extra Large (TSS.503W-L1)	108-160 cm (42 ½" - 63")
Cleaning/disinfection	
This sling can be wiped down with the following types of cleaning/disinfection agents: Please refer to the "Instructions for use" for more information.	Ethanol, Isopropanol, Hydrogen peroxide, Quaternary ammonium, and Bleach.