



Contemporary bedside furniture that supports patients and caregivers in essential daily activity

DURABLE AND MODERN BEDSIDE FURNITURE

Overbed table & bedside cabinet



High quality furniture for essential daily requirements

Modern and elegant with clean lines and rounded profiles, the bedside furniture is designed and built for the people who use it most.

The Overbed Table has a two-piece tilting table top, allowing the patient to position part of the table at multiple angles for activities such as reading. The book rest is included as standard. The smaller side table section remains flat, providing a practical surface for items such as drinks.

The single hand operation, designed for a smooth variable height adjustment, raises the table from 735 mm to 1100 mm and automatically locks into position providing a stable surface at all times. If the bed is raised inadvertently, the table top will also rise to avoid accidental entrapment.



The Bedside Cabinet has two large drawers easily accessible by the patient whilst in bed.

Practical and flexible

The Bedside Cabinet includes, as standard, bag or towel hooks on each side which are within easy reach for everyday patient use, a lockable top drawer for privacy and an integrated bottle holder rail in the bottom drawer which holds up to three 1.5 litre bottles.

Optional perimeter rails help to provide extra stability for items on the top surface and also assist caregivers when positioning the cabinet and table.



Strong, durable and easy clean

The drawer lining of the Bedside Cabinet is smooth and built with water resistant plastic. Both drawers can be simply removed, without tools, for full access and ease of cleaning.

All Bedside Cabinet models are supplied with four easy to access twin casters, two of which have brakes. The low profile, twin wheel castors on the Overbed Table provide high mobility for easy and accurate positioning around the bed.

Flat surfaces are made from high quality smooth and non-absorbent HPL (high pressure laminate). HPL is renowned for its easy-to-clean qualities, in addition to its durability in daily use and resistance to chemicals, moisture, stains and dirt.



Colours and options



Single colour

Light Grey offers neutral aesthetic.



Colour highlights

Selected colour highlighted on tilt table top side section, bedside cabinet shelf and top surface.



Wood finish

Choice of wood style finishes.



AH Blue



Lilac



Light Blue




Light Grey



White Oak




Copper Beech



North Walnut

Overbed tables are available in 2 models:



February 2019. Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. ® and ™ are trademarks belonging to the Arjo group of companies.
© Arjo, 2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure injury and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

Arjo AB · Hans Michelsengatan 10 · 211 20 Malmö · Sweden · +46 10 335 4500

www.arjo.com