

HEIGHT ADJUSTABLE BATH CHAIR

# Calypso



Calypso is a height adjustable bath chair that enables transporting, transferring in and out of a bathtub and bathing of active to semi-active residents, making a one-to-one assisted bathing process possible.

## ARJO MOBILITY GALLERY



A



B



C

### **Designed to support the bathing workflow**

The height adjustable Calypso bath chair enables a single caregiver to safely and comfortably transfer a resident or patient from the bedside, chair or wheelchair to the bath, shower or toilet, and back again, without any manual lifting. Thus, enabling a more ergonomic working position throughout the care process.

### **Enables bathing in seated or reclined position**

Calypso features individually foldaway back and hand rests. With the backrest folded away the resident can, subject to mobility level, lean back in the tub for a better immersion. The foldaway back and hand rest facilitates transfer on and off the bath chair.

### **Ease of maneuverability**

Calypso bath chair features both positioning and transfer handles with easy grip to facilitate maneuverability during placement and transportation of the bath chair.

# Key features

Calypso comes in two versions, hydraulic or electric.

## BACKREST / HAND RESTS

For resident comfort and safety

## TRANSFER HANDLES

Used during transportation of the resident

## SEAT

Available in two versions: the Classic seat and a version in soft, easy to clean Polyurethane material

## TOILET INTEGRATION

To facilitate toileting. Calypso can be positioned directly over the toilet

## HEIGHT ADJUSTABLE

For an ergonomic work position

## POSITIONING HANDLES

For bed/bath side maneuvering and positioning of the bath chair

## HAND CONTROL

To adjust the height of equipment

## BED PAN AND HOLDER (OPTIONAL)

Enabling in-chair toileting needs

## FOOT SCALE (OPTIONAL)

For easy weighing of resident

## FOOT BRAKE

Manually applied on all castors



PRODUCT INFORMATION	CALYPSO ELECTRIC PUR	CALYPSO ELECTRIC CLASSIC	CALYPSO HYDRAULIC
Max. patient weight, SWL	136 kg (300 lbs.)	136 kg (300 lbs.)	136 kg (300 lbs.)
Total weight of lift (lift + scale)	60 kg (131 lbs.)	56 kg (122 lbs.)	47 kg (103 lbs.)
Total width	641 mm (25 1/4")	641 mm (25 1/4")	631 mm (24 7/8")
Total width seat	472 mm (18 5/8")	479 mm (18 7/8")	479 mm (18 7/8")
Total length	933 mm (36 3/4")	933 mm (36 3/4")	933 mm (36 3/4")
Seat height, lowest position	450 mm (17 3/4")	430 mm (16 7/8")	470 mm (18 1/2")
Seat height, highest position	1015 mm (40")	995 mm (39 1/8")	970 mm (38 1/4")
Min. height of equipment	1130 mm (44 1/2")	1130 mm (44 1/2")	1105 mm (43 1/2")
Max. height of equipment	1695 mm (66 3/4")	1695 mm (66 3/4")	1605 mm (63 1/4")
Max. chassis width	641 mm (25 1/4")	641 mm (25 1/4")	631 mm (24 7/8")
Min. chassis width between legs	459 mm (18 1/8")	459 mm (18 1/8")	484 mm (4")
Max. chassis length	933 mm (36 3/4")	933 mm (36 3/4")	933 mm (36 3/4")
Lift stroke range	565 mm (22 1/2")	565 mm (22 1/2")	500 mm (19 5/8")

ACCESSORIES	
Safety belt	To ensure resident safety. Calypso comes with one safety belt.
Batteries and charger 24V	Quick-replaceable battery system with two batteries, ensuring availability for use at all times.
Bed pan and holder (optional)	To facilitate toileting.
Scale (optional)	Foot-operated scale to weigh resident.

## SPACE REQUIREMENTS

- Blue area shows the minimum working area required for the staff to be able to use the mechanical aids in an ergonomic way from one side.
- Light blue area shows required extension of working area to facilitate activities from either side to provide adequate access.

For more info on installation requirements visit Architect portal at [www.arjo.com](http://www.arjo.com)

