

## SHOWER TROLLEY

# Carevo®



Featuring  
patented  
Flexi Zone and  
Ergo-Access  
areas

The Carevo® shower trolley is designed to provide a comfortable shower experience for patients and residents who need to remain in a lying position, while enabling the caregiver to work efficiently and ergonomically.

Designed with Emma in mind



### Designed for patient comfort

The soft surface material and patented Flexi Zone are designed to enhance patient or resident comfort throughout the hygiene routine.

### Universal control panel

Clearly marked control panels located on both sides and at the head of the trolley, allowing easy access for the carer.

### Promote patient security

Comfort handles provide support for your resident or patient and a greater sense of security while promoting their participation.

### Facilitates improved working posture


Height adjustability and patented Ergo-Access areas enable the caregiver to get closer to the resident and maintain a comfortable working posture, thus reducing static physical overload.\*

\*Source: Knibbe NE, Knibbe JJ & Waaijer E (2013)

# Key features



## SPACE REQUIREMENT



Blue area denotes the required working space to facilitate hygiene and transfer activities while providing adequate access for the resident, equipment and assisting carer.

## PRODUCT INFORMATION

Safe working load (SWL) patient + water	200 kg (440 lb)
Max patient weight	182 kg (400 lb)
Weight (incl mattress)	82 kg (181 lb)
Total length	2005 mm (79 in)
Width (stretcher)	586 mm (23 in)
Total width (incl. side supports)	759 mm (29 ¾ in) – 888 mm (35 in)
Length of chassis	1211 mm (47 ¾ in)
Width of chassis	700 mm (27 ½ in)
Height, top position floor to stretcher	955 mm (37 ½ in)
Height, low position floor to stretcher	605 mm (23 ¾ in)
Stroke range	350 mm (13 ¾ in)

## ACCESSORIES

- Head pillow
- Wedge pillow
- Mattress
- Drainage hose
- Battery charger and battery
- Surface disinfectant

For more detailed information visit [www.arjo.com](http://www.arjo.com)