

BARIATRIC SHOWER COMMODE CHAIR

Carmina



Designed for assisted showering, toileting and hygiene care for plus size patients and residents, Carmina offers a solution for transferring from the bedside, chair or wheelchair to the shower or toilet, and back again.

ARJO MOBILITY GALLERY



A



B



C

A solution for regular showering

Carmina is a mobile shower chair that supports the hygiene workflow, providing secure transport to and from the shower room and enabling hygiene care in a seated position.

Toileting

Carmina is designed to support the patient or resident during toileting routines by providing a commode function when access to the toilet is not possible.

Comfort for patient or resident and caregiver

Designed to accommodate the various body shapes and weight distributions of plus size patients, with comfort and dignity in mind. The adjustable backrest with self-adjusting angle allows a comfortable position to be achieved. Foldable armrests allow good access for the caregiver, whilst the additional push poles support maneuverability.

Key features



ARMREST
The armrests fold up, giving access for caregiver and enables lateral transfers.

ADJUSTABLE BACKREST
Easily adjustable backrest has four positions and self-adjusting angle.

PUSH POLES
The push poles assist the caregiver when manoeuvring the chair.

TIP-UP SEAT
The tip-up seat is designed to make cleaning and disinfection easier and provides access to the commode container.

SELF-ADJUSTING COMMODE CONTAINER HOLDER
Container holder lifts up, self-adjusting to different toilet heights.

ADJUSTABLE FOOTREST
Pedal-operated footrest.

PRODUCT INFORMATION

Maximum Patient Weight, Safe Working Load	320 kg (705 lb)
Total weight, excl. packaging (Carmina Basic / Carmina)	37 / 75 kg (81 / 165 lb)
Maximum total weight (lift + patient)	357 / 395 kg (785 / 869 lb)
Overall width of product	960 / 997 mm (37 3/4 / 39 1/4 inches)
Overall depth of product	760 / 980 mm (29 7/8 / 38 5/8 inches)
Overall height of product	1265 mm (49 3/4 inches)
Seat depths (4 positions)	420 to 525 mm (16 1/2 to 20 5/8 inches)
Seat width	915 mm (36 inches)
Turning diameter	1300 mm (51 1/8 inches)
Ambient temperature range	+10 C to +40 C (+50 F to +104 F)
Water temperature range	+10 C to +40 C (+50 F to +104 F)
Relative humidity range	30% to 75%
Operating force brakes	<300 N

ALLOWED COMBINATIONS

Carmina shower and commode chair integrates with the following Arjo products	Maxi Sky 1000, Maxi Sky 600, Tenor, Maxi Move
--	---

SPACE REQUIREMENTS

- Blue area shows the minimum working area required for the staff to be able to use the mechanical aids in an ergonomic way from one side.
- Light blue area shows required extension of working area to facilitate activities from either side to provide adequate access.

For more info on installation requirements visit Architect portal at www.arjo.com

