

Arjo disposable slide sheets



Arjo disposable slide sheets are patient/resident-specific tubular slides and flat sheets with built-in handles.

ARJO MOBILITY GALLERY



C



D



E

© Arjo 2019

C. Able to partially bear weight on at least one leg, uses a wheelchair for moving a distance.

D. Unable to bear weight through feet. Is able to sit if well supported.

E. Difficulty maintaining sitting, hence might remain lying/partially lying in bed, and may only sit in a special supportive chair.

Please contact Arjo for further information.

- **Cross-contamination prevention**

As a patient-specific product, Arjo disposable slide sheets offer a solution for busy care environments concerned about reducing cross-contamination. Every sliding aid comes individually wrapped in a heat-sealed plastic bag. The product is not designed to be laundered.

- **Prevent caregiver injuries by reducing physical effort**

Assisted in-bed repositioning tasks expose caregivers and patients at risk of injury. Arjo disposable slide sheets are made from a high-strength, latex-free, low-friction material enabling multi-directional movement within the bed for a variety of repositioning tasks. Use of slide sheets for in-bed repositioning has shown to reduce caregiver exposure to excessive push and pull forces in addition to reduce injury to patients associated with shear and friction when repositioning without slide sheets.¹

Key features

All Arjo disposable slide sheets are full bed length, which eliminates the need for additional tubes and flat sheets

PATIENT IDENTIFICATION TAG

Each Arjo disposable slide sheet has a tag to write the patient’s name on and a green water-soluble label, which will dissolve if the sheet is laundered, offering you a visible indication when not to use the sheet.



STRONG BUILT-IN HANDLES

The built-in reinforced handles offer something to hold during patient transfers.



DISPOSABLE PULL STRAPS

The disposable pull straps can be attached to the sheets for lateral transfers avoiding caregivers having to overreach, thereby reducing the risk of injury.

PRODUCT INFORMATION		
Flites	Width	Length
Arjo disposable slide sheet	760 mm / 30 in	1960 mm / 77 in
Arjo disposable slide sheet – XL	1140 mm / 45 in	1960 mm / 77 in
Arjo disposable slide sheet – XXL	1520 mm / 60 in	1960 mm / 77 in
Arjo disposable pull straps for slide sheets	40 mm / 1 ½ in	800 mm / 31 ½ in

Tubes	Width (when tube is flat)	Length
Arjo disposable slide tube – half bed	980 mm / 39 in	980 mm / 39 in
Arjo disposable slide tube	760 mm / 30 in	1960 mm / 77 in
Arjo disposable slide tube – XL	1155 mm / 45 ½ in	1960 mm / 77 in

1. Baptiste A, et al. (April 2006) “Friction-Reducing Devices for Lateral Patient Transfers: A Clinical Evaluation”. Article in AAOHN Journal. DOI: 10.1177/216507990605400407. Source: PubMed.