

# Arjo Slide Sheets



Arjo slide sheets with handles is a sliding aid that assists the caregiver with in-bed repositioning and transfers with dependent patients or residents.

## ARJO MOBILITY GALLERY



C



D



E

© Arjo 2019

**C.** Able to partially bear weight on at least one leg, uses a wheelchair for moving a distance.

**D.** Unable to bear weight through feet. Is able to sit if well supported.

**E.** Difficulty maintaining sitting, hence might remain lying/partially lying in bed, and may only sit in a special supportive chair.

Please contact Arjo for further information.

- **A simple and cost-effective solution**  
The Arjo slide sheets are easy to access and are an economical solution because they do not require a lift to perform a turn or lateral transfer.
- **Prevent caregiver injuries by reducing physical effort**  
Done manually, lateral transfer creates a risk of injury, most likely to the caregiver's shoulders and lower back. Friction reducing devices reduce the pushing and the pulling force from tasks like in-bed repositioning and transfers with dependent patients or residents.<sup>4</sup>
- **Help prevent pressure injuries**  
Pressure injuries prevalence and incidence is high in many healthcare facilities and they have a significant economic impact for the organizations, and also affect a patient's quality of life, morbidity, and mortality. This type of injury can be prevented with mobility.<sup>1-3</sup>
- **Improve patient safety**  
Using a sliding aid while performing an in-bed repositioning or a transfer will improve patient safety.<sup>4</sup>

# Key features

## PULL STRAPS (OPTIONAL)

Optional pull straps can be attached to the sheets for lateral transfers, avoiding caregivers needing to overreach.

## STRONG BUILT-IN HANDLES

Each sheet has built-in reinforced handles for the caregiver to hold. The low-friction device helps the caregiver manage a range of in-bed repositioning and lateral transfers with dependent patients.

## PATIENT COMFORT

Arjo slide sheets are manufactured in a soft fabric, which can be fitted beneath the patient promoting comfort and reducing shear and friction during in-bed repositioning and transfers.<sup>4</sup>



## FLAT SHEET DESIGN

The flat sheet design with handles provides a low-friction surface to help the caregiver manage a range of in-bed repositioning and lateral transfers with dependent patients.

## FLEXIBILITY

The sheets are available in six sizes that cover most patient sizes for the most common applications.

## PRODUCT INFORMATION

	Width		Length	
Arjo slide sheet	730 mm	28 ¾ in	1730 mm	68 in
Arjo slide sheet – long	840 mm	33 in	2000 mm	78 ¾ in
Arjo slide sheet – XL	1095 mm	43 in	1730 mm	68 in
Arjo slide sheet – long XL	1095 mm	43 in	2000 mm	78 ¾ in
Arjo slide sheet – XXL	1460 mm	57 ½ in	1730 mm	68 in
Arjo slide sheet – long XXL	1460 mm	57 ½ in	2000 mm	78 ¾ in
Arjo pull straps for slide sheets	–	–	660 mm	26 in

1. Dealey C, et al. J Wound Care. 2012; 21(6):261-266.
2. Brem H, et al. Am J Surg. 2010; 200:473-477.
3. Spear M. Plastic Surgical Nursing. 2013; 33(3):147-149.
4. Baptiste A, et al. Friction-reducing devices for lateral transfers: A clinical evaluation. AAOHN Journal 2006; 54(4):173-180.