

Maxi Twin[®] Compact



Maxi Twin Compact is a smaller version of Maxi Twin floor lifter, designed to access limited spaces in care facilities.

ARJO MOBILITY GALLERY



D

E

© Arjo 2019

D. Unable to bear weight through feet. Is able to sit if well supported.
E. Difficulty maintaining sitting, hence might remain lying/partially lying in bed, and may only sit in a special supportive chair.
Please contact Arjo for further information.

- **Small size**
The smaller size compared with Maxi Twin enables access to areas where space is limited.
- **Floor pick up**
The high lifting stroke enables lifting from the floor.
- **Manoeuvrability**
The twin-mast design distributes the patient's or resident's weight evenly.
- **Precision transfers**
Integrated with the spreader bar is a powered dynamic positioning system (PDPS) that helps in repositioning the patient in the sling, minimising manual handling.
- **Transfer into car**
Following assessment of the environment and car type, Maxi Twin Compact can be used for transfers into and out of a car.

Key features

TWIN-MAST DESIGN

The twin-mast design distributes the patient's weight evenly. The open mast design gives legroom between the two masts.

ERGONOMIC HANDLING

The Maxi Twin mast design provides convenient positioning handles for the use of caregivers.

SMALL SIZE

The lifter's small footprint (length, height and width) enables access to areas where space is limited.



POWERED DPS

The powered dynamic positioning system (PDPS) minimises manual handling and improves precision when positioning the patient.



Lifter can be used for transfers into and out of a car.

PRODUCT INFORMATION


Max safe working load	160 kg	350 lb
Lifter weight (with spreader bar and battery)	48 kg	103 lb
Max lifting height	1375 mm	54 1/8 in
Min lifting height	290 mm	11 3/8 in
Lifting stroke	1080 mm	42 1/2 in
Min total storage height	1125 mm	44 1/4 in
External width legs closed	660 mm	26 in
External width legs open	1105 mm	43 1/2 in
Internal width leg closed	520 mm	20 1/2 in
Internal width leg open	965 mm	38 in
Leg clearance (floor to top of legs)	115 mm	4 1/2 in
Chassis clearance (floor to bottom of chassis)	30 mm	1 1/4 in
Lifter length	1100 mm	43 in
Turning diameter	1310 mm	51 1/2 in
Lifter protection class	IP24	
Hand control protection class	IPX7	
Battery	24V, 2.5 Ah NiMH	
Battery charge indicator		
Emergency stop and system failure override		
Automatic safety cut-out if lowered onto an obstacle		
Low-friction castors, the rear two with brakes		

ACCESSORIES

Wide range of special purpose slings and flites
Arjo disinfectants

SPACE REQUIREMENT

The space requirements in the illustration are based on a Maxi Move passive floor lifter. Please note that Maxi Twin Compact has an even smaller footprint.



- Blue area shows the minimum working area required to enable the caregiver to use the floor lifter ergonomically from one side.
- Light blue area shows extension of working area required to facilitate activities from either side to provide adequate access for the patient, floor lifter and caregiver.