

CLEANING AND DISINFECTION PANEL

MultiClean



MultiClean is a wall-mounted disinfectant mixing and dispensing system which can be used to clean and disinfect bathing and showering equipment, such as our range of fixed or height-adjustable tubs, hygiene lifters, hygiene chairs and shower trolleys.

Easy to use

Instructions on the front of the panel to guide the correct procedure. Selection of either the disinfection or rinse function is marked on the control panel.

Safe and efficient

The unit has a lockable door to restrict unauthorized access to the main water isolation valve and the Arjo cleaner/disinfectant. Using Arjo range of disinfectants, the concentration is preset by the user, ensuring that the correct dilution ratio is delivered.

Versatile

MultiClean can be used to disinfect the air spa system on existing products as well as the air spa system on baths featuring the "Twist Lock" disinfection port.

PRODUCT INFORMATION

Weight, excluding disinfectant	11 kg (24 1/2 lbs) standard 13 kg (29 lbs) with pump
Weight, including disinfectant	14 kg (31 lbs) standard 16 kg (35 1/2 lbs) with pump
Dimensions (W x D x H)	484 x 200 x 578 mm (19 x 8 x 22 3/4 in)
Inlet coupling (cold water)	15 mm (1/2 in)
Overflow	15 mm (1/2 in) standard 19 mm (3/4 in) with pump
Min flow rate – inlet	6 l/min (1.6 Gal/min)
Min flow rate – outlet	5 l/min (1.3 Gal/min)
Min dynamic pressure	2 bar (29 Psi) standard 0.3 bar (2.9 Psi) with pump
Max static pressure	6 bar (87 Psi)
Ambient temperature	-20 to 70 C (-4 to 158 F)
Relative humidity	10–80% at 20 C (68 F)