

SPECIFICATION SHEET

Prioma® Birthing Bed



Fully integrated patient and nurse control panels



Large under tray



Ergonomically designed, sturdy and reliable leg stirrups



Retracting backrest

Reliability and comfort for one of life's greatest moments.

The Prioma Birthing Bed is an electrically operated delivery bed designed to provide the ultimate comfort and support to the mother, as well as efficient and reliable functionality for the delivery team throughout all stages of the birthing process.

- Safe Working Load (SWL) 200 Kg
- 150 mm twin central locking castors
- Linak actuators
- ABS removable backrest panel
- High pressure compact laminate seat & foot panels
- Head end drop down side rails
- Integral bed extension - 580 mm
- Battery back up
- Manual CPR backrest quick release
- Removable lightweight head & foot panels
- Trendelenburg tilt 14°, reverse Trendelenburg tilt 12°

Standard accessories

A full range of attachable accessories come as standard to optimise positioning for delivery and clinical procedures.

- Under tray
- Grab handles
- Footrest plate
- Upper and lower leg supports
- Leg stirrups
- Mattress sections

Optional features

- Under bed light - set in permanent on-mode
- 125 mm twin central locking castors
- Linak nurse control handset

Optional accessories

- Storage trolley - powder coated finish
- Self help pole - includes strap and handle, available in stainless steel or powder coated finish
- IV pole - stainless steel finish
- "C" size oxygen bottle holder - clamps to head or foot end of bed
- Extra set of mattress covers
- Dust cover



Optional storage trolley

Basic Dimensions

Overall length - bed configuration	2130 mm
Overall length - obstetric configuration	1580 mm
Bed extension length	580 mm
Mattress surface length - between panel faces	2010 mm
Overall width	1000 mm
Deck height range - 150mm twin central locking castors	485 mm to 850 mm +/- 5 mm
Trendelenburg and reverse Trendelenburg tilt - electronic	0 to 14° (12° reverse Trendelenburg)
Backrest angle adjustment - electronic	0 to 70°
Mattress size seat and backrest section	850 mm (w) x 1350 mm (l) x 125 mm (h)
Mattress size extension bolster	850 mm (w) x 750 mm (l) x 170 mm (h)

Product Weight (approximate)

Bed without accessories	160 Kg
Mattress sections	14 Kg
Safe Working Load (SWL)	200 Kg
Maximum patient weight	135 Kg

Electrical Data

Power in - 100 - 240V ~ 50/60Hz Max 3.9A
Duty rating - intermittent 10%, max 2 min / 18 min
Electric Shock Protection - Class 1 Type B
Liquid ingress protection - IPX4
Internally powered equipment - battery backup only - 2 x 12V series connected batteries, sealed chargeable lead/Acid Gel 1.2 amp hours

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind. www.arjo.com

Arjo AB, Hans Michelsengatan 10, SE-211 20 Malmö, Sweden. Phone: +46 (0) 10 335 45 00

www.arjo.com

ARJO

Arjo.A00566.2.0.INT.EN