

WIPEABLE CLIP SLING FOR ARJO STANDING AND RAISING AIDS*

Arjo Wipeable Standing Clip Sling



With no laundering requirement our Wipeable Standing Clip Sling is designed to be able to be wiped clean between patients or residents; hence helping to reduce the cross infection risks associated with shared patient handling devices^{1,2}.



Design benefits

A: Completely wipeable - Avoid the laundry process, i.e. no lost or missing slings, the sling is always available. The body, the straps and attachment points of the sling are made of non-porous material and hence completely wipeable.

B: Welded seams - No stitching on the seams of the sling is designed to reduce potential of fluid ingress and contamination of the sling.

C: Smooth patient contact surface - The inside of the sling in contact with the patients, is free from seams, straps and other components, designed to create a smooth surface towards the patients body.

D: The belt is fixed to the sling - The belt and sling are attached to prevent the two being separated and the belt being lost.

E: Colour coded sizes - All slings are size colour-coded for ease when choosing the right sling. There are four sizes available: S, M, L and XL.

*Compatible lifters:

Sara® Steady, Sara® Flex and Sara® 3000

PRODUCT SPECIFICATIONS

General	
*Compatible lifters	Sara® Stedy, Sara® Flex and Sara® 3000
SWL	200 kg / 440 lbs
Size (Article Number)	Waist size in cm (in)
Small (TSS.500W-L1)	64-84 cm (25" - 33")
Medium (TSS.501W-L1)	76-104 cm (30" - 41")
Large (TSS.502W-L1)	88-124 cm (34 ½" - 49")
Extra Large (TSS.503W-L1)	108-160 cm (42 ½" - 63")
Cleaning/disinfection	
This sling can be wiped down with the following types of cleaning/disinfection agents: Please refer to the "Instructions for use" for more information.	Ethanol, Isopropanol, Hydrogen peroxide, Quaternary ammonium, and Bleach.

Arjo.A001611.0.INT.EN November 2019. Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.
20P_Bro-202007-0130-ANZ