



From
care facility
to home

CARING FOR VETERANS

Government product catalog

October 2025

VETERANS MOBILITY GALLERY & VA MSST	3
PATIENT HANDLING SOLUTIONS	
Sit-to-stand lifts	4
Floor lifts	5
Ceiling lifts	6
SLING SOLUTIONS	
In-bed & lateral transfer slings	7
Seated transfer slings	8
Seated transfers, standing & walking slings	9
LATERAL TRANSFERS SOLUTIONS	
Sliding sheets	10
Air-assisted transfer mattresses	11
HYGIENE SOLUTIONS	
Bathing systems	12
Bath lifts	13
Shower systems	14
VENOUS THROMBOEMBOLISM (VTE) PREVENTION	
Deep Vein Thrombosis solutions	16
MEDICAL BEDS	17
THERAPEUTIC SUPPORT SURFACES	
Non-powered foam mattresses	18
Powered mattresses	19
THERAPEUTIC SEATING SOLUTIONS	20

Disclosure:

These Arjo solutions are to be used as a guideline. Individual patient assessment should be carried out to determine appropriateness of use. It is not to intended to replace clinical judgment.

The Veteran Mobility Gallery

The Arjo Veterans Mobility Gallery™ serves as a strategic tool, adapted to support the Veterans Affairs which aligns to the VA Mobility Screening and Solution Tool. It is categorized into five distinct levels to aid in the understanding and support of veterans' mobility needs. This initiative helps in identifying and communicating the specific mobility assistance requirements of veterans, thereby enhancing care outcomes.

This gallery is a homage to the courage, resilience, and diverse needs of our veterans from all branches of the military – Marine Corps, Air Force, Navy, Army, and the Coast Guard – illustrated through the personalized mobility journeys of five representative characters.

The Veterans Administration Mobility Screening and Solutions Tool (VA MSST)

The VA MSST has four levels. These levels help determine the appropriate use of powered and non-powered equipment based on a patient's independence.¹

Level 1: Requires the most assistance, often needing powered equipment.

Level 2: Requires moderate assistance, may use a mix of powered and non-powered equipment.

Level 3: Requires minimal assistance, primarily non-powered equipment.

Level 4: Independent, no equipment needed.

Reference

1. Reliability and validity of the Veterans Administration Mobility Screening and Solutions Tool | BMC Health Services Research | Full Text (biomedcentral.com).



A - Albert

Ambulatory but may use a walking cane/crutch for support.



B - Barbara

Ambulatory to some degree but may use a walking frame/crutches for support.



C - Carl

Able to partially bear weight on at least one leg, uses a wheelchair for mobilizing a distance.



D - Doris

Unable to bear weight through her feet. Is able to sit if well supported.



E - Emma

Difficulty maintaining sitting hence might remain lying/partially lying in bed, may only sit in a special supportive chair.

Patient handling solutions

Sit-to-stand lifts



B



MSST

Level 3N



HC

HOME CARE APPROVED

Sara® Steady non-powered sit-to-stand

- Enables a single caregiver to assist patients and residents to perform sit-to-stand transfers throughout the day
- Weight limit: 400 lbs

Product description	Product #	Contracts	
Sara Steady	NTB2000	FSS ECAT NX EQ	36F79725D0231 SPE2DH-25-D-0017 36C10G21D0037
Sara Steady Compact	NTB3000	OPEN MARKET	



C



MSST

Level 2, Level 3P



HC

HOME CARE APPROVED

Sara® Flex powered sit-to-stand

- Powered sit-to-stand device that enables a single caregiver to assist a patient or resident from a seated to a standing position
- Weight limit: 440 lbs

Product description	Product #	Contracts	
Sara Flex	HEB0000-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Sara Flex with scale	HEB1000-01		



C



MSST

Level 2, Level 3P

Walker ambulatory support aid

- Walker provides the support for taking that first step – for both patient or resident and physiotherapist
- Weight limit: 300 lbs

Product description	Product #	Contracts	
Lift Walker	GCB4100-031	OPEN MARKET	
Lift Walker, battery powered	GCB4103-031	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017



C



MSST

Level 2, Level 3P



HC

HOME CARE APPROVED

Sara® Plus powered sit-to-stand

- A powered sit-to-stand designed for active transfers, balance, stepping and gait training
- Weight limit: 420 lbs

Product description	Product #	Contracts	
Sara Plus	HEP0001-US	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Sara Plus with scale	HEP1001-US		

Patient handling solutions

Floor lifts



Maxi Move® floor lift (Item shown KMCLUN-D)

- Maxi Move floor lift is designed to enable a single caregiver to manage demanding everyday patient or resident transfer and repositioning tasks. It is a versatile solution that can be adapted using a variety of spreader bars to accommodate patient transfer needs
- Available with or without powered Dynamic Positioning System (DPS) spreader bar
- Weight limit: 500 lbs

Product description	Product #	Contracts	
Lift, scale, extra low castors* and no spreader bar	KMCEUN	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Lift, scale, extra low castors* and powered DPS	KMCEUN-D		
Lift, extra low castors* and powered DPS	KMCEXN-D		
Lift, scale, low castors* and manual DPS	KMCLUN		
Lift, scale, low castors* and powered DPS	KMCLUN-D		
Lift, low castors* and manual DPS	KMCLXN		
Lift, low castors* and powered DPS	KMCLXN-D		
Lift, scale, standard castors and manual DPS	KMCSUN	OPEN MARKET	
Lift, standard castors and medium manual DPS	KMCXN		
Lift, standard castors and medium powered DPS	KMCSXN-D		
Lift, scale, standard castors and powered DPS	KMCSUN-D	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017

*Castors heights:
Standard castors - 4.5"
Low castors - 4"
Extra low castors - 2.25"



Maxi Twin® floor lift (Item shown KTBB2ESX2WW)

- The open design of Maxi Twin can enhance safe working practices for caregivers when transferring dependent patients or residents. Maxi Twin Compact can be used for transfers into and out of a vehicle
- Available with or without powered Dynamic Positioning System (DPS)
- Select between a 2-point spreader bar, 4-point spreader bar
- Weight limit: 400 lbs

Product description	Product #	Contracts	
Lift with 2-point spreader bar – loop sling	KTBB2ESX2WW	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Lift with 4-point large spreader bar – loop sling	KTBB2FSX2WW		
Lift with manual DPS spreader bar – clip sling	KTBB4ASX2WW	OPEN MARKET	
Lift, scale & powered DPS spreader bar – clip sling	KTBB4BSU2WW	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Lift & powered DPS spreader bar – clip sling	KTBB4BSX2WW		



Tenor® bariatric floor lift

- Tenor bariatric mobile floor lift is specifically designed for lifting and transferring plus-size patients or residents
- Weight limit: 704 lbs


Product description	Product #	Contracts	
Tenor	KHA1000-US	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Tenor with scale	KHA1010-US		

Patient handling solutions

Ceiling lifts



Notice: Ceiling lifts require professional installation and track systems. Please discuss with your Arjo Representative for more details.



Maxi Sky® 2 ceiling lift (Item shown MS000-02-01)

- Maxi Sky 2 ceiling lift delivers a complete and intuitive patient handling solution assisting in safeguarding the wellbeing of caregivers, patients and residents
- Available with or without powered Dynamic Positioning System (powered DPS)
- Optional Enhanced Chagrining System (ECS)
- Weight limit: 600 lbs

Product description	Product #	Contracts	
Motor, 2 function*, ECS and Gate-US	MS230-01-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-001
Motor, 2 function*, ECS and Gate-PDPS-U	MS230-02-01		
Motor, 2 function*, ECS	MS200-01-01		
Motor, 4 function*, ECS	MS400-01-01		
Motor, 4 function*, ECS and powered DPS	MS400-02-01		
Motor and 4 function*	MS410-01-01		
Motor, 2 function*, ECS, MB & powered DPS	MS200-02-01		
Motor, 2 function*, C-Stat, Reg. HC & MB	MS210-01-01		
Motor, 2 function*, MB, C-Stat & powered DPS	MS210-02-01		
Motor, 4 function*, C-Stat and powered DPS	MS410-02-01		

*2 function includes up/down movement of motor, the four function includes up/down and left/right movement of motor



Maxi Sky® 2 Plus bariatric ceiling lift (Item shown MS205-10-01)

- Maxi Sky 2 Plus has a 4-point telescopic spreader bar that automatically adjusts to the patient's size and position. It can be used in combination with a wide range of bariatric slings to lift, transfer, turn and reposition plus-size patients
- Weight limit: 1,000 lbs

Product description	Product #	Contracts	
Lift, set of 2 motors	MS000-01-01	OPEN MARKET	
Lift, set of 2 motors, Infection Control and 2 function*	MS205-10-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017

*2 function includes up/down movement of motor, the four function includes up/down and left/right movement of motor



Maxi Sky® 2 Infection Control ceiling lift (Item shown MS205-01-01)

- To support an optimized ceiling lift cleaning and disinfection process, the Maxi Sky 2 Infection Control offers non-porous components: lifting strap, emergency cord and spreader bar
- Available with or without powered Dynamic Positioning System (powered DPS)
- Optional Enhanced Charging System (ECS)
- Weight limit: 600 lbs

Product description	Product #	Contracts	
Motor, Infection Control, 2 function* and ECS	MS205-01-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Motor, Infection Control, 4 function* and ECS	MS405-01-01		
Motor, Infection Control, 2 function*, ECS and powered DPS	MS205-02-01		
Motor, Infection Control, 4 function*, ECS, MB and powered DPS	MS405-02-01		
Motor, Infection Control, 4 function*, C-Stat, MB and powered DPS	MS415-02-01		
Motor, Infection Control and 2 function*	MS215-01-01		
Motor, Infection Control, 2 function* and powered DPS	MS215-02-01		
Motor, Infection Control and 4 function*	MS415-01-01		
Motor, Infection Control, 2 function*, MB, ECS and GATE-US	MS235-01-01		
Motor, Infection Control, 2 function*, ECS and GATE-IC-PDPS-US	MS235-02-01		

*2 function includes up/down movement of motor, the four function includes up/down and left/right movement of motor



Maxi Sky® 440 portable ceiling lift

- The Maxi Sky 440 is an easy-to-operate portable ceiling lift system, which allows a single caregiver to perform transfers of patients, avoiding the stress and strain associated with manual lifting
- Weight limit: 440 lbs

Product description	Product #	Contracts	
Lift, 2-point portable cassette	LE00009	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017

Sling solutions

In-bed & lateral transfer slings



Repositioning slings

Product description	Product #	Size	Contracts	
Loop repositioning sling (SWL 600 lbs)	MAA6000-L1	One size	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Loop comfort repositioning sling - regular size (SWL 600 lbs)	MAA7000R-GB			
Loop comfort repositioning sling - wide size (SWL 600 lbs)	MAA7000W-GB			
Disposable loop repositioning sling (SWL 600 lbs)	AHD001-L1			



Mortuary transfer device

Product description	Product #	Contracts	
Mortuary bariatric spreader bar kit with 2x72" & 2x47" bands (SWL 600 lbs)	700.05395	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Mortuary standard spreader bar	700.05380		
Mortuary standard large spreader bar	700.05385		
Mortuary spreader bar kit with 2x57" & 2x47" bands (SWL 300 lbs)	700.05390		
Scale adapter for 4-point spreader bar	700.05420		

WHEN ORDERING:

To place a sling order, use the **product #** - your selected **sling size** - then **L1** as shown in the example below.

EXAMPLE:

When ordering a **seated transfer sling**, general purpose, clip-padded legs, size **XL**

MAA4000M-XL-L1

Product # Size L1

Sling solutions

Seated transfer slings



General purpose slings

Product description	Product #	Size	Contracts
Clip sling with padded legs (SWL 600 lbs)	MAA4000M	S,LL,XXL	OPEN MARKET
Loop sling with padded legs (SWL 600 lbs)	MLA2000	S,M,L,XL	
Loop in situ* sling with padded legs (SWL 418 lbs)	MLAAS2000M	S, M, L, XL	
Disposable 4-point clip sling, TAA compliant (SWL 600 lbs)	MFA1000M	S,M, L,LL,XL, XXL	
Disposable loop sling (SWL 600 lbs)	MFA2100	S, M, L, XL, XXL	

*Can be left under the patient between transfers if clinical assessment allows.



Bathing slings

Product description	Product #	Size	Contracts
Clip mesh bathing sling with padded legs (SWL 600 lbs)	MAA4061M	XS, S, M, L,LL, XL	OPEN MARKET



Toileting slings

Product description	Product #	Size	Contracts
Clip toilet sling with padded legs (SWL 600 lbs)?	MAA4031M	XS, S, M, L, LL, XL, XXL	OPEN MARKET
Loop toilet sling with padded legs (SWL 600 lbs)	MLA4031	S, M, L, XL	



Bariatric slings

Product description	Product #	Size	Contracts
Bariatric loop comfort sling (SWL 1,000 lbs)	MAA8000	M, L, XL, XXL	OPEN MARKET
Bariatric loop sling (SWL 1,000 lbs)	MAA8010		
Bariatric loop hammock sling with divided legs (SWL 1,000 lbs)	MAA8020		
Bariatric loop hammock sling (SWL 1,000 lbs)	MAA8030		

WHEN ORDERING:

To place a sling order, use the **product #** - your selected **sling size** - then **L1** as shown in the example below.

EXAMPLE:


When ordering a **seated transfer sling , general purpose, clip-padded legs, size XL**

MAA4000M-XL-L1

Product # Size L1

Sling solutions

Seated transfer, standing & walking slings



Amputee slings

Product description	Product #	Size	Contracts
Amputee clip sling (SWL 600 lbs)	MAA2090M	XS, S, M, L, LL, XL	OPEN MARKET
Amputee right leg clip sling (SWL 600 lbs)	MAA4070M	M, L, LL, XL	
Amputee left leg clip sling (SWL 600 lbs)	MAA4080M	M, L, LL, XL	
Amputee clip bathing sling (SWL 600 lbs)	MAA7090M	M, L, LL, XL	
Amputee loop hammock sling (SWL 600 lbs)	MLA7000	S, M, L, XL	
Amputee clip toilet sling (SWL 600 lbs)	MAA2091M	M, L, XL	

Sara Flex® slings

Product description	Product #	Size	Contracts	
Sara Flex standing sling (SWL 440 lbs)	TSS.500	S	OPEN MARKET	
Sara Flex standing sling (SWL 440 lbs)	TSS.501	M		
Sara Flex standing sling (SWL 440 lbs)	TSS.502	L		
Amputee clip bathing sling (SWL 600 lbs)	TSS.503	XL		
Amputee loop hammock sling (SWL 600 lbs)	TSS.504	XXL		
Sara Flex seated transfer sling (SWL 440 lbs)	TSS.511	S, M, L, LL, XL, XXL		
Sara Flex disposable standing sling (SWL 440 lbs)	MFA3000-L1	One size	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017



Ambulation slings

Product description	Product #	Size	Contracts
Loop walking jacket (SWL 418 lbs)	MAA5000	S, M, L, XL	OPEN MARKET
Loop walking sling (SWL 440 lbs)	TEM	S, M, L	

WHEN ORDERING:

To place a sling order, use the **product #** - your selected **sling size** - then **L1** as shown in the example below.

EXAMPLE:

When ordering a **seated transfer sling**, **general purpose**, **clip-padded legs**, **size XL**

MAA4000M-XL-L1

Product # Size L1

Lateral transfer solutions

Sliding sheets



Washable slide sheets

- Arjo slide sheets are made of a premium low friction fabric that can be laundered and re-used
- The sewn-in handles gives a comfortable grip for pulling and the Arjo pull straps can easily be attached to the handles when extended reach is necessary

Product description	Product #	Size	Contracts	
68" X 28 3/4" Arjo slide sheet	NSA0500-INT4	N/A	FSS	36F79725D0231
Extra large 68" X 43" Arjo slide sheet	NSAV5260-INT4	XL	ECAT	SPE2DH-25-D-0017
Long extra large 78 3/4" X 43" Arjo slide sheet	NSAV5320-INT3	XL Long	OPEN MARKET	
68" X 57 1/2" extra, extra large Arjo slide sheet	NSAV5000-INT4	XXL	FSS	36F79725D0231
			ECAT	SPE2DH-25-D-0017
78 3/4" X 57 1/2" extra, extra large, long Arjo slide sheet	NSAV5330-INT3	XXL Long	OPEN MARKET	
Accessories				
Washable 26" pull straps - set of 8	NSA0510-INT4		OPEN MARKET	



Disposable slide sheets and tubes

- Arjo disposable slide sheets for single patient use are designed to accompany a patient throughout the hospital stay. They are available in different designs and sizes both as flat slide sheets with handles and slide tubes

Product description	Product #	Size	Contracts	
30"X 77" Arjo slide sheet, box of 20	NPA0500	N/A	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
45"X 77" extra large Arjo slide sheet, box of 20	NPAV5260	XL		
60"X 77" extra, extra large Arjo slide sheet, box of 20	NPAV5000	XXL		
30"X 77" Arjo slide tube, box of 20	NPA0700	N/A		
39"X 39" half bed Arjo slide tube, box of 50	NPA0900	N/A	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
45"X 77" extra large Arjo slide sheet, box of 20	NPA0750	XL	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Accessories				
30"X 77" disposable pull straps, box of 50	NPA0510		FSS ECAT	36F79725D0231 SPE2DH-25-D-0017



Washable slide tubes

- Arjo slide tubes are made of a premium low friction fabric that can be laundered and re-used. The products have two surfaces that have been sewn together into a tubular design, which means that only one product is needed rather than two flat slide sheets to achieve the same friction reduction

Product description	Product #	Size	Contracts	
23" X 17" Arjo slide tube small	NSA0600-INT4	S	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
22 13/16" X 26 3/8" Arjo slide tube Radiology	NSAV5060	N/A	OPEN MARKET	
77" X 17" Arjo slide tube long	NSA0700-INT4	N/A	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
35 1/2" X 82 11/16" Arjo slide tube extra long	NSA0800-INT4	XL Long		
35 1/2" X 43 3/8" Arjo slide tube - half bed	NSA0900-INT4	N/A		

Lateral transfer solutions

Air-assisted transfer mattresses



AirPal® disposable air-assisted transfer mattress

- AirPal supports the patient on a bed of air – an elevating effect created by air flowing through the perforated underside of the AirPal transfer mattress, lifting the patient away from the bed surface
- Disposable for single patient use to help reduce the risk of infection and lower the laundry burden
- Two safety belts to secure the patient onto AirPal
- Weight limit: 1,000 lbs

Product description	Product #	Contracts	
AirPal SPS long disposable transfer mattress - 10 per case	AP-028SPSLH-DR	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
AirPal pad full SPS Long disposable transfer mattress 34"x78" - 10 per case	AP-034SPS-DR		
AirPal pad short disposable transfer mattress 34"x42" - 10 per case	AP-034SPSH-DR		
AirPal pad full sps long disposable transfer mattress 39"x78" - 10 per case	AP-039SPS-DR		
AirPal pad short disposable transfer mattress sps 39"x42" - 10 per case	AP-039SPSH-DR		
AirPal disposable sling mattress - 39"	AP-039SPSS-DR		
AirPal pad full straps - 50"x78" - 5 per case	AP-050SPS-DR		
AirPal disposable sling mattress - 50"	AP-050SPSS-DR		



AirPal® launderable air-assisted transfer mattress

- AirPal supports the patient on a bed of air – an elevating effect created by air flowing through the perforated underside of the AirPal transfer mattress, lifting the patient away from the bed surface
- A launderable transfer mattress, made of a durable nylon fabric. It can be wipe cleaned or laundered at 158°F (70°C)
- Two safety belts to secure the patient onto AirPal
- Weight limit: 1,200 lbs

Product description	Product #	Contracts	
Long launderable transfer mattress 34"x78"	AP-034N-DR	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Long launderable transfer mattress 39"x78"	AP-039N-DR		
Long launderable transfer mattress 50"x78"	AP-050N-DR		



AirPal air supply

- For use with all AirPal transfer mattresses
- Air supply 1100W, 120V/60Hz
- Standard air filter
- The air supply unit has an extended hose, making it useful during MRIs and Radiolucent diagnostic procedures

Product description	Product #	Contracts
Air supply – 7" diameter x 12 1/2" height	MAS000001-US	OPEN MARKET
Extended air hose	8664411	
Transport cart	MCA003	

Hygiene solutions

Bathing systems



Century® sit bath

- The Century sit bath is a fixed-height keyhole shaped tub to aid good immersion. This bath is a space-saving alternative for traditional bathing. It is designed to enable transfers and the keyhole-shape gives the caregiver easier access to the patient during bathing sessions
- Mobility categories Albert, Barbara and Carl may use Century with an Alenti or Calypso bath lift

Product description	Product #	Contracts	
Sit bath	AK61611-US	Whirlpool bath contract	36C10G23D0030
Bath with hydromassage and disinfection	AK62611-US	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Bath with HydroSound® and disinfection	AK63611-US	Whirlpool bath contract	36C10G23D0030



Parker® bath

- Featuring a distinctive swinging door, Parker is designed to enable independent to assisted transfers into and out of the tub. Parker is height-adjustable and goes from sitting to reclining, allowing the caregiver and resident to choose their preferred position to enhance the bathing experience
- The Doris mobility category can use Parker with a ceiling lift

Product description	Product #	Contracts	
Bath with air spa	AL27010-US	ECAT	SPE2DH-25-D-0017
Bath with air spa and Sound & Vision® system	AL27110-US	Whirlpool bath contract	36C10G23D0030



Rhapsody® bath

- Rhapsody offers a range of bath that are fully compatible with our range of bath lifters, giving you the flexibility to choose the right solution for your resident or patient. All baths are fully height-adjustable, ensuring you are always able to work in a safe and comfortable position and can enjoy a positive one-to-one interaction throughout your bathing task

Product description	Product #	Contracts	
Rhapsody electronic P220 Hydromassage with integrated liquids & disinfection, single shower, 1-level autofill & emergency lowering	AR32212US1111	OPEN MARKET	
Rhapsody programmable P300 Hydromassage with integrated liquids & automatic disinfection, dual showers, 3 level autofill & emergency lowering	AR32312US0111		
Rhapsody P220 Hydrosound® with liquids & disinfection, dual shower handles, autofill, emergency lowering	AR33212US0011		
Rhapsody P300 Hydrosound with liquids & disinfection, dual shower handles, disinfection lockout contact time, autofill, emergency lowering	AR33312US0011	Whirlpool bath contract	36C10G23D0030
Rhapsody P300 Hydromassage with Sound & Vision®, with liquids & disinfection, dual shower handles, automatic disinfection, autofill & emergency lowering	AR35312US01110	OPEN MARKET	

Hygiene solutions

Bath lifts



Miranti® height-adjustable bath trolley

- Fully compatible with a range of Arjo recumbent baths, Miranti is designed for the needs of patients and residents with the mobility level of Doris and Emma, who need to remain in a lying position throughout the bathing routine

Product description	Product #	Contracts	
Battery powered	CEB8203-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Battery powered with scale	CEB8303-01	OPEN MARKET	
Battery powered bath trolley without scale	CEB8203-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Battery powered bath trolley with scale	CEB8303-01		
Power Drive bath trolley without scale	CEB8103-01		
Power Drive bath trolley with scale	CEB8003-01		



Alenti® hygiene chair

- Alenti enables a single caregiver to collect the patient or resident from their bedside or wheelchair and transfer them directly into the bath and back again. Alenti can also be used for showering, toileting and other hygiene tasks, reducing the number of required transfers for daily hygiene care routines
- Mobility category of Albert, Barbara and Carl may use Alenti to enter and exit the recumbent sit bath

Product description	Product #	Contracts	
Battery powered hygiene chair without scale	CDB8103-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Battery powered hygiene chair with scale	CDB8153-01		



Calypso® hygiene chair

- Calypso is a height adjustable bath chair that enables transporting, transferring in and out of a bathtub and bathing of active to semi-active residents, making a one-to-one assisted bathing process possible
- Mobility category of Albert, Barbara and Carl may use Calypso to enter and exit the recumbent or sit bath

Product description	Product #	Contracts	
Battery powered, classic vinyl padded seat hygiene chair	CDB6003-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Battery powered, PUR seat hygiene chair without scale	CDB7003-01		
Battery powered, PUR seat hygiene chair with scale	CDB7023-01		



Bolero® height-adjustable bath trolley

- Bolero is a height adjustable bath trolley that can be used for transporting, transferring in and out of the bathtub and bathing dependent residents in reclined or fully lying position
- Mobility category of Albert, Barbara and Carl may use Bolero to enter and exit the recumbent sit bath

Product description	Product #	Contracts	
Battery powered bath trolley without scale	CEB6003-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Battery powered bath trolley with scale	CEB6023-01	ECAT	SPE2DH-25-D-0017

Hygiene solutions

Shower systems



Carendo® hygiene showering chair

- The Carendo multipurpose hygiene chair facilitates an efficient hygiene and showering routine within a single transfer, while allowing a comfortable working posture for the caregiver throughout the process
- Carendo is designed to enable a single caregiver to perform the full hygiene routine, including dressing and undressing, toileting, showering and other hygiene tasks in a single transfer

Product description	Product #	Contracts	
Hygiene showering chair (standard accessories seat cushion, head pillow and safety belt)	BIB2003-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Carendo kit (BIB2003-01, back cushion, bedpan holder and bedpan with lid)	BIB2005-01	FSS ECAT NXEQ	36F79725D0231 SPE2DH-25-D-0017 36C10G22D0034



Carino® hygiene showering chair

- Designed to reduce the number of transfers required for assisted showering and hygiene, Carino enables a single caregiver to perform everyday hygiene tasks such as washing, toileting and foot care

Product description	Product #	Contracts	
Hygiene showering chair (standard accessories head cushion and safety belt and bed pan holder)	BOC1003-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Carino kit (BOC1003-01, back cushion, seat cushion, bedpan with lid and toilet adapter)	BOC1010-01		



Concerto® and Basic shower trolleys

- Concerto and Basic shower trolley help to facilitate efficient showering and hygiene care routines for dependent patients or residents who need to remain in a lying position
- Promotes correct working posture

Product description	Product #	Contracts	
Shower Trolley Basic (does not include: straight steer, leveling control or tilting function)	BAB5000-WW	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Concerto 1900 mm - hydraulic foot pedal	BAB1000-WW		
Concerto 1900 mm - Electric	BAB1100-US		
Concerto 1600 mm - hydraulic foot pedal	BAB2000-WW	OPEN MARKET	
Concerto 2250 mm - hydraulic foot pedal	BAB3000-WW	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Concerto mattress	8131422-014	OPEN MARKET	
Concerto flush sink stainless steel-P Trap (P Trap is for wall drain)	80015-53		



Carevo® shower trolley

- The Carevo shower trolley is designed to provide a comfortable shower experience for patients and residents who need to remain in a lying position, while enabling the caregiver to work efficiently and ergonomically

Product description	Product #	Contracts	
Shower trolley with adjustable back raising, electric 1900 mm, includes nurse controls and drainage hose	BAC1111-01	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017

Hygiene solutions

Shower systems



Carmina® bariatric hygiene showering chair

- Designed for assisted showering, toileting and hygiene care for plus size patients and residents, Carmina offers a solution for transferring from the bedside, chair or wheelchair to the shower or toilet, and back again
- Carmina is suitable for plus size patients and residents with the mobility level of Albert, Barbara or Carl

Product description	Product #	Contracts	
Hygiene showering chair with foot support	BLB2000-013	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017



Shower Panel

- The Shower and Disinfection Panel provides a controlled delivery system for assisted showering, as well as for the cleaning and disinfection of shower equipment. Assisted showering with secure warm water supply. The wall mounted shower panel is a part of the shower system when used with shower chairs or trolleys. It features a thermostatic mixer with an anti-scalding protection. During showering the water can be closely managed, delivering a gentle water stream at an optimum temperature
- Supports infection prevention with an integrated disinfection system which automatically dilutes the concentrated disinfectant for an easy and secure cleaning and disinfection process

Product description	Product #	Contracts	
Universal Shower Panel with disinfection.	BCA1204-04	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Shower Panel with Flush Tank - Must be combined with Shower Cabinet with toilet or flush sink (080015/080090 series) for Shower Trolley	BCA1214-04		

Venous thromboembolism prevention

Deep Vein Thrombosis solutions



Flowtron® ACS900 pump

- Reduce the risk of VTE in your facility with the Flowtron ACS900 pump, featuring SmartSense™ 2 Automatic Garment Recognition and Compliance Monitoring technology.

Product description	Product #	Contracts	
Flowtron ACS900 pump - standard length tubing	526000-01-P	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Flowtron ACS900 pump - OR length tubing	526000-02-P		



Flowtron® DVT garments

- Built on a uniform (single-chamber bladder) compression philosophy, Flowtron DVT is an extensive range of foot, calf and thigh garments addressing a wide variety of clinical needs and patient types.

Product description	Product #	Contracts	
Small calf garment (10 pairs/box) - fits up to 14" calf circumference	DVT5	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Flowtron standard calf garment (pair) - fits up to 17" calf circumference	DVT10		
Flowtron large calf garment (pair) - fits up to 23" calf circumference	DVT20		
Flowtron standard thigh garment (pair) - fits up to 28" thigh circumference (Case of 10)	DVT30		
Flowtron large thigh garment (pair) - fits up to 35" thigh circumference (Case of 10)	DVT40		
Flowtron bariatric calf garment (pair) - fits up to 32" calf circumference (Case of 10)	DVT60L		
Flowtron standard foot garment - fits women size up to size 9, men up to size 7 (Case of 10)	FG100		
Flowtron large foot garment - fits women size 9½ and larger, men size 7½ and larger (Case of 10)	FG200		
Flowtron calf garment with orange binding (pair) - fits up to 17" calf circumference	L501-M		
Flowtron thigh garment (pair) - fits up to 28" thigh circumference (Case of 10)	L503-M		
Flowtron standard sequential calf garment - fits up to 17" calf circumference (Case of 10)	TRP10		
Flowtron large sequential calf garment - fits up to 23" calf circumference (Case of 10)	TRP20		
Flowtron standard sequential thigh garment - fits up to 28" thigh circumference (Case of 10)	TRP30		
Flowtron large sequential thigh garment - fits up to 35" thigh circumference (Case of 10)	TRP40		
Flowtron bariatric calf garment - fits up to 32" calf circumference (Case of 10)	TRP60L		

Medical beds



Enterprise® 9000X hospital bed

The Enterprise 9000X hospital bed is designed for intensive care and high-dependency patients, offering advanced safety & comfort for both patients & caregivers.

- Enhanced patient and caregiver safety
- Compliance with international standards
- Robust, intuitive, and easy-to-use design
- Quality-focused care solutions
- Weight limit: 550 lbs

A proven choice for hospitals worldwide.

Product description	Product #	Contracts	
Bed frame 9800 — 2-way nurse call with environmental controls, SafeSet & EMR, no handset, standard back rest (curved deck for Bio-Contour)	9814DF101BCUSB	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
IndiGo® stand-alone retrofit unit for Enterprise	INDIJAA		



Citadel® Plus bariatric bed

- The Citadel Plus Bariatric Care System has been designed to designed for dignified bariatric care
- Engineered for high-demand hospital environments, Citadel Plus supports patient recovery while enhancing safety and efficiency. Universal platform design helps transition between various Arjo beds addressing wide range of needs in one familiar platform
- Flexible sizing that adjustable width & length accommodate various patient sizes & ease transport through doorways and elevators
- Weight limit: 1,000 lbs

Product description	Product #	Contracts	
Citadel Plus Bariatric Bed frame — 2-way nurse call, environmental controls, handset, EMR, SafeSet, 150 mm twin castor, foot end, DIN, foot pedal & Power Drive	FX811B3B4AMABB	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017



Citadel® Patient Care System

The Citadel Patient Care System is built on a universal platform, designed to seamlessly transition from non-integrated to fully integrated care.

Integrated for better outcomes

- One-button CPR to activate CPR mode on both mattress and frame instantly
- Automatic pressure adjustment mattress adapts as the head of bed is raised
- Turning is restricted when the head of bed exceeds 30° or when side rails are lowered for enhanced safety
- Citadel's paired bed & mattress are designed together to help address entrapment risks
- Weight limit: 595 lbs

Product description	Product #	Contracts	
Citadel C200 bed frame — 2-way nurse call, environmental controls, EMR, SafeSet, DIN rail & Indigo Drive assist	CX812A3K1AMADO	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
C200 control unit	MATT0038	OPEN MARKET	
C200 mattress	MATT0041		



Total Lift™ Bed

- Designed for early mobilization in critical care that supports upright positioning & progressive weight-bearing with minimal caregiver effort
- Early rehab support that moves patients from lying to full standing to combat immobility-related deconditioning
- Efficient operation that enables safe mobilization using only one or two caregivers
- Weight limit: 425 lbs

Product description	Product #	Contracts	
Total Lift Bed V5 — Trendelenburg, mattress, weight bearing, Tilt Quick Release & Power Drive	VGTLB425V5MWPD	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017

Therapeutic support surfaces


Powered mattresses



AtmosAir® 9000X mattress system

- Open pressurized system redistributes air to maintain consistent pressure levels
- Durable premium cover that last 2.5 times longer than standard fabrics
- Cover is chemical-resistant, antimicrobial & abrasion-resistant
- Infection control with welded seams & fully welded carry handles, included as standard
- Heelguard, sloped heel for added heel protection
- Firm foam outside parameter for safer exit out of the bed
- Weight limit 550 lbs

Product description	Product #	Contracts	
34"x78"mattress with premium cover	633037US	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
35"x80" mattress with premium cover	633038US		
35"x84" mattress with premium cover	633039US		



AtmosAir Velaris®

- A flexible, cost-effective system that replaces most powered surfaces with one adaptable platform -ideal for low to high-risk patients (excluding highest risk patients)
- Dual therapy modes:
 - With pump, alternating pressure & full offloading via AltoVac®
 - Without pump, reactive pressure redistribution with ARM® for consistent comfort
- Dementia Design Accredited product with **1A rating**
- Weight limit 550 lbs

Product description	Product #	Contracts	
78"x34"mattress with premium cover	633026EN	OPEN MARKET	
80"x35"mattress with premium cover	633027EN		
84"x35"mattress with premium cover	633028EN		



Auralis® alternating pressure mattress system

- Advanced support for high-acuity patients at risk of pressure injury
- Dual therapy modes, Self Set Technology™ adjusts pressure in active & reactive modes for personalized care
- Premium cover design with durable, water-resistant & infection-control with welded seams and protective flap
- Engineered for comfort & safety to patients with limited mobility & compromised skin integrity
- Optional Skin IQ® cover for optimal microclimate management for high-at-risk patients
- Weight limit 550 lbs

Product description	Product #	Contracts	
80"x34" - 110 mattress overlay	636T01S-P	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
80"x34" - 200 mattress replacement system	636M02S-P		
Pump	636003US-P		
Auralis seat cushion	636C01S-P	OPEN MARKET	
Skin IQ complete kit: coverlet, fan & power supply. Case of 4	SIQCPST04US	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Skin IQ 365 (reusable) mattress coverlet only	SIQCORUS1UN		



Auralis® Plus mattress system

- Advanced pressure care for plus-size patients up to 1,000 lbs
- Dual therapy with self-adjusting alternating & reactive modes with optional microclimate management via Skin IQ & Heelguard technology
- Dignity & consistency for bariatric patients that aligned with Arjo's therapeutic surfaces for a seamless care environment
- Width-adjustable to support efficient transfers & transport on Citadel Plus bed frame
- Durable, water-resistant & infection-control premium cover
- Weight limit 1,000 lbs

Product description	Product #	Contracts	
Auralis Plus mattress	636B02-P	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Pump	636003US-P		
Optional microclimate management: Skin IQ 1000 mattress kit of coverlet, fan & power supply	SIQCPBAS1US		

Therapeutic support surfaces


Powered mattresses



AtmosAir Velaris® Plus Flex

- AtmosAir Velaris Plus Flex is designed for bariatric patients. It accommodates a 'One Surface Strategy' as it can replace all but the most advanced automatic surfaces with one easy to use, flexible mattress system
- Product features removable bolsters to help during transport
- Velaris pump (works on standard & bariatric Velaris) can be added for active therapy & full offloading
- Note: AtmosAir Velaris is suitable for all care settings and is a Dementia Design Accredited Product (DSDC) with a **1A rating**
- Weight limit 1,000 lbs

Product description	Product #	Contracts
84"x 48"mattress	633903EN	OPEN MARKET
Left side bolster	633935	
Foot bolster	633939	
Right side bolster	633937	
Mattress base/core mattress only for Velaris Plus Flex with side clips for bolster add on	633915	




AtmosAir Velaris® Seat Cushion

- The AtmosAir Velaris Seat Cushion is part of Arjo's adaptable alternating pressure system designed for pressure injury prevention across a wide range of care settings
- Standard size dimensions: 17x17x4" (432x432x100 mm)

Product description	Product #	Contracts
17x 17x4 AtmosAir Velaris standard seat cushion (non-powered)	633018	OPEN MARKET
18x 18x4 AtmosAir Velaris large seat cushion (non-powered)	633019	


Therapeutic seating solutions



Sorrento 2™ motorized chair

- The Sorrento 2 lends itself to changing user needs in multi-user environments, while keeping the patient/resident in mind. It provides postural support and pressure redistribution. The Sorrento 2 tilt-in-space chair is clinically designed for the patient
- Weight limit: 440 lbs


Product description	Product #	Contracts	
Sorrento 2 motorized chair	PN1550	ECAT	SPE2DH-25-D-0017
Optional removable 1" arm covers	PN1573	OPEN MARKET	
Optional leg rest wedge	PN447	ECAT	SPE2DH-25-D-0017
Optional lateral supports	PN432	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Optional leg lateral supports	PN1577	OPEN MARKET	



Phoenix 2™ motorized chair

- The Phoenix 2 is the most supportive and functional chair and has been designed to accommodate spinal conditions that require maximum support for highly dependent patients.
- Weight limit: 440 lbs


Product description	Product #	Contracts	
Phoenix 2 motorized chair	PN1556	ECAT	SPE2DH-25-D-0017
Optional removable arm covers	PN1573	OPEN MARKET	
Optional leg lateral supports	PN1577	OPEN MARKET	
Optional leg rest wedge	PN447	ECAT	SPE2DH-25-D-0017
Optional lateral support	PN432	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017



Sorrento Bariatric™ chair

- The Sorrento Bariatric chair is a tilt-in-space bariatric chair on wheels.
- Weight limit: 650 lbs

Product description	Product #	Contracts	
Chair with remote control powered	PN543	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
1" padded arm covers (25 mm)	PN493	ECAT	SPE2DH-25-D-0017
Removable seat cover	PN479		
Seat belt	PN426	OPEN MARKET	
Leg rest wedge	PN449	ECAT	SPE2DH-25-D-0017
Lateral supports	PN433		




Atlanta 2™ chair

- A robust and durable chair, Atlanta 2 provides a safe and comfortable seat for those with vigorous involuntary movements, while designed to reduce the risk of falling from the chair.
- Weight limit: 350 lbs

Product description	Product #	Contracts	
Atlanta 2 chair-manual with footstool	PN1867	ECAT	SPE2DH-25-D-0017
Atlanta 2 chair-manual with height-adjustable footstool	PN1970	OPEN MARKET	
Atlanta 2 chair-motorized with height-adjustable footstool	PN2165		
Optional horseshoe pillow	PN450		
Optional concave pillow	PN451		
Optional head pillow	PN452		
Optional 1" padded arm covers	PN489	ECAT	SPE2DH-25-D-0017
Optional 2" padded arm covers	PN490		
Optional unpadded arm covers	PN484		

Therapeutic seating solutions



Milano™ chair

- Designed for accessibility, the Milano is perfect for multiple users and can easily be adjusted to fit different patients and residents.
- Weight limit: 350 lbs

Product description	Product #	Contracts	
Chair with manual tilt/recline chair	PN1009	ECAT	SPE2DH-25-D-0017
Motorized chair	PN1375		
Optional lateral support	PN432	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017
Optional horseshoe pillow	PN450	OPEN MARKET	
Optional concave pillow	PN451		
Optional head pillow	PN452		
Optional table tray	PN809	ECAT	SPE2DH-25-D-0017



Sydney GoFlat™ seat

- The Sydney GoFlat seat has been designed to enable the safe and efficient transfer of critical care patients from lying in bed to a fully seated, clinically optimized position offering care professionals new options in early patient mobilization.
- Weight limit: 440 lbs

Product description	Product #	Contracts	
Chair	PN1453	FSS ECAT	36F79725D0231 SPE2DH-25-D-0017



PATIENT HANDLING SOLUTIONS

Wide range of lifts for safe transfer of patients



SLING SOLUTIONS

Wide range of slings for care of patients



LATERAL TRANSFERS SOLUTIONS

Wide range of lateral transfer solutions for safe transfer of patients



HYGIENE SOLUTIONS

Ergonomic, adaptable bathing and showering systems for all levels of mobility



VENOUS THROMBOEMBOLISM PREVENTION

Active compression systems (IPC) for the prevention of DVT



MEDICAL BEDS

Medical and extended care beds, treatment couches, and functional accessories



THERAPEUTIC SUPPORT SURFACES

Therapeutic support surfaces and systems for prevention and treatment



THERAPEUTIC SEATING SOLUTIONS

Therapeutic support chairs that help reduce pressure injuries, prevent sliding, falls and safely increase early mobilization

At Arjo, we believe that empowering movement within healthcare environments is essential to quality care. Our products and solutions are designed to promote a safe and dignified experience through patient handling, medical beds, personal hygiene, disinfection, diagnostics, and the prevention of pressure injuries and venous thromboembolism. With over 6500 people worldwide and 65 years caring for patients and healthcare professionals, we are committed to driving healthier outcomes for people facing mobility challenges.

Arjo Inc., 2000 Center Drive · Hoffman Estates, Illinois 60192 · 1-800-323-1245 · www.arjo.us

