

Evenda Design that connects

A long term-care bed with a feeling of home



The Evenda® bed range supports the needs of long-term care residents and caregivers in a homelike design. The bed range includes three models: Moveno, Comfio and Sereno, which provide customizable solutions depending

on individual resident challenges and facility preferences. Textile paddings, headboards, footboards, and side panels make it possible to design a bed solution that feels homey and comfortable.



Evenda Moveno

Configured to promote independence and out-of-bed mobility, the Evenda Moveno is suited to a more mobile resident:

Footboard with grab bar to support the resident when moving around the bed space

Swing-out Mobility Assist Handle to help the resident get in and out of bed

In closed position the Mobility Assist Handle provides a good support for the resident to reposition in bed

Under-bed light can help guide the resident when getting up at night



Evenda Comfio

Configured to support safe and efficient care for immobile residents:

4 large castors with central brake to facilitate in-bed transportation

Choice of split, telescoping side rails or 3/4 folding rails

Footboard with accessory bar. Mattress pumps can be attached so hooks do not protrude into the bed space

Battery back-up so the bed can be operated even if temporarily unplugged from power source



Evenda Sereno

Designed to provide a safe and recognizable environment for residents living with cognitive challenges:

Simplified hand control for resident. Hand control for caregiver at the foot of the bed

Wider frame and a low height of only 6.7" (17 cm) to provide a sense of security to help reduce the fear of falling out of the bed

Head and footboard design with color contrasts to help make the bed more recognizable

Under-bed light can help guide the resident when getting up at night



Every detail connects to the bigger picture



- 1 **Auto-regression bed articulation** to help reduce resident migration down the bed
- 2 **Mobility Assist Handle with two positions** that can also be folded away
- 3 **Range of hand controls** including simplified hand control for the resident
- 4 **Available in a wide version**
- 5 **Range of headboard & footboard designs** to customize to needs and preferences
- 6 **Optional manual CPR release**
- 7 **Beds with 6 and 8 castors have a low height of only 6.7" (17 cm)** to accommodate residents very prone to falls

Designed to suit different care needs

Evenda can be configured to meet your facility preferences

● Standard

○ Optional



		MOVENO	COMFIO	SERENO
CASTERS	6 casters (17 cm low height)	○		○
	8 casters (17 cm low height)	●		●
	4 casters (12.2"/30.99 cm low height - 125 mm casters) (13.2"/33.53 cm low height - 150 mm casters)		●	
CASTER COVERS INCLUDED		●		●
HEADBOARDS AND FOOTBOARDS See page 4 for all design options	Style A	○	○	○
	Style B	○	○	○
	Style C	●	●	●
	Style D	○	○	○
SIDE RAILS	no side rails	○		○
	50:50 split, telescoping side rails	○	○	○
	40:60 split, telescoping side rails	○	○	○
	3/4 folding side rails	○	○	○
MOBILITY ASSIST HANDLE		○		○
CONTROLS	Standard keypad hand control	●	●	●
	Keypad hand control with tilt	○	○	○
	Simplified hand control	○	○	○
BATTERY BACK-UP		○	●	○
FRAME WIDTH	35.4" (90 cm)	●	●	○
	42" (107 cm)	○		○
CPR RELEASE		○	●	○
UNDER BED LIGHT		●	●	●

Key features

SYSTEM SPECIFICATIONS			
	Moveno	Comfo	Sereno
Overall standard length	82.5"/2100 mm	84.5"/2150 mm	82.5"/2100 mm
Overall width	39"/990 mm	41"/1040 mm	47.5"/1210 mm
Mattress length standard	78"-78.5"/1980-2000 mm	78"-78.5"/1980-2000 mm	78.5"/2000 mm
Mattress width	35.5"/900 mm	35.5"/900 mm	42"/1070 mm
Overall extended length	Position 1: 86.5"/2200 mm, Position 2: 90.5"/2300 mm		
Mattress extended length	Position 1: 82.5"/2100 mm, Position 2: 86.5"/2200 mm		
Overall length short	80"/2030 mm	80"/2030 mm	80"/2030 mm
Top height position	29.5"/750 mm	31.5"/800 or 32.5"/825 mm	29.5"/750 mm
Low height position	6.7"/170 mm	12"/310 or 13"/335 mm	6.7"/170 mm

OVERALL DIMENSIONS	
Backrest length	31.5"/800 mm
Seat section length	9.6"/245 mm
Calf section length	19.6"/500 mm
Thigh section length	11"/285 mm
Safe working load/max. patient weight	507 lbs/430 lbs (230 kg/195 kg)
Head down tilt	12°
Foot down tilt	12°
Backrest angle	70°
Backrest pause	30°, 45°
Thigh angle	40°
Calf angle	0-18°
Ingress protection	IP24
Shock protection	CLASS 1 OR CLASS 2, TYPE B
Power in	100-240V 50/60 Hz, max. 3,9 A
Safety standard	IEC 60601-2-52

ACCESSORIES***
Lifting pole
IV poles
DIN rail with hooks
Plastic deck covers
Battery back-up
Head end wall spacer
Mobility assist handle
Personal item pocket
Urine bottle holder
Side rail cover, telescoping sides
Padded cover for the headboard
Flexible hand control holder

EVENDA IS AVAILABLE IN DIFFERENT DESIGNS*
A Modern and easy to clean**
B Dual color for increased contrast
C Accessory rail for mattress pumps or other accessories
D Grab rail to support the resident

* Must be selected at time of order
 ** Also available with low footboard
 *** Can be retrofitted



At Arjo, we believe that empowering movement within healthcare environments is essential to quality care. Our products and solutions are designed to promote a safe and dignified experience through patient handling, medical beds, personal hygiene, disinfection, diagnostics, and the prevention of pressure injuries and venous thromboembolism. With over 6,000 people worldwide and 60 years caring for patients and healthcare professionals, we are committed to driving healthier outcomes for people facing mobility challenges.

Arjo Inc., 2349 West Lake Street · Addison, Illinois 60101 · 1-800-323-1245 · www.arjo.us

