




UNDERPADS

Dri-Flo[®]

Incontinence protection pads designed to provide air permeability and absorbency which may contribute to microclimate management. Tri-dimensional matting creates air-permeable durability.

Key features

Dri-Flo Underpads are constructed using a three-layer bonding process that integrates the top, core and backsheet of the underpad. This process produces a highly absorbent, breathable underpad that allows patients to receive the full benefits of underlying microclimate management therapies.¹



FEATURE	BENEFIT
Soft, permeable top layer	Allows fluid to pass to the core of the underpad. Keeps moisture away from the patient's skin.
Tri-dimensional matting	Provides tensile strength and durability.
Super absorbent polymers	Draws fluid away from patient's skin and encapsulates it within the core. Absorbs secondary voids reducing associated odors.
Breathable back sheet	Allows underlying microclimate therapies to function properly. ¹

Dri-Flo Underpads:

- Provide air permeability which may aid in keeping the patient's skin dry and comfortable.
- Maintain structural integrity due to their tensile strength. This may decrease the number of underpad changes required.

PRODUCT CODE	DESCRIPTION AND SIZE	CASE PACKAGING
7046	Large Dri-Flo Underpads 23" x 35" (58,4 cm x 88,9 cm)	6 bags of 12 pads (72 pads per case)
47001	Extra Large Dri-Flo Underpads 32" x 36" (81,3 cm x 91,4 cm)	6 bags of 8 pads (48 pads per case)

1. Arjo Data on File. Evaluation of Skin IQ: Cleveland Clinic - 100018399. September 2011.

March 2019. Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. ® and ™ are trademarks belonging to the Arjo group of companies.
© Arjo

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

Arjo UK Ltd - Head Office · Houghton Hall Business Park · Houghton Regis · LU5 5XF
Phone: +44 (0) 8457 342000 · Private customer: +44 (0) 3457 342000 · Fax: +44 (0) 1582 745745 · Email: sales.admin@arjo.com

Arjo UK Ltd - Service and Technical Centre · St Catherine Street · Gloucester · GL1 2SL
Phone: +44 (0) 8456 114114 · Private customer: +44 (0) 3456 114114 · Service fax: +44 (0) 1452 525207 · Email: ukservice@arjo.com

Arjo UK Ltd - Belfast Office · Units 16-18 · 19 Heron Road · Belfast, BT3 9LE · Northern Ireland
Phone: +44 (0) 28 9050 2000 · Fax: +44 (0) 28 9050 2001 · Email: enquiries.belfast@arjo.com

Arjo (Ireland) Ltd · EA House · Damastown Industrial Park · Mulhuddart · Dublin 15 · D15 XWR3 · Republic of Ireland
Phone: +353 (0)1 809 8960 · Fax: +353 (0)1 809 8971 · Email: enquiries.dublin@arjo.com

www.arjo.co.uk

