

# Overbed Table HT110

## Full top, non-tilting



Lightweight, easy to manoeuvre to reduce moving and handling risks. With a slight upwards pressure the table-top can be raised or lowered which then locks into position automatically providing a stable surface.

With infection control in mind, the table has a smooth non-absorbent surfaces making it extremely easy to wipe clean.

Available colours:

- Blue
- Grey

## Specifications

Specifications		
	Maximum height	1115mm
	Minimum height	770mm
	Top size	900mm x 395mm
	Top Material	ABS Plastic
	Column Material	Mild steel frame
	Castors	50mm
	Safe working load	25kg
	Product Weight	20kg

Features		
	Easy Lift Safety, Fully Mechanical Height Adjustment.	•
	Locks in Compression no Release Required for Raising	•

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure injury and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.