



A CALMING SENSORY EXPERIENCE

Nordic Sensi[®] Chair

Care version

Difficulties with communicating and performing everyday tasks can take a toll on elderly people with dementia. In time, residents may develop disruptive behaviours such as aggression, restlessness, confusion and agitation, all of which may disturb other residents and require additional attention from carers.

Helping people with dementia relax through sensory stimulation with automated rocking, soothing music and tactile stimulation can have a beneficial effect on their behaviour and quality of life. The fully-automated Nordic Sensi Chair offers users a unique sensory experience.^{1,2,3}

Addressing individual needs in the care setting

Pre-set rocking programmes combine smooth, quiet motion and soothing music with tactile stimulation to promote physical and mental relaxation, or a refreshing and stimulating experience, depending on the needs of the resident. Each setting is also easily customisable with just the touch of a button.

A choice of optional accessories could add to the user experience. Residents may find that the overall pressure from the weighted blanket combined with the support of the U-shaped cushion and additional positioning accessories provide a sense of security and "cocooning" effect.



POSITIONING KIT



U-SHAPED CUSHION



WEIGHTED FIBRE BLANKET



WEIGHTED BLANKET



TOP PROTECTION WITH AND WITHOUT POCKETS

CUSTOM POSITIONING

Adjust both the angle of the backrest and the rocking motion. The legrest can also be repositioned for additional comfort.

SWING-OUT ARMRESTS

Enables easy access while providing the resident with a sense of security. The armrest covers are available in burgundy, kiwi, fudge and mauve.*

EXTERNAL AUDIO INPUT

Enables personalised music to be played.

CONTRAST

Contrasting colours may help guide the residents to the seat section.



BUILT-IN AUDIO SYSTEM WITH MUSICURE®

According to the Danish Alzheimer's association, music can have a positive and motivating effect on residents living with dementia. MusicCure has been specially composed to aid residents' relaxation.⁴ In addition, the subwoofers in the lower back offer a gentle tactile stimulation.

REMOVABLE INCONTINENCE COVER

Easy to clean and disinfect. Available in burgundy, kiwi, fudge and mauve.*

WHEELED PEDESTAL

Can be used with Arjo patient lift solutions, while allowing for easy manoeuvring and portability.

SOFT IMITATION LEATHER

Easy to clean and disinfect. Available in black or white.

*Colour options may vary by country.

1. Rocking Chair Therapy - a non-pharmacological approach to comfort and wellbeing - Arjo whitepaper 2019
2. Wellness Nordic Relax Chair user feedback - Arjo whitepaper - 2019
3. Novel therapeutic chair benefits residents with dementia - Arjo whitepaper 2019
4. <https://www.alzheimer.dk/viden-om-demens/behandling-pleje-og-omsorg/naar-musik-er-medicin/>

Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. ® and ™ are trademarks belonging to the Arjo group of companies.
© Arjo, 2022

Nordic Sensi® Chair and MusiCure® copyright, trademarks and logos are the intellectual property of Wellness Nordic A/S and MusiCure in Healthcare ApS.

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

Arjo UK Ltd - Head Office · Houghton Hall Business Park · Houghton Regis · LU5 5XF
Phone (local rate): +44 (0) 3457 342000 · Fax: +44 (0) 1582 745745 · Email: customerservice@arjo.com

Arjo UK and Ireland - Service and Technical Centre · Regus Gloucester Docks · North Warehouse · Gloucester · GL1 2EP
Phone (local rate): +44 (0) 3456 114114 · Email: uk.service@arjo.com

Arjo UK Ltd - Belfast Office · 19 Heron Road · Belfast, BT3 9LE · Northern Ireland
Phone: +44 (0) 28 9050 2000 · Fax: +44 (0) 28 9050 2001 · Email: enquiries.belfast@arjo.com

Arjo (Ireland) Ltd · EA House · Damastown Industrial Park · Mulhuddart · Dublin 15 · D15 XWR3 · Republic of Ireland
Phone: +353 (0)1 809 8960 · Fax: +353 (0)1 809 8971 · Email: enquiries.dublin@arjo.com

www.arjo.co.uk

