



High end
customisable
floor bed

ARJO AGED CARE BED

Sorrento Floor Bed

A premium aged care bed offering high level quality and style

Uncompromising style with premium care features

The care bed is central to the daily wellbeing of the resident and is designed to provide a secure environment to optimise therapeutic rest and recuperation in style.



Fully customisable to suit any decor

The Sorrento Floor Bed is designed to reduce the risk of harm from falls whilst allowing you to provide a home-style bed for your residents and peace of mind for their families and caregivers.

The Sorrento Floor Bed has all the features you expect from a home style care bed with a range of advanced features. This floor bed showcases a perfect balance of design, functionality and comfort.



Sorrento Floor Bed

4-Section fully electric mattress platform



Premium Tuscany Panel in Bottega Oak featured



WALL BUMPER
(Standard)



BI-CASTER BRAKE
Four 75mm braking castors
and four free-wheeling
(Standard)




**BUILT-IN BASE
EXTENSION**
An extra 195mm to the
mattress base length
(Standard)



**ADVANCED CONTROL
PANEL WITH CPR**
(Standard)



**LIGHTWEIGHT HAND
CONTROL**
(Standard)



IV POLE
(Optional)



EGRESS RAIL
(Optional)



**THREE QUARTER
LENGTH RAIL**
(Optional)



SELF HELP POLE
(Optional)



Floor Height



Semi Chair



Full Chair

Sorrento Floor Bed your style

With a range wide range of design options, the Sorrento Floor Bed can be tailored to suit any project, decor or style.

Choose from an extensive range of panel options, colours, accessories, furniture and therapeutic surfaces that suit the environment and clinical needs of the resident.

1 Choose your head and foot panel profile



Tuscany - 01 *★Premium option*



Stockholm - 02



Adelaide - 03



Monaco - 04

2 Choose your material colour

Available in both vertical or horizontal grain

Woodmatt - Smooth textured



01 Nordiac Oak
02 Angora Oak
03 Castino Beech
04 Tasmanian Oak
05 Antico Oak



06 Natio Walnut
07 Havana Oak
08 Empire Oak
09 Perugian Walnut
23 Costal Oak



24 Prime Oak
25 Florentine Walnut
26 Bottega Oak
27 Estella Oak

Matt - Smooth



10 Gesso Lini
11 Tuross Oak
12 Natural Oak
13 Marni Lini
14 Soft Walnut



15 Maison Oak
16 Rocco Lini
17 Artisan Oak
18 Belgian Oak
19 Cava Lini



20 Truffle Lini
21 Shannon Oak
22 Castagna Relativo
28 Beech Select



Premium modern
- Tuscany in Bottega Oak (Premium Panel)



Classic natural
- Stockholm in Beech Select (Standard Panel)


3 Choose your accessories and furniture




Egress rail
LAU-HGR-L-BR- Left
LAU-HGR-R-BR- Right



Three quarter length rail
LAU-FSR-CL-BR



IV pole
LAU-IV110



Self Help pole
LAU-SP002



Sorrento Floor Bed
assist bar (pair)
G420-88ASB

Full range of bed side lockers, cabinets, wardrobes & over bed table tables available on request

★★ Available in various colours



Tallboy, bedside cabinet, and wardrobe - with finger pull**



Tallboy, bedside cabinet, and wardrobe - with handles**



U Base over bed table**
Thermo-laminated finish with raised edge



U Base over bed table**
Laminated top with ABS edging

Product specifications

	SINGLE	KING SINGLE
Overall Length 211.5cm	•	•
Overall Length Extended 231cm	•	•
Overall Width 92cm	•	-
Overall Width 109cm	-	•
Overall Width (rails fitted) 102.5cm	•	-
Overall Width (rails fitted) 120cm	-	•
Mattress Surface Length (between panels) 207cm	•	-
Mattress Surface Length (between panels) 202.5cm	-	•
Bed Platform Height Range 16.0cm to 69.0cm (+/-5mm)	•	•
Back Rest Angle Adjustment 0° to 70°	•	•
Knee Bend Angle Adjustment 0° to 35°	•	•
Safe working load 250kg	•	•
Maximum Patient Weight 215kg	•	•
Electrically adjusted height range	•	•
Trendelenberg head and foot down tilt	•	•
Four sectioned electrically operated mattress platform with independent adjustment of backrest and knee section	•	•
Four soft contoured mattress retainer handles	•	•
Integrated foot end mattress base extension. Extends the bed by 195mm	•	•
OPTIONAL ACCESSORIES		
Egress folding side rails	•	•
Three quarter length folding side rails	•	•
Sorrento FB Assist Bar	•	•
USB capability	•	•
Range of head and foot panels for customisation	•	•

Designed and manufactured in compliance with International Standard IEC 60601-2-52

**For further information please contact 1800 072 040
or email: customerservice-au@arjo.com**

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. Our products and solutions ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, effective prevention of pressure injuries and venous thromboembolism and helping professionals across care environments raise the standard of safe and dignified care. Everything we do, we do with people in mind.

Regional Head Office · Arjo Australia Pty Ltd · Level 3 Building B, 11 Talavera Road · Macquarie Park NSW 2113 · Australia · 1800 072 040

www.arjo.com.au

arjo

Mattress options



ComfortFLx

A high quality multi layered hospital grade foam mattress, developed to immerse the patient in the mattress and provide pressure redistribution.



Air Lo 8 king single alternating

Mattress replacement alternating system designed for king single medical beds with heel guard technology.



AtmosAir Velaris

Hybrid mattress replacement system. Step up to alternating therapy as required by addition of the Velaris pump.



Alpha Active 4

Alpha Active 4 Mattress Replacement System for high risk of pressure injury on single bed size.